

# Processing powder data

JPW

Spring 2002

## Abstract

The real aim of this project was for the current author to have a go at literate programming. As a side effect it is hoped that a documented set of data reduction routines for ID31 will be produced. These routines should not only be correct, but also include documentation to inform the user what they do and how they do it, in excruciating detail. They should fully replace the programs used for BM16, and hopefully be understandable and extensible enough to accomodate any future developments.

## Contents

|          |   |           |
|----------|---|-----------|
| <b>1</b> | <b>Introduction</b>                               | <b>1</b>  |
| <b>2</b> | <b>Specifications</b>                             | <b>1</b>  |
| <b>3</b> | <b>Program structure</b>                          | <b>2</b>  |
| <b>4</b> | <b>Reading SPEC files</b>                         | <b>2</b>  |
| 4.1      | Some variables to hang on to . . . . .            | 3         |
| 4.2      | File opening . . . . .                            | 3         |
| 4.3      | Finding a specific scan . . . . .                 | 4         |
| 4.4      | Reading and interpreting file headers . . . . .   | 5         |
| 4.5      | Reading lines of data . . . . .                   | 10        |
| 4.6      | A module for all of that specfile stuff . . . . . | 10        |
| 4.7      | Test program . . . . .                            | 11        |
| 4.8      | Utility programs . . . . .                        | 13        |
| <b>5</b> | <b>Rebinning the data</b>                         | <b>16</b> |
| 5.1      | Useful variables . . . . .                        | 16        |
| 5.2      | Some functions for binning . . . . .              | 19        |
| 5.3      | Initialisation . . . . .                          | 22        |
| 5.4      | The binning algorithm . . . . .                   | 23        |
| 5.5      | A rebinning module . . . . .                      | 27        |
| 5.6      | Converting to another unit . . . . .              | 29        |
| 5.7      | Test programs . . . . .                           | 31        |
| 5.8      | Median filtering . . . . .                        | 40        |
| <b>6</b> | <b>Summation and normalisation</b>                | <b>41</b> |
| 6.1      | Combining data from multiple detectors . . . . .  | 41        |
| 6.2      | Taking the median of the channels . . . . .       | 51        |
| 6.3      | Saving the ascan array . . . . .                  | 52        |
| 6.4      | Superzapem . . . . .                              | 53        |
| 6.5      | Combining data from different scans . . . . .     | 54        |
| 6.6      | Pattern scaling . . . . .                         | 54        |
| 6.7      | A module for the summation stuff . . . . .        | 55        |

|           |  |            |
|-----------|--|------------|
| <b>7</b>  | <b>Detector offset calibration</b>               | <b>56</b>  |
| <b>8</b>  | <b>Writing output files</b>                      | <b>60</b>  |
| 8.1       | Various file formats . . . . .                   | 60         |
| 8.2       | Diagnostic plots . . . . .                       | 64         |
| 8.3       | epf and inp formats . . . . .                    | 68         |
| 8.4       | Logfiles . . . . .                               | 69         |
| 8.5       | Utils . . . . .                                  | 70         |
| 8.6       | Output files module . . . . .                    | 71         |
| <b>9</b>  | <b>Driver programs</b>                           | <b>72</b>  |
| 9.1       | Specifying user options . . . . .                | 72         |
| 9.2       | Get monitor col . . . . .                        | 81         |
| 9.3       | helpmsg . . . . .                                | 87         |
| 9.4       | Something like binit . . . . .                   | 88         |
| 9.5       | Something like binem . . . . .                   | 90         |
| 9.6       | Checking scans for consistency . . . . .         | 92         |
| 9.7       | Merging inp files . . . . .                      | 95         |
| 9.8       | AOB . . . . .                                    | 97         |
| <b>10</b> | <b>Tidying up</b>                                | <b>97</b>  |
| <b>11</b> | <b>Peak Fitting</b>                              | <b>98</b>  |
| 11.1      | Introduction . . . . .                           | 98         |
| 11.2      | Peak function . . . . .                          | 98         |
| 11.3      | Background function . . . . .                    | 99         |
| 11.4      | Defining the problem . . . . .                   | 100        |
| 11.5      | Marquardt Damping . . . . .                      | 101        |
| 11.6      | Modules . . . . .                                | 101        |
| 11.7      | Forming the least squares matrix . . . . .       | 102        |
| 11.8      | Parameter estimation . . . . .                   | 104        |
| 11.9      | Matrix inversion . . . . .                       | 105        |
| 11.10     | Driver program . . . . .                         | 106        |
| <b>12</b> | <b>Sifit</b>                                     | <b>107</b> |
| <b>13</b> | <b>Code building</b>                             | <b>110</b> |
| <b>14</b> | <b>Indices</b>                                   | <b>111</b> |
| <b>15</b> | <b>Appendix 1</b>                                | <b>114</b> |
| <b>16</b> | <b>One Page Guide</b>                            | <b>134</b> |
| <b>17</b> | <b>Pseudo-Random number generator</b>            | <b>136</b> |
| <b>18</b> | <b>Sorting algorithm from netlib for medians</b> | <b>136</b> |

## 1 Introduction

Powder diffraction data were traditionally recorded by moving a detector to a particular angle, recording the number of counts produced in a given time, and then moving on to measure at the next angle of interest. The resulting data are then a set of angular positions and counts, and perhaps esds (or statistical weights). Conventionally these data are then analysed by programs which insist on constant angular step sizes. Collecting data in this format by step scanning has a disadvantage, the detector arm has to keep stopping and starting, which is fairly time consuming. Also, using more than one detector makes it very

difficult to ensure data from different channels will differ by an integer number of steps. The method used at ESRF is to scan the diffractometer at constant speed and record the counts arriving and angular positions as a function of time and then to take that (large) dataset and process it into a conventional format for the user. A series of scans can also be combined together, allowing greater flexibility in experimental design. This document is intended to illuminate the mysteries of the data processing carried out.

## 2 Specifications

The program(s) will need to do the following things:

- Put the data on a constant step sized two theta scale
- Combine data from multiple scans (and SPEC files even)
- Check scans for consistency ( $R_{merge}$ ) when combining them together
- Determine detector offsets
- Determine detector efficiency corrections
- Apply any other corrections, for example flat plate stuff
- Produce two dimensional datasets versus external variables (eg temperatures)
- Identify and deal with problems (glitches, zingers, shclose etc etc)
- Demonstrate some sort of correctness
- Carry out the tasks efficiently (be fast)
- Doubtless other things currently undreamed of... (be vaguely extensible)

The data arrives in files, called SPEC files, which contain motor positions, counts, times and perhaps temperatures. Getting from this information to the final constant step sized single histogram involves some extra information. The angular offset and efficiency of each detector need to be known or determined. The counts need to be rebinned onto a commensurate two theta scale and corrected for detector efficiency and summed together. Finally some error bars for the data points need to be determined and an overall scaling can be applied to the histogram.

## 3 Program structure

Having one large program to do all of these jobs is likely to result in something which is as fast as possible, but might lack flexibility in the future. The intention is to produce a series of smaller programs and subroutines which will work together and can also be combined into monolithic beasts. (FIXME) This document is initially going to be a bit of a mess, but the hope is that by being forced to improve and think about this as a readable account of how the binning is done, the programs will benefit. We can delete that sentence once we're finished!

Fortran 90 is being used throughout, with no apologies from the author. Knowledge of that language would be handy for anyone modifying these routines in the future. Ideally there will be enough text and comments in here to make things fairly self explanatory. The code should be standard conforming (at least as far as the compilers can enforce the standard). Porting to fortran 77 probably just means doing something about the allocatable arrays and removing some syntax. Adding some fortran 95 bits would be desirable if we can get hold of a compiler that generates parallel code. The code is in something of a condensed format in places, with multiple statements on a single line. Given the amount of explanation in the accompanying text it seemed reasonable to save space. A global find and replace of ; with carriage return plus six spaces would probably fix this. Real numbers are universally kind=8 (double precision) and integers are single precision. Results have been found to be numerically identical on both windows

and linux based systems, suggesting that any floating point problems remain fairly subtle. (FIXME) This section wants an overview of the programs from the users perspective. Or maybe a user manual section for the programs in a userio section?

We'll begin with a section about SPEC files and some routines to handle them. Then move onto rebinning, describing what is actually recorded by the instrument during an experiment and how to convert that into something we are more used to. This will culminate in a program which can take a single scan out of a SPEC file and write out binned files for each of the detectors separately.

Once the data are on convenient  $2\theta$  scales they can be added together (more easily). Combining the data from different channels and different scans means just summing them together, with appropriate corrections for detector efficiencies (and other nastiness). The determination of the detector efficiency is logically carried out here, along with some checks for compatibility of the scans.

Some statistics, like  $R_{merge}$  will be introduced along with the calculation of the statistical weights for the final histograms. Unfortunately the  $R_{merge}$  calculation for a particular scan requires you to already know how the total looks, so that needs access to both individual scans and the summed total - can be done on two passes of the SPEC file or by processing the scans separately and then together.

The only thing obviously missing is the calculation of the detector offsets. Due to the design of the instrument these should never change - so this gets a section to itself and a bit of a standalone method and program.

Peak fitting is so ubiquitous that some routines might be added to the package to do this - mainly for strain scanning experiments, otherwise it gets a bit too involved.

## 4 Reading SPEC files

Fortran code for doing this probably exists already somewhere at ESRF, but we didn't find it. Something is knocked together here which is specific to BM16/ID31 files. It would probably make a lot more sense to have used the specfile library which is designed for c programmers. Probably this could be wrapped up somehow but it looked a bit daunting. For now, some fortran specific stuff is used, but we might come back to doing it properly later. A cursory examination of the specfile library suggests it tries to index all of the scans in the file before doing anything useful, which is something we want to avoid. When someone comes along with a 100MB file containing a thousand scans that might take a long time, so we aim to traverse a specfile in a single pass, collecting the information we want as we go along. File IO is likely to be the time critical step in the whole operation although some profiling should be done if things aren't fast to begin with.

### 4.1 Some variables to hang on to

Some shared information amongst the various specfile routines is needed, so we'll make a list of what we want here.

```

< SPECVARS 3 > ≡
    integer(kind=4) :: iunit, iscan, ispecerr, ncolumns
    integer(kind=4), parameter :: LINELENGTH=512, WORDLENGTH=100
    integer(kind=4), parameter :: NWORDS=60, NLINES=50
    integer(kind=4), dimension(NLINES) :: headerwords
    real(kind=8) :: scanstart, scanend
    real(kind=8), dimension(NWORDS, NLINES) :: wordvalues
    real(kind=8), dimension(NWORDS) :: Q
    character(len=6) :: lnfmt
    character(len=LINELENGTH) :: filnam, line
    character(len=LINELENGTH) :: specsfilename, specsdate
    character(len=WORDLENGTH) :: scantype
    character(len=WORDLENGTH), dimension(NWORDS, NLINES) :: words
    character(len=WORDLENGTH), dimension(NWORDS) :: columnlabels
    integer(kind=4) :: i0, i9, id, isp, im
    character(len=WORDLENGTH) :: FIRSTDET, LASTDET, TWOTTH
    character(len=WORDLENGTH) :: MONITORCOL="Monitor"
    logical topofscan
    data lnfmt /'(a512)'/
    data iunit /15/

```

◇

Fragment referenced in 11a.

Those variables are intended to contain the following information: **iunit** is the fortran unit number which the currently opened file is attached to, **iscan** is the current scan, **ispecerr** indicates something is going wrong. **ncolumns** is the number of columns of data. The parameters **LINELENGTH** and **WORDLENGTH** specify the length of some character strings and **NWORDS** and **NLINES** dimension arrays for us. **LINELENGTH** is the maximum length of a line (distance between carriage returns) which we allow. **specsfilename** will hold the name that SPEC was using for the name of the file and **specsdate** will hold the date the file inwas created on (when the experiment was originally carried out). This information can be inserted into datafile titles if nothing else is provided. We think (but are not sure) that longer lines will be truncated if read into a string which is too short.

The integer array **headerwords** makes a note of how many words are found on the #O lines and **NLINES** is the maximum number of #O lines we can handle. Those words will be stored in the array **words**. **wordvalues** will hold the numbers corresponding to the words held in **words**. **filnam** is the name of the SPEC file and **line** is the contents of the most recently read line. The scan type (turboscan, ascan, whatever) will be stored in **scantype**. Finally the column labels (from the #L) line are stored in **columnlabels**.

## 4.2 File opening

A file is needed before anything can be done with it so we start by trying to get hold of one. At this point we'll need to know the filename which the user is going to have to specify if it wasn't already supplied. Calling the routine when a file is already open causes that file to be closed and the current one to be opened (so only one file will ever be open at once from here). Normally the filename will be supplied from the command line, and the main program should ensure it is already filled in for us.

```

⟨ getfile 4 ⟩ ≡
    subroutine getfile
    logical :: od
    ! See if there is already a file open, if there is then close it
    inquire(unit=iunit,opened=od)
    if(od)close(iunit)
    ! filnam must already be filled in, or we'll try to open garbage
1    open(unit=iunit,file=filnam,form='formatted',                &
    & access='sequential', status='old',iostat=ispecerr)
    if(ispecerr.ne.0) then
    write(*,'(3a,i6)') 'Problem with file ',filnam(1:len_trim(filnam)),&
    & ' iostat=',ispecerr
    write(*,'(a)') 'Please supply a valid SPEC file name'
    read(*,lnfmt)filnam
    goto 1
    endif
    topofscan=.false.
    ! Supply defaults in case dumb specfile has no #L
    ncolums=14 ; columnlabels(1)='2_theta'; columnlabels(2)='Epoch'
    columnlabels(3)='Seconds'; columnlabels(4)='MA0'
    columnlabels(5)='MA1'; columnlabels(6)='MA2'
    columnlabels(7)='MA3'; columnlabels(8)='MA4'
    columnlabels(9)='MA5'; columnlabels(10)='MA6'
    columnlabels(11)='MA7'; columnlabels(12)='MA8'
    columnlabels(13)='Monitor'; columnlabels(14)='Fluo det'
    i0=ichar('0');i9=ichar('9')
    id=ichar('.');isp=ichar(' ');im=ichar('-')
    return
    end subroutine getfile

```

Fragment referenced in 11a.

Default values are provided for the number of data columns, **ncolums**, and their labels, in **columnlabels**. The values **i0**, **i9**, **id**, **is** and **im** are used in the **rdnums** routine, and initialised here as we expect this will always be called before we try to read any numbers.

This will need to be modified or replaced if it sits behind a gui, or perhaps it should just send a unit number to **useroptions**?

### 4.3 Finding a specific scan

Having gotten a file opened we'll want to read it. Usually we'll go after a particular scan so here's a routine for doing just that, provided the scans are in order in the SPEC file. It positions the file at the top of the data beginning of the scan. Any header information will already be read in by the **readheader** routine, so there is no need to call **readheader** as well as **findscan**. (That would actually cause problems). There is an implicit assumption that scans will always be in the file in ascending order and that we will always access them in ascending order. If this is not the case we can add some kind of SPEC file preprocessor which puts the scans in order, which is more suited to some kind of script.

```

⟨ findscan 5 ⟩ ≡
    subroutine findscan(n)
    integer(kind=4),intent(in)::n
    character(len=7) :: rd ! enough space for yes, no and unknown
    inquire(unit=iunit,read=rd) ! Can we read the file?
    if(rd(1:3).ne.'YES') call getfile ! If not go open it
    if(iscan.eq.n .and. topofscan) return
    if(iscan.gt.n) rewind(n) !
2    call readheader ! in case top of file
    if(iscan.ne.n)then ! work through file to find #
1    read(iunit,lnfmt,err=10,end=10)line
        if(line(1:1).eq.'#') goto 2
        go to 1
    else
    ! Scan found - check scantype elsewhere - not a job for specfile module
    return
    endif
10   write(*,*)
        write(*,'(a,i5,a)')'Scan ',n,' not found in file '//
        & filnam(1:len_trim(filnam))
        ispecerr=-2
    return
    end subroutine findscan

```

Fragment referenced in 11a.

The routine returns sets **ispecerr** to -2 if the scan is not found in the file. This is to distinguish from the end of a scan, which is currently caught by setting **ispecerr** to -1. Note that the routine uses the **readheader** routine to interpret any header lines found, rather than trying to read them itself. The logical **topofscan** is set by **readheader** when it gets to a data line after a set of header lines. These complications are just to prevent us ever needing to read the same line twice.

#### 4.4 Reading and interpreting file headers

It will be necessary to read and intelligently interpret the various bits of header information which go with each scan. At the very least we need to check the kind of scan and the header information for each of the columns. If the detector offsets ever get stored in the specfile this is a good time to dump them on a temp.res file or something. Probably strain scanning experiments will want various motor positions recorded somewhere, so we just pull everything possible out here.

For interpreting the header we need some information about the meaning of the various # labels in the header. This has been alluded to earlier, but here is a definitive list.

- #F filename
- #E A number - not sure what it is!
- #D The date
- #C The creator of the file
- #On Index of the labels which will correspond to following #P cards
- #Sn Scan number and information
- #G Diffractometer geometry for spec fourc mode (see link below)
- #Q Four circle variables. h,k,l = q[0,1,2], wavelength=q[3]
- #Pn values corresponding to motor names in #O
- #N number of columns in output data

- #L names of columns in output data (should be #N columns)
- blank, the data themselves
- blank line - end of a scan?
- scan ended by control-c ? Did I imagine this ?

For the rest of #Q see: [http://www.certif.com/spec\\_manual/fourc\\_4.9.html](http://www.certif.com/spec_manual/fourc_4.9.html)  
 The **readheader** routine interprets these lines.

```

⟨ readheader 6 ⟩ ≡
  subroutine readheader
  character(len=1) :: letter
  character(len=128) :: motor
  integer(kind=4) :: i,j
  real(kind=8) :: a
  integer(kind=4),parameter :: four=4
  if(line(1:1).eq.'#') goto 2 ! if already on header, don't skip
1  read(iunit,lnfmt,end=100)line
2  if(line(1:1).eq.'#') then; letter=line(2:2); topofscan=.true.
    select case((letter))
      case('O')
        read(line(3:3),'(i1)')i
        call split(line(4:LINELENGTH),words(:,i+1),WORDLENGTH,NWORDS,j)
        headerwords(i+1)=j
      case('P')
        read(line(3:3),'(i1)')i
        call rdnums(four,NWORDS,wordvalues(:,i+1))
      case('Q')
! Copy the Q from the header into our specfile module
        read(line(3:LINELENGTH),*,end=1, err=1) Q
        case('N')
          read(line(3:len_trim(line)),*,end=1)ncolumns
        case('L') ! signals end of header, bug out here !
          call split(line(3:LINELENGTH),columnlabels,WORDLENGTH,NWORDS,i)
          if(ncolumns.ne.i)
& write(*,'(a,i5)')'error reading header for scan',iscan
        case('S')
          read(line(3:len_trim(line)),*,end=1)iscan,scantype,motor,
& scanstart, scanend
        case('F')
          specsfilename=line(3:len_trim(line)) ! get original filename
        case('D')
          read(line(3:len_trim(line)),'(a256)',end=1)specsdate
        case('C'); continue ! comments
        case('G'); continue ! no idea - always zero
        case('E'); continue ! epoch - we don't care what it was for now.
        case default
        end select
      else ! Not a # line, so assume end of header
        read(line,*,err=1,end=1)a ! catch blank lines in header
        topofscan=.true.          !!! Must be able to read a number !
        return                    !!! escapes from routine here !!!!!
    endif
    goto 1
100 ispecerr=-1; return; end subroutine readheader

```

Fragment referenced in 11a.



A short routine to split a line into a series of words was needed in the **readheader** routine and might generally come in useful, so here it is. Note that the gap between words has to be at least two spaces to allow motor names which contain a space.

```

< split 7a > ≡
      subroutine split(instring,outstrings,lenout,n,i)
      integer(kind=4),intent(in) :: lenout,n
      character(len=*),intent(in) :: instring
      character(len=lenout),dimension(n),intent(out) :: outstrings
      integer(kind=4),intent(out) :: i
      integer(kind=4) :: j,k,l
      j=1; k=1
      do i=1,len_trim(instring) ! hope len > len_trim or array overstep
      if(instring(i:i+1).ne.' ')then
        outstrings(j)(k:k)=instring(i:i)
        if(k.ne.1 .and. k.lt.lenout) k=k+1
        if(instring(i:i).ne.' ' .and.k.eq.1)k=k+1
      else
        if(k.gt.1)then ; do l=k,lenout
          outstrings(j)(l:l)=' ' ! blank pad end of string
        enddo ; j=j+1; k=1; if(j.gt.n)exit ; endif
      endif
      enddo
      ! Blank pad last string if necessary
      if(k.gt.1)then
        do l=k,lenout;outstrings(j)(l:l)=' ';enddo
      endif
      i=j; return; end subroutine split

```

Fragment referenced in 11a.

A quite involded routine which was needed in **readheader**, to convert a string containing a space separated list of numbers into an array. This same routine is used from the **getdata** routine. It is designed to be significantly faster than a simplistic fortran **read(line,\*) vars**, at least on some platforms. Briefly, it parses along the string, ignoring any leading white space and then inteprets the characters it finds into and array of numbers. The only characters a line is allowed to contain (currently) are therefore 1,2,3,4,5,6,7,8,9,0,., and -. Anything else causes **ispecerr** to be set to -1, which normally signals the end of the file (FIXME - add E formats?).

```

< rdnums 7b > ≡
      subroutine rdnums(ic,n,values)
      ! Placed here in a subroutine in case of formatting or error handling problems
      integer(kind=4),intent(in) :: n, ic
      real(kind=8),dimension(n),intent(inout) :: values
      if(line(ic:ic).eq.'#')then; ispecerr=-1; return; endif
      read(line(ic:len_trim(line)),*,err=10,end=20)values(1:ncolumns)
      goto 100
10    write(*,*)'Error reading line, looking for ',ncolumns,' values'
      write(*,*)line(ic:len_trim(line))
      write(*,*)values
20    continue
100   return
      end subroutine rdnums

```

Fragment referenced in 11a.

$\langle junk1\ 8 \rangle \equiv$

```

integer(kind=4) :: i, j, k, l; real(kind=8) :: t
integer(kind=4) :: ii0,ii9,iid,iisp,iim
logical :: pastdot, neg , innum
k=1; t=0.1d0; l=0 ; values=0.0d0
ii0=ichar('0');ii9=ichar('9');iisp=ichar(' ');iid=ichar('.');iim=ichar('-')
innum=.false. ; pastdot=.false. ; neg=.false.
do i=ic,LINELENGTH ! strip leading spaces
  if(ichar(line(i:i)).ne.iisp)then; l=i; exit ; endif; enddo
if(l.eq.0)then; ispecerr=-1 ; return ; endif
do i=1,LINELENGTH ! instring is always line
  j=ichar(line(i:i))
  if(j.le.ii9 .and. j.gt.ii0) then ! digit
    innum=.true.
    if(pastdot)then
      values(k)=values(k) + real(j-ii0,8)*t ; t=t/10.0d0
    else
      values(k)=values(k)*10.0d0 + real(j-ii0,8)
    endif; cycle ;endif  !!! next char !!!
  if(j.eq.iisp)then ! space
    if(.not.innum)cycle  !!! next char !!!
    if(neg)values(k)=-values(k)
    if(k.lt.n) then
      k=k+1
    else ; return;
    endif
    t=0.1d0; pastdot=.false.; neg=.false.; innum=.false.
    cycle;
  endif  !!! next char !!!
  if(j.eq.ii0)then ; innum=.true.
    if(pastdot) then ; t=t/10.0d0 ; else
      values(k)=values(k)*10.0d0 ; endif ; cycle; endif
  if(j.eq.iid)then; innum=.true.; pastdot=.TRUE. ; cycle ; endif
  if(j.eq.iim)then; innum=.true.; neg=.TRUE. ; cycle ; endif
  if(line(ic:ic).eq.'#')then; ispecerr=-1; return; endif
  if(line(i:i).eq.'e' .or. line(i:i).eq.'E')then
    values=0.0d0 ! give up and use fortran read !
!   write(*,*)'About to call fortran read'
    read(line(ic:len_trim(line)),err=10,end=20)values(1:ncolumns)
    goto 100 ! bug out after successful 'e' format read
  endif
! error reading line
  write(*,*)'Could not interpret character... line was:'
  write(*, '(a'))line(1:len_trim(line))
  write(*, '(a,i4,a,i4)')line(ic:ic),j,'at position',ic
!   write(*,*)'ii0',ii0,'ichar(0)=' ,ichar('0'),'ii9','=',ichar('9')
10  continue
enddo
20  continue
100  return; end subroutine rdnums ◇

```

Fragment never referenced.

⟨junk 9a⟩ ≡

```

10      write(*,'(a)') 'rdnums couldn't interpret the line'
      write(*,*) 'ic=',ic,'n=',n,'j=',j
      write(*,'(a)') line(ic:len_trim(line))
      values=0.0d0
      ispecerr=-1;
      return ! Give error message if bad character
! end of line reached - OK as we were able to get line from file
20      return ! with hopefully first few values filled in
      enddo; return; end subroutine rdnums

```

Fragment never referenced.

This horribly complicated routine might later get replaced with something else. It gave a factor of more than 2 (55 seconds to 18 seconds) for reading a 50MB file. To make head or tail of what it is doing you would probably want to expand all the semicolons to put the code on separate lines. It has been tested and it thought to be correct however - it's fairly straightforward to watch under a debugger to see if it reads a string correctly. Note that **line** should already be filled in, this routine reads whatever is in that variable.

Obtaining the information from the file headers, indexed by keyword requires a small function which sifts through the headers which have been read in and then gets the corresponding value. This is to be called after the header is read in, if required. It will be able to supply the initial *2θ* position and the positions of any motors of interest (like X trans, Y trans, Chi etc).

⟨getheadervalue 9b⟩ ≡

```

real(kind=8) function getheadervalue(string)
character(len=*),intent(in) :: string
integer(kind=4) :: i, j, k
k=len(string); getheadervalue=0.0
do i=1,NWORDS; do j=1,NLINES
  if(string(1:k).eq.words(i,j)(1:k))then
    getheadervalue=wordvalues(i,j); return
  endif
enddo; enddo
return; end function getheadervalue

```

Fragment referenced in 11a.

Usage is then simply **variable = getheadervalue('X translation')**. We might add a logfile option to binem to dump out selected motors - will need to do something quite general about spaces in motor names here, they cannot be supplied as command line arguments without quoting.

Understanding the contents of the data columns when they are read in means we need to be able to ask for the column label of a particular counter. The next routine goes through the **columnlabels** array which was filled in during **readheader** to work out which column a particular label corresponds to. It returns -1 if the column is not there (use for example to work out if we have recorded a particular counter value).

```

< whichcolumn 10a > ≡
    integer(kind=4) function whichcolumn(string)
    ! Interprets the #L line information
    character(len=*),intent(in) :: string
    integer(kind=4) :: i, j
    j=len(string) ; whichcolumn=-1
    do i=1,NWORDS
        if(string(1:j).eq.columnlabels(i)(1:j))then
            whichcolumn=i ; return
        endif
    enddo
    return; end function whichcolumn

```

Fragment referenced in 11a.

## 4.5 Reading lines of data

At last we come to the most interesting thing, which is to read in a line of data. Arguments to a subroutine are used to pass the information back. The calling routine is going to have to specify how many numbers to read and supply an array to put them into. All of the glitch elimination, zinger elimination, detection of closed shutters and whatnot is going to happen elsewhere - here we only send back some numbers and set **ispecerr** to a non zero value if something went wrong. The caller should check this value and any data returned from here. A futuristic situation would be to have a version of this routine which can pass back the data as it arrives from the instrumentation, along some kind of unix pipe thing. A quick test on the amber cluster (linux computers at ESRF) showed that if a file is appended with **cat && file** then a fortran program does indeed see the extra lines even if they appear after it starts running. However, if the top of the file is modified (so the fortran program would lose its place) then it seems to work off of a copy of the file in the state it was in when it was opened. These details can be worked out as and when we try for an online binning program.

```

< getdata 10b > ≡
    subroutine getdata(a,n)
    ! A is an array, dimensioned by the number of data items to be expected
    integer(kind=4),intent(in) :: n
    real(kind=8),dimension(n),intent(out) :: a
    integer(kind=4), parameter :: one = 1
    if(topofscan)then; topofscan=.false. ; else
        read(iunit,lnfmt,err=10,end=10)line           ! goes via line
        topofscan=.false.
    endif
    call rdnums(one,n,a)                               ! reads from line
    if(line(1:1).eq.'#')then; ispecerr=-1 ;goto 100 ; endif
    go to 100
10    ispecerr=-1 ! End of file
100   return ; end subroutine getdata

```

Fragment referenced in 11a.

The **getdata** routine is just calling the **rdnums** routine to interpret the line of numbers. This assumes all numbers can sensibly be read into a floating point format (so integers are converted to reals).

## 4.6 A module for all of that specfile stuff

Just placing all the variables and subroutines together in a module for easy access elsewhere. Apart from trying out a fortran 90 feature the author was curious about, this is supposed to be an aid to structured programming.

```

⟨ specfiles 11a ⟩ ≡
    module specfiles
    ⟨ SPECVARS 3 ⟩
    contains
    ⟨ getfile 4 ⟩
    ⟨ findscan 5 ⟩
    ⟨ getdata 10b ⟩
    ⟨ readheader 6 ⟩
    ⟨ getheadervalue 9b ⟩
    ⟨ whichcolumn 10a ⟩
    ⟨ rdnums 7b ⟩
    ⟨ split 7a ⟩
    end module specfiles

```

◇

Fragment referenced in 11b, 12, 13, 14b, 15, 16, 31a, 32, 40a, 57b, 86b, 89, 91.

So in theory any programmer can come along in the future and if they add the line **use specfiles** to their programs then they'll have access to all those goodies. The file **specfiles.f90** should be compilable into a library (type thing).

```

"specfiles.f90" 11b ≡
    ⟨ specfiles 11a ⟩

```

◇

## 4.7 Test program

In order to check the routines in this section are working we have a simple program which gets a file name from the command line (or prompts) and dumps out header and column information on stdout.

```

"testspec.f90" 12≡
  < specfiles 11a>
    program testspec
    use specfiles
    integer(kind=4) :: i,j
    real(kind=8),allocatable :: data(:), prevdata(:)
    real :: time1, time2
    call getarg(1,filnam)
    call getfile
    call cpu_time(time1)
    call readheader ! Get's the header at the top of the file
2  write(*,'(a5,i5)') 'Scan ',iscan
    write(*,'(a)') 'Header words and values'
    do j=1,NLINES
      if(headerwords(j).gt.0)then
        do i=1,headerwords(j)
          write(*,'(a20,f16.8)') words(i,j),wordvalues(i,j)
        enddo
      endif ! headerwords(j).gt.0
    enddo
    write(*,'(a,i5)') 'Number of columns = ', ncolumns
    if(.not.allocated(data))allocate(data(ncolumns))
    if(.not.allocated(prevdata))allocate(prevdata(ncolumns))
    if(size(data).ne.ncolumns)then
      deallocate(data); deallocate(prevdata)
      allocate(data(ncolumns));allocate(prevdata(ncolumns))
    endif
    j=0
1  prevdata=data ! Start of loop through reading data
    call getdata(data,ncolumns)
    j=j+1
    if(ispecerr.eq.0)goto 1 ! loop to end of scan
    ispecerr=0 ! OK, we reached the end of the scan
    write(*,'(a)') 'Last data line in scan'
    do i=1,ncolumns
      write(*,'(a20,f16.8)') columnlabels(i),prevdata(i)
    enddo
    write(*,'(a,i10)') 'Number of data points = ',j
    call findscan(iscan+1) ! get the next scan (assumes sequential)
    if(ispecerr.eq.0)goto 2 ! process next scan
    deallocate(data); deallocate(prevdata)
    call cpu_time(time2)
    write(*,*)
    write(*,'(a,f8.4,a)') 'Time to read file = ',time2-time1,' seconds'
    write(*,'(a)') 'Thats all folks!'
    end program testspec

```

◇

This program could form the basis of a rebinning program - all it needs to do as extra work is to determine the user options for which scan to bin, which stepsize etc and then pass the data to the rebinning, summing and output routines. It could also be modified to write out summary information about the scan (total counts) or to pull out information about temperatures and so on.

During testing it was noted that there is a factor 5 in the speed that the program runs if the line of data is actually read in, or just checked to see if it the second character on the line is a digit. This suggests we might try to optimise the rdnums subroutine.

(FIXME) Some more work will be necessary (somewhere) to catch corrupted specfiles, and allow files to be read where scans are not in order.

## 4.8 Utility programs

Some useful programs for extracting data from SPEC files are given here. Plotit just pulls out a column of data from the SPEC file and creates an mtv format file for display using the plotmtv program.

```
"plotit.f90" 13≡
  < specfiles 11a>
    program plotit
    use specfiles
    integer(kind=4) :: j
    integer(kind=4) :: nscan, ncol
    character(LEN=WORDLENGTH) :: label
    real(kind=8),allocatable :: data(:), prevdata(:)
    real :: time1, time2
  < cmdline 14a>
    call cpu_time(time1)
    if(.not.allocated(data))allocate(data(ncolumns))
    if(.not.allocated(prevdata))allocate(prevdata(ncolumns))
    if(size(data).ne.ncolumns)then
      deallocate(data); deallocate(prevdata)
      allocate(data(ncolumns));allocate(prevdata(ncolumns))
    endif
    j=0
    write(*,'(a)')'$ DATA = CURVE2D'
    write(*,'(3a)')'% xlabel = "',columnlabels(1),'"'
    write(*,'(3a)')'% ylabel = "',columnlabels(ncol),'"'
    write(*,'(a)')'% linetype = 1 markertype = 2'
    write(*,'(a)')'% linecolor = 3 markercolor = 4'
1    prevdata=data          ! Start of loop through reading data
    call getdata(data,ncolumns)
    j=j+1
    if(ispecerr.eq.0)then
      write(*,*)data(1),data(ncol)
      goto 1
    endif
    deallocate(data); deallocate(prevdata)
    call cpu_time(time2)
    write(*,*)'$ END'
!    write(*,'(a,f8.4,a)')'Time to read file = ',time2-time1,' seconds'
!    write(*,'(a)')'Thats all folks!'
100  end program plotit                                ◇
```

The command line stuff is generally the same from some of the other programs below, so we separate it off here.

```

< cmdline 14a > ≡
    call getarg(1,filnam)
    call getfile
    call getarg(2,line)
    read(line,*,err=10,end=10)nscan
    call findscan(nscan)
    call getarg(3,label)
    read(label,*,err=11,end=11)ncol
    goto 20
! allow labels as well as numbers for column headings
11  ncol=whichcolumn(label(1:len_trim(label)))
    if(ncol.lt.1)then
        write(0,'(a,a)')'Couldn't find column labelled ',
&          label(1:len_trim(label))
        write(0,'(a,i5)')'Number of columns = ', ncolumns
        do j=1,ncolumns
            write(0,'(i5,1x,a)')j,columnlabels(j)
        enddo
        goto 10
    endif
    goto 20
10  write(0,*)'Probs with command line'
    goto 100
20  continue
    ◇

```

Fragment referenced in 13, 15, 16.

To get a list of the columns present in a particular scan of a particular file the columns program was written. In practice it's easier to use a **grep #L filename** from the unix prompt, the only advantage here is to get the columnnumbers. This is superseded by the modified plotit which takes a columnlabel in text as an argument.

```

"columns.f90" 14b≡
< specfiles 11a >
    program columns
    use specfiles
    integer(kind=4) :: i,j
    integer(kind=4) :: nscan, ncol
    real(kind=8),allocatable :: data(:), prevdata(:)
    real :: time1, time2
    call getarg(1,filnam)
    call getfile
    call getarg(2,line)
    read(line,*,err=10,end=10)nscan
    goto 20
10  write(*,*)'Probs with command line'
    write(*,*)'Usage: columns filename scan'
    goto 100
20  call findscan(nscan)
    if(ispecerr.eq.-2)goto 100
    write(*,'(a,i5)')'Scan number = ',iscan
    write(*,'(a,i5)')'Number of columns = ', ncolumns
    do j=1,ncolumns
        write(*,'(i5,1x,a)')j,columnlabels(j)
    enddo
100  end program columns
    ◇

```



For extracting data ready for peak fitting some error bars are required. The **c2xye** program does just that, and also tags on some header information for later plotting in plotmtv.

```
"c2xye.f90" 15≡
  < specfiles 11a>
    program c2xye
    use specfiles
    integer(kind=4) :: i,j
    integer(kind=4) :: nscan, ncol
    character(LEN=WORDLENGTH) :: label
    real(kind=8),allocatable :: data(:), prevdata(:)
    real :: time1, time2
  < cmdline 14a>
    if(.not.allocated(data))allocate(data(ncolumns))
    if(.not.allocated(prevdata))allocate(prevdata(ncolumns))
    write(*,'(a)')'$ DATA=CURVE2D'
    write(*,'(a)')'% linetype=0 markertype=2'
    write(*,'(3a)') '% xlabel="",columnlabels(1),'"'
    write(*,'(3a)') '% ylabel="",columnlabels(ncol),'"'
    write(*,'(a)') '% title="A fit"'
1    prevdata=data          ! Start of loop through reading data
    call getdata(data,ncolumns)
    if(ispecerr.eq.0)then
        write(*,*)data(1),data(ncol),sqrt(data(ncol)+1.) ! default weight
        goto 1
    endif
    deallocate(data); deallocate(prevdata)
    call cpu_time(time2)
100  end program c2xye ◇
```

Surprisingly there didn't seem to be anything for plotting a mesh scan in 3D, so we made something here. (mesh scan means scanning one motor through a range for a range of positions of a second motor, it's very useful for aligning things on a beamline).

```

"plotmesh.f90" 16≡
  < specfiles 11a>
    program plotmesh
    use specfiles
    integer(kind=4) :: i,j
    integer(kind=4) :: nscan, ncol
    character(LEN=WORDLENGTH) :: label
    real(kind=8),allocatable :: data(:), prevdata(:)
    real :: time1, time2, olddata2
    logical inc
  < cmdline 14a>
    call cpu_time(time1)
    if(.not.allocated(data))allocate(data(ncolumns))
    if(.not.allocated(prevdata))allocate(prevdata(ncolumns))
    if(size(data).ne.ncolumns)then
      deallocate(data); deallocate(prevdata)
      allocate(data(ncolumns));allocate(prevdata(ncolumns))
    endif
    j=0
    write(*,'(a)')'$ DATA = CURVE3D'
    write(*,'(3a)')'% xlabel = "',columnlabels(1),'"'
    write(*,'(3a)')'% ylabel = "',columnlabels(2),'"'
    write(*,'(3a)')'% zlabel = "',columnlabels(ncol),'"'
    inc=.true.
    call getdata(data,ncolumns)
    j=j+1
    prevdata=data
    write(*,*)data(1),data(2),data(ncol)
    call getdata(data,ncolumns)
    j=j+1
    write(*,*)data(1),data(2),data(ncol)
    olddata2=data(2)
1    prevdata=data                ! Start of loop through reading data
    call getdata(data,ncolumns)
    j=j+1
    if(ispecerr.eq.0)then
      if(data(2)-olddata2.gt.1e-6)write(*,*)
      write(*,*)data(1),data(2),data(ncol)
      olddata2=data(2)
      goto 1
    endif
    deallocate(data); deallocate(prevdata)
    call cpu_time(time2)
    write(*,*)'$ END'
100  end program plotmesh

```

◇

## 5 Rebinning the data

A collection of routines for taking raw counts and putting them into bins are collected together in a module here. A couple of test programs are included at the end of this section. This module is supposed to deal with all of the things we need that are to do with bins.

### 5.1 Useful variables

If we can get the data structures right then we might hope that the algorithms will become self evident (FIXME find insert original quote). There are a few bits of information we need before we are able to do

any rebinning at all, along with some space to actually store the data itself. In order to get the angular position of each detector we'll need to know it's offset:

```

< OFFSET 17 > ≡
    integer(kind=4), parameter :: NCHAN = 9  ! ID31 will have nine channels
    integer(kind=4) :: NCHANNEL = 9  ! ID31 will have nine channels
    real(kind=8), dimension(NCHAN) :: offset, mult, multerr
    logical :: tempres
    integer(kind=4) logexdet(NCHAN)
    integer(kind=4) :: iexrc ! region count
    integer(kind=4), allocatable, dimension(:) :: iexarray
    real(kind=8), allocatable, dimension(:, :) :: exarray
    data logexdet /9*0/ ! whole array initialised to false
    data iexrc /0/      ! no excluded regions by default

```

Fragment referenced in 28.

The **logexdet** stuff is for excluding channels. If it is set to zero then the channel is used as normal, if it is one the channel is completely thrown away, if it is two then the channel is excluded in particular regions. **iexrc** holds the number of lines read into

**NCHAN** is a parameter for the number of channels on the instrument and **offset** stores their  $2\theta$  offsets. The actual numbers for the offsets are initialised in the **initialiserebin** routine.

```

⟨ offsetdefaults 18 ⟩ ≡
! Germanium numbers
  offset(1) = 7.88643510d0
  offset(2) = 5.91013889d0
  offset(3) = 3.89405184d0
  offset(4) = 1.97036951d0
  offset(5) = 0.00000000d0
  offset(6) = -2.12832415d0
  offset(7) = -4.03585040d0
  offset(8) = -6.00076222d0
  offset(9) = -8.03007953d0

! Silicons latest numbers
  offset(1) = 8.0634515d0
  offset(2) = 5.8865159d0
  offset(3) = 3.9594961d0
  offset(4) = 2.0986688d0
  offset(5) = 0.0000000d0
  offset(6) = -1.9488783d0
  offset(7) = -3.9966086d0
  offset(8) = -6.0474594d0
  offset(9) = -8.0536348d0

! Some new numbers for silicon - there seems to be some drift?
  offset(1) = 8.05624406d0
  offset(2) = 5.88332391d0
  offset(3) = 3.95657893d0
  offset(4) = 2.09530059d0
  offset(5) = 0.00000000d0
  offset(6) = -1.94883199d0
  offset(7) = -3.99982480d0
  offset(8) = -6.04725698d0
  offset(9) = -8.05585677d0

! some more new numbers for silicon - there is drift for sure

  offset(1) = 8.02703050d0
  offset(2) = 5.88348041d0
  offset(3) = 3.95722668d0
  offset(4) = 2.09585757d0
  offset(5) = 0.00000000d0
  offset(6) = -1.94681946d0
  offset(7) = -3.99878112d0
  offset(8) = -6.04566287d0
  offset(9) = -8.05515342d0

```

◇

Fragment referenced in 23.

Those numbers came from a nac dataset collected with the 311 monochromator at about 30keV. It was one of the sharpest datasets JPW could find, although they can be updated with the **id31offsets** program and stored in the temp.res file if something better comes along.

The offsets and efficiencies were previously stored in a file called temp.res. This is an example temp.res file, the first column is the two theta offset, the second is the efficiency and the third is the esd on the efficiency. No idea what the fourth is all about (it is called Tau somewhere in the older programs, perhaps related to flat plate corrections?).

7.88736900,1.47159656,0.00318100,0.0

```

5.91065300,0.84581036,0.00268899,0.0
3.89374900,1.00674892,0.00283546,0.0
1.97086700,0.93762646,0.00276469,0.0
0.00000000,0.92496591,0.00275221,0.0
-2.12814400,0.86504408,0.00263751,0.0
-4.03525200,0.87438534,0.00263823,0.0
-6.00101500,0.88839695,0.00260034,0.0
-8.02961200,1.18542541,0.00276954,0.0

```

Note that these are the original "BM16" offsets - the ones hard wired into the program have now been updated. We read these files in section 5.3, where the initialisation is carried out. In practice it might be wise to make it fairly easy to modify the offsets or efficiencies independently of each other. They will be compiled into the program in the current version, although perhaps this could change and they could be placed in some file - temp.res would be a good place, with a flag on the efficiency esd's to indicate they are still to be determined (make them negative for example?).

The data format to hold a binned scan is going to need to be flexible enough to handle varying numbers of datapoints. Therefore we will want some dynamic memory allocation, with dimensions of twice the number of channels as we have so that there is space to carry along a separate monitor spectrum along for each channel. Also we make a note of the stepsize used, **step**, the range of  $2\theta$  values which the array is holding (**tthlow** and **tthhigh**) and the total number of points in **npts**. **aminstep** is a default minimum stepsize to prevent mistakes (like zero) causing the program unnecessary difficulties. The default values will hopefully suffice for most applications, but will be modifiable user options.

```

< SCAN 19 > ≡
    real(kind=8), allocatable :: ascan(:, :)
    real(kind=8) :: step, tthlow, tthhigh, aminstep
! if requested by the user
    real(kind=8) :: user_step = 0.003d0
    real(kind=8) :: user_tthlow = -30.0d0
    real(kind=8) :: user_tthhigh = 160.0d0
    integer(kind=4) :: npts
    data tthlow, tthhigh, step, aminstep /-30.0d0, 160.0d0, 0.003d0, &
    & 0.0002d0 /

```

Fragment referenced in 28.

That should be all we need, all the information for reading SPEC files is in the preceeding section, and will be used as needed here.

Memory usage by the program will be mainly in this **ascan** array. Given that the smallest stepsize allowed was going to be  $5 \times 10^{-5}$  (the precision of the diffractometer) and the maximum range for a pattern is  $360^\circ$ . For nine channels the array could have dimensions of  $18 \times \frac{360}{5 \times 10^{-5}}$ , which is  $1.296 \times 10^8$ . Since a double precision number occupies 64 bits, or 8 bytes of memory this array would use  $8 \times 1.296 \times 10^8 / (1024)^2 = 988\text{MB}$ . That is clearly too much for most machines, so we will bear this in mind and either reduce the two theta range to be used or increase the step size. Only about  $190^\circ$  should ever be needed, with a step of  $0.0002^\circ$ , 0.72 arcsec, we use about 130MB, which is not unreasonable. In practice, Bragg peaks are not actually that narrow. Given that the accuracy of the diffractometer is about 1 arcsec, this is around  $3 \times 10^{-4}^\circ$ , so binning into steps smaller than that would seem unwise. Since the peaks tend to get narrower as energy increases and the pattern is compressed, the effect of reduced angular range should reduce the memory requirement. It will always be possible to process large datasets in a series of angular ranges which can then be joined together, but it is not expected that this will be necessary for the current diffractometer.

The contents of the ascan array (and other large binning arrays) is described graphically in figure 1.

## 5.2 Some functions for binning

The bins are to have  $2\theta$  centers at tthlow, tthlow + step, tthlow + 2×step, etc. Functions for working out which bin a particular  $2\theta$  value falls into along with the limits and centre of the bin are given here. An




Figure 1: Illustration of the large arrays in use in the programs and their intended memory contents.

integer is which used for accessing the ascan array is used to label each bin, these will start at tthlow and end at tthhigh, returning -1 for a  $2\theta$  value which is out of range. Figure 2 will hopefully make it clearer what these functions are doing.

**ibin** returns the integer bin number for a particular bin.

```

< ibin 20a > ≡
integer(kind=4) function ibin(tth)
real(kind=8), intent(in) :: tth
real(kind=8)::x
if(tth.ge.tthlow .and. tth.le.tthhigh) then
  x=(tth-tthlow)/step+0.5d0
  ibin=int(x)
else ; ibin=-1 ; endif      ! Out of range for this ascan array
return ; end function ibin

```

Fragment referenced in 28.

The upper  $2\theta$  limit for a bin is to be given by a function tthhb:

```

< tthhb 20b > ≡
real(kind=8) function tthhb(n)
integer(kind=4), intent(in) :: n
integer(kind=4)nplusone
nplusone=n+1
tthhb=tthlb(nplusone)
return; end function tthhb

```

Fragment referenced in 28.




Figure 2: Illustration of the bin definitions. **ibin** calculates an integer to be assigned to any particular bin and the other functions find the limits and center of a bin labelled by **ibin**.

The lower  $2\theta$  limit for a bin is to be given by a function **tthlb**:

```

< tthlb 21a > ≡
  real(kind=8) function tthlb(n)
  integer(kind=4), intent(in) :: n
  tthlb = (real(n,8)-0.5d0)*step + tthlow
  return; end function tthlb

```

◇

Fragment referenced in 28.

The center of the bin is given by a function **bincen**:

```

< bincen 21b > ≡
  real(kind=8) function bincen(n)
  integer(kind=4), intent(in) :: n
  bincen = real(n,8)*step + tthlow
  return; end function bincen

```

◇

Fragment referenced in 28.

Checking the parameters before initialising the ascan array should avoid a lot of problems and also allow us to ensure that there is always a bin centered on  $2\theta = 0$ . We check **tthlow** and **tthhigh** and adjust them to be the edges of the bins we assume they were aiming at, given a bin centred at 0, and also we check the stepsize is physically reasonable. The diffractometer is sensitive to angular displacements of  $5 \times 10^{-5}$ , so any binsize smaller than that would be a bit silly.

```

⟨ checkrebinpars 22 ⟩ ≡
  subroutine checkrebinpars
    real(kind=8) :: x
    if ((units.ne.'T').and. (.not. wavelength_set) ) return
    if(step.lt.aminstep .and. (units.eq.'T'))then
      write(*,'(a,G12.4)')'Step is a bit small, resetting to ',aminstep
      step=aminstep
      user_step = step
    endif
    x=tthlow/step
    tthlow=real(int(x),8)*step
    x=(tthhigh-tthlow)/step
    npts=int(x)
    tthhigh=tthhb(npts)
    if (units .eq. 'T') then
      write(*,'(3(G12.5,a),G12.5)')tthlow,' < tth < ',tthhigh,      &
      &' step=',step,' npts=',npts
    endif
    if (units.eq.'Q') then
      write(*,'(3(G12.5,a),G12.5)')tthlow,' < 2pi/d < ',tthhigh,      &
      &' step=',step,' npts=',npts
    endif
    if (units.eq.'R') then
      write(*,'(3(G12.5,a),G12.5)')tthlow,' < Q^2 < ',tthhigh,      &
      &' step=',step,' npts=',npts
    endif
    return
  end subroutine checkrebinpars
  ◇

```

Fragment referenced in 28.

### 5.3 Initialisation

Some of the data will need to be initialised before any rebinning can be carried out. Traditionally a file called "temp.res" is used to hold the detector offsets. We will also need to decide upon the scan range and stepsize and allocate an array to hold the data.



```

< initialiserebin 23 > ≡
    subroutine initialiserebin
    integer(kind=4) :: ierr, i
    real(kind=8) :: junk
    !      real :: time1, time2
    ! constants to machine precision
    four_pi = 8.0d0*asin(1.0d0)
    pi_over_360 = asin(1.0d0)/180.0d0
    ! default values
    < offsetdefaults 18 >
    open(unit=16,file='temp.res',form='FORMATTED',
    & access='SEQUENTIAL', status='OLD', iostat=ierr)
    if(ierr.eq.0)then
        do i=1,nchannel
    ! Should clarify policy on whether to read/use these?
        read(16,*,err=10,end=12)offset(i),mult(i),multerr(i),junk
        enddo
    ! Report temp.res found and read in
        write(*,'(a)')'temp.res file found and read in'
        tempres=.true.
        goto 11
    endif
10  write(*,'(a)')'Not able to read temp.res file,'
    write(*,'(a)')'You need this file for the detector calibration'
    write(*,'(a)')'Please copy this file to the current directory'
    write(*,'(a)')' eg: " cp ~/temp.res ." '
    stop
12  write(*,'(a)')'Reached end of temp.res file early?'
11  close(16)
    ! two theta low, two theta high, step and that npts is correct
    if ((units.eq.'T') .or. wavelength_set) call unit_lims
    return ; end subroutine initialiserebin

    subroutine rebinallocate
    implicit none
    integer ierr
    if( allocated(ascan) ) deallocate(ascan) ! can reset
    allocate(ascan(nchannel*2,npts),stat=ierr)
    if(ierr.ne.0)stop 'Memory allocation error'
    !      call cpu_time(time1)
    ascan=0.0d0 ! Clear any junk
    !      call cpu_time(time2)
    !      write(*,*)'Time taken to zero array =',time2-time1,'/s'
    return ; end subroutine rebinallocate
    ◇

```

Fragment referenced in 28.

The commented out lines are used to find out if the program is thrashing the hard disk. If the time taken to zero the array is more than a small fraction of a second the program will run very slowly. This happens if it is running on a computer which has run out of RAM and is using the hard disk for virtual memory.

## 5.4 The binning algorithm

Finally we get to the guts of the thing, the actual rebinning algorithm. This will eventually take a line of entries from the SPEC file and assign them to bins. Some explanation of SPEC files is required here. Figure 3 shows a graphical representation of the data in the SPEC file. The starting  $2\theta$  position is given




Figure 3: Illustration of the SPEC file contents, the initial two theta value is in the file header, subsequent lines record the current two theta value and counts arriving since the last line, indicated by the arrows.




Figure 4: Illustration of the binning procedure. Counts are assigned to bins assuming they arrived uniformly during the sampling time.

in the file header and then datalines contain the current value of  $2\theta$  and the number of counts which have arrived at the detector since the last dataline.

For processing a single line of data we will need to know the value of  $2\theta$  from the previous line (or header) and the current  $2\theta$  value. The counts are then assumed to have been arriving uniformly between those two values and are apportioned to bins accordingly. The bins are designated by their borders and the convention will be taken that there is always a bin centred on zero. The algorithm will proceed as follows:

1. Obtain the lower and upper of the  $2\theta$  detector positions,  $2\theta_l$  and  $2\theta_u$
2. Identify the bin into which this  $2\theta_l$  falls and note the upper limit of the bin  $2\theta_b$
3. If the upper detector position is also within this bin put all the counts in this bin and exit
4. If the upper position is in the next bin apportion the counts to the two bins so the lower bins gets the fraction  $\frac{2\theta_b - 2\theta_l}{2\theta_u - 2\theta_l}$  and put the rest in the upper bin and exit.
5. Otherwise put the fraction  $\frac{2\theta_b - 2\theta_l}{2\theta_u - 2\theta_l}$  in the lower bin, the corresponding fraction in the topmost bin and split the rest of the counts amongst the rest of the bins. exit.

Figure 4 is intended to illuminate this process, in case the text wasn't clear.

The following subroutine takes an upper and lower  $2\theta$ , a number of counts, and a channel number and places these counts into the appropriate bin.

```

< bin 25 > ≡
      subroutine bin(tth1,tth2,cts,ichan)
      ! tthh & tthl are the high & low two thetas from the SPEC file
      ! cts is the number of counts and ichan is the column to use in ascan
      real(kind=8), intent(in) :: tth1, tth2, cts
      integer(kind=4), intent(in) :: ichan
      integer(kind=4) :: ibl, ibh, j
      real(kind=8) :: frac, tthh, tthl
      !   tthl=min(tth1,tth2)                                ! Ensures ascending data
      !   tthh=max(tth1,tth2)    ! strange bug on Fe304 exp??
      if(tth1 .ge. tth2)then
        tthl=tth2; tthh=tth1
      else
        tthl=tth1; tthh=tth2
      endif
      ibl = ibin(tthl)                                ! Bin of lower point
      ibh = ibin(tthh)                                ! Bin of upper point
      if ((ibl.le.0).or.(ibh.lt.ibl).or.(ibh.gt.npts))return ! error
      if ( ibh .eq. ibl) then                          ! All in one bin
        ascan(ichan,ibl) = ascan(ichan,ibl) + cts
        return
      elseif ( ibh .eq. ibl+1) then                    ! Spread over two bins
        frac=(tthh-ibl)/(tthh-tthl)                  ! Fraction in bin 1
        ascan(ichan,ibl)=ascan(ichan,ibl) + cts*frac
        ascan(ichan,ibh)=ascan(ichan,ibh) + cts*(1.0d0-frac)
        return
      else
        frac=(tthh-ibl)/(tthh-tthl)                  ! Fraction in bin 1
        ascan(ichan,ibl)=ascan(ichan,ibl) + cts*frac
        frac=(tthh-tthl)/(tthh-tthl)                  ! Fraction in last bin
        ascan(ichan,ibh)=ascan(ichan,ibh) + cts*frac
        frac= step/(tthh-tthl)                        ! Fraction in middle bins
        do j = ibl+1,ibh-1,1
          ascan(ichan,j)=ascan(ichan,j) + cts*frac
        enddo
        return
      endif
    end subroutine bin

```

◇

Fragment referenced in 28.

Finally we create a routine that takes a line of data from a SPEC files and uses the rebinning code above to process the entire line of entries. The detector offset must be added to the  $2\theta$  value for each detector before it can be binned. The specification of the contents of the ascan array is decided here, column 1 will correspond to detector 1, column 2 is a monitor column for detector 1 and so on. Users of the routine must supply a high and low  $2\theta$  value along with a set of counts and a monitor column.

```

⟨ processline 26 ⟩ ≡
  subroutine processline(tth1,tth2,cts,n,ctsmon)
  real(kind=8), intent(in) :: tth1, tth2, ctsmon
  integer(kind=4), intent(in) :: n
  real(kind=8), intent(in), dimension(n) :: cts
  integer(kind=4) :: i, ich
  real(kind=8) :: tthh, tthl
  do i = 1,n
    if(logexdet(i).eq.1)cycle ! skip if excluded completely
    if(userandomstart.and.rstchan(i).and.(tth2.lt.randomval))cycle
    if(userandomend.and.renchan(i).and.(tth2.gt.randomvalend))cycle
    tthl = tth1 - offset(i)
    tthh = tth2 - offset(i)
    if(logexdet(i).eq.2)then
! Is this an excluded region for this channel
      if(led(i,tthl,tthh))cycle
    endif
    tthl = convert_unit_function( tthl )
    tthh = convert_unit_function( tthh )
    ich = 2*i-1
    call bin(tthl,tthh,cts(i),ich)
    ich = 2*i
    call bin(tthl,tthh,ctsmon,ich)
    sumtotal(i)=sumtotal(i)+cts(i)
  enddo
  sumtotalmon=sumtotalmon+ctsmon
  return
end subroutine processline

```

Fragment referenced in 28.

Checking if a particular channel is excluded in a particular two theta interval needs the function **led**. It just runs through a stored version of the excluded region file to return true or false.

⟨ *led* 27 ⟩ ≡

```

logical function led(ic,xl,xh)
use specfiles
integer(kind=4),intent(in) :: ic
real(kind=8),intent(in) :: xl, xh
integer(kind=4)i
if(.not.allocated(iexarray).or. .not. allocated(exarray)) then
! oh dear, function should never have been called
  write(*,'(a)')'Bug in program, bailing out, sorry!'
  write(*,'(a)')'Please mail a bug report to wright@esrf.fr'
  stop
endif
led=.false.
do i=1,iexrc
  if((iexarray(i)+1).eq.ic)then ! found the right line
    if(int(exarray(i,3)).lt.1)then
      if(xl.gt.exarray(i,1) .and.xl.lt.exarray(i,2))led=.true.
      if(xh.gt.exarray(i,1) .and.xh.lt.exarray(i,2))led=.true.
    else
      if(xl.gt.exarray(i,1) .and.xl.lt.exarray(i,2) .and. &
& iscan.eq.int(exarray(i,3)))led=.true.
      if(xh.gt.exarray(i,1) .and.xh.lt.exarray(i,2) .and. &
& iscan.eq.int(exarray(i,3)))led=.true.
    endif
  endif
endif
enddo
return
end function led

```

◇

Fragment referenced in 28.

## 5.5 A rebinning module

The data and functions created so far should be all that we need to carry out the rebinning process. They are collected together for in module so they can be used elsewhere.

```

< rebin 28 > ≡
    module rebin
    < OFFSET 17 >
    < SCAN 19 >
        real(kind=8) :: sumtotal(nchan), sumtotalmon, minmon=1.0d0
        real(kind=8) :: winhigh,winlow,winhighread,winlowread
        integer(kind=4) :: wincol
        real(kind=8) :: minrenormsig=5.0
        real(kind=8) :: randomstart = 0, randomval=0
        real(kind=8) :: randomend = 0, randomvalend=0
        character(len=90) :: wincnt
        logical :: winlog=.false.
        logical :: userandomstart = .false.
        logical :: rstchan( nchan ) = .false.
        logical :: userandomend = .false.
        logical :: renchan( nchan ) = .false.
        logical :: wavelength_set = .false.
        real(kind=8) :: wavelength=0.0d0
        ! T = Two theta
        ! Q = 4 * pi * sin( two_theta * pi / 360.0 ) / wavelength
        ! R = [q squared] = Q * Q
        character(len=1) :: units = 'T'
        real(kind=8) :: four_pi, pi_over_360
    contains
    < ibin 20a >
    < tthhb 20b >
    < tthlb 21a >
    < bincen 21b >
    < checkrebinpars 22 >
    < initialiserebin 23 >
    < bin 25 >
    < processline 26 >
    < led 27 >
    < window 37 >
    < convertunitfunction 30 >
    < random 136a >
    end module rebin

```

◇

Fragment referenced in 31a, 32, 40a, 57b, 86b, 89, 91.

Now we need some information from the user on the command line. Namely, the unit to use (two theta remains as a default) and the wavelength as an option

$\langle \text{unitoptions 29} \rangle \equiv$

```

! These will end up in the useroptions module
!       case('wvl')  call setwavelength(string)
!
      subroutine setwavelength(s)
      use rebin
      character(len=*), intent(in) :: s
      if(s(1:5).eq.'wvln=')read(s(6:len_trim(s)),*,err=1)wavelength
      wavelength_set = .true.
      write(*,*)'Using wavelength of ',wavelength
      return
1     write(*,*)'Error reading wavelength',s
      stop
      end subroutine setwavelength
!
!       case('uni')  call setunits(string)
      subroutine setunits(s)
      use rebin
      character(len=*), intent(in) :: s
!       12345678
      if( s(1:8) .eq. 'units=Q2' ) then
          units = 'R'
          write(*,*)'Binning into Q^2 = 16.pi^2.sin^2(theta)/wavelength^2'
          return
      endif
!       1234567
      if( s(1:7) .eq. 'units=Q' ) then
          units = 'Q'
          write(*,*)'Binning into Q = 4.pi.sin(theta)/wavelength'
          return
      endif

      write(*,*)'Did not understand your option',s
      write(*,*)'Use nothing for twotheta, units=Q or units=Q2'
      stop
      end subroutine setunits

```

◇

Fragment referenced in 87a.

This module now contains everything we will need for rebinning our data.

## 5.6 Converting to another unit

March 2009, some people want to bin into constant steps in  $Q$ ,

$$Q = \frac{4\pi \sin \theta}{\lambda}$$

We get the wavelength,  $\lambda$  either from the command line, or from the header of the scans  $Q(4)$  fortran or  $Q[3]$  in C.

For completeness we'll do  $Q$  squared at the same time. Might be interesting for plotting data on a log scale (thermal factors go as  $q$  squared).

D-spacing offers an exciting possibility for dividing by zero, so it is skipped for now.

$\langle \text{convertunitfunction } 30 \rangle \equiv$

```

real(kind=8) function convert_unit_function( x )
real(kind=8), intent(in) :: x
real(kind=8) :: qt
select case((units))
  case('T') ! Two theta, do nothing
    convert_unit_function = x
    return
  case('Q') ! Q - constants from initialise_rebin
    convert_unit_function = four_pi*sin(x*pi_over_360)/wavelength
    return
  case('R') ! Q^2
    qt = four_pi * sin( x * pi_over_360 )/wavelength
! use a signed quantity
!   write(*,*)qt,qt*qt,SIGN( qt * qt , qt ),x
    convert_unit_function = sign( qt * qt , qt )
    return
  case default
    convert_unit_function = x
    return
end select
return
end function convert_unit_function

!
!
! To be called once per scan - picks up wavelength from command line
subroutine convert_unit_setupQ( Q )
! from rebin module :::: logical :: wavelength_set = .false.
real(kind=8), dimension(60), intent(in) :: Q
if ( (units .eq. 'T') .or. wavelength_set) then
  return
endif
if( Q(4) .gt. 1.0E-15 ) then
  wavelength = Q(4)
! Assume always the same for all scans:
  write(*,*)
  write(*,*)'Got wavelength',wavelength,'from spec file Q[3]'
  wavelength_set = .true.
endif
! Check for errors
! let them have Angstroms as meters (1e-11)
if ((wavelength .lt. 1.0E-15) .and. (units .ne. 'T') ) then
  write(*,*)'Your wavelength is a bit small ',wavelength
  write(*,*)'Giving up, try putting wvln=1.234 on command line'
  stop
endif
!   write(*,*)'calling setupQ from convert_unit_setupQ'
call unit_lims()
return
end subroutine convert_unit_setupQ

subroutine unit_lims()
real(kind=8) :: stmp
! apply the limits according to if the user changed them
tthlow = convert_unit_function( user_tthlow )
tthhigh = convert_unit_function( user_tthhigh )
! set the step size to be right at 30 degrees twotheta
step = convert_unit_function( 30.0d0+user_step )
step = step - convert_unit_function( 30.0d0 )
! round it to be some printable representation
write(*,*)
! round to something nicely printable
if(units.ne.'T')then

```



## 5.7 Test programs

A quick test of the functions for defining the bins is a program to print out a series of  $2\theta$  values at the start and ends of the scan and around zero. Then also to print a series of  $2\theta$  values and the bins which they are assigned to. Testbins1 checks the initialisation routines and fires off some numbers to see that they end up in the right bins.

```
"testbins1.f90" 31a≡
  < specfiles 11a>
  < rebin 28>
  < report 31b>
      program testbins1
      use rebin ! pull in the module defined above
      integer(kind=4) :: i, j, n
      real(kind=8) :: tth, steplocal
      character *80 string
      step=0.1d0 ; tthlow=-10.023487d0 ; tthhigh=60.0129384d0
      call getarg(1,string)
      read(string,*,err=1, end=1)step
      call getarg(2,string)
      read(string,*,err=1, end=1)tthlow
      call getarg(3,string)
      read(string,*,err=1, end=1)tthhigh
1      write(*,'(a)')'Calling initialise rebin'
      call initialiserebin
      call report
      write(*,'(a)')'Calling initialise rebin'
      call initialiserebin
      call report
! Now for some numbers to bin
      write(*,'(a)')'Two Theta, IB, LOW, HIGH, CENTRE'
      steplocal=0.0242349876
      n=nint((tthhigh-tthlow)/steplocal)-1 ! random hard to bin numbers?
      do i=1,n
          tth=(tthhigh-tthlow)*float(i)/float(n) + tthlow
          j=ibin(tth)
          write(*,'(f15.9,I9,3f12.4)')tth,j,tthlb(j),tthhb(j),bincen(j)
      enddo
      end program testbins1
```

◇

This brief reporting routine outputs the status of some variables after initialising things.

```
< report 31b> ≡
      subroutine report
      use rebin
      integer(kind=4) :: i, j, k
1000 format('low ',f12.5,' high ',f8.5,' step ',f6.4,' npts ',i5)
1001 format('bin ',i5,' low ',f12.5,' cen ',f12.5,' high ',f12.5)
      write(*,1000)tthlow,tthhigh,step,npts
      do i=1,4
          j=npts+2-i
          k=i-2
          write(*,1001)j,tthlb(j),bincen(j),tthhb(j)
          write(*,1001)k,tthlb(k),bincen(k),tthhb(k)
      enddo
      i=ibin(0.0d0)
      write(*,1001)i,tthlb(i),bincen(i),tthhb(i)
      end subroutine report
```

◇

Fragment referenced in 31a.

The program was tested with Compaq Visual Fortran (6.6A) with options options for full warnings and standards checking. No problems were detected and the output suggests the functions for assigning  $2\theta$  values to bins are correct, so any remaining bugs are fairly subtle. Part of the output is as follows:

```
Calling initialise rebin
temp.res file found and read in
-10.000    < tth <    60.050    step= 0.10000    npts=          700
low   -10.00000 high 60.05000 step 0.1000 npts 700
bin   701 low    60.05000 cen    60.10000 high    60.15000
bin   -1 low   -10.15000 cen   -10.10000 high   -10.05000
bin   700 low    59.95000 cen    60.00000 high    60.05000
bin    0 low   -10.05000 cen   -10.00000 high   -9.95000
bin   699 low    59.85000 cen    59.90000 high    59.95000
bin    1 low   -9.95000 cen   -9.90000 high   -9.85000
bin   698 low    59.75000 cen    59.80000 high    59.85000
bin    2 low   -9.85000 cen   -9.80000 high   -9.75000
bin   100 low   -0.05000 cen    0.00000 high    0.05000
...
    tth          bin      low limit  high limit  centre
    ...
    59.443821391      694      59.3500    59.4500    59.4000
    59.468068536      695      59.4500    59.5500    59.5000
    ...
```

We can check the correctness of the rebinning algorithm by generating  $2\theta$  values and counts in the range of a hypothetical scan and then asking for them to be rebinned. A gross check for consistency is that the total number of counts in the rebinned scan must be the same as the total number of counts supplied for binning.

```
"testbins2.f90" 32≡
  < specfiles 11a >
  < rebin 28 >
    program testbins2
    use rebin ! pull in the module defined above
    integer(kind=4) :: i
    real(kind=8) :: tth1, tth2, x, counts, s, c, t, y
    ! initialiserebin sorts out tthlow and tthhigh to be bin edges
    step=0.001d0 ; tthlow=-1.000561344d0 ; tthhigh=11.023479875d0
    npts=nint((tthhigh-tthlow)/step)+1
    call initialiserebin
    s=0.0d0 ; c=0.0d0 ; t=0.0d0 ; y=0.0d0
    do i=1,1000000
      call random_number(x)          ! Returns x in range 0. <= x < 1.
      tth1 = x*10.0d0                ! From 0 to 10 degrees
      tth2 = tth1 + (x-0.5)/100.0d0  ! Random offset
      counts = ( x + 0.1 ) * 1000.0d0 ! Random counts in range
      y=counts-c; t=s+y; c=(t-s)-y ; s=t ! Kahan's summation formula
      call bin(tth1,tth2,counts,1)    ! put everything in channel 1
    enddo
    x=ascan(1,1); c=0.0d0 ; t=0.0d0 ; y=0.0d0
    do i=2,npts
      y=ascan(1,i)-c;t=x+y;c=(t-x)-y; x=t ! Kahan's summation formula
    enddo
    write(*,*)'Binned cts=',x,' from ',s,' with f90 sum of bins=',    &
    & sum(ascan(1,:))
    end program testbins2
```

◇

The program proceeds with the following output, which suggests our algorithms are working. Note that during testing it was noted that we drop things in the sixth decimal place using (kind=4) real numbers (single precision), which motivated a change to making everything double precision. For the purposes of powder diffraction experiments the decimal places below ought to be enough, although we should be good to 15 with kind=8, suggesting something is a bit weird here <sup>1</sup>. Using Kahan's summation formula didn't help the precision much, but nevertheless, losing 0.05 of a count in a total of 600 million is probably accurate enough for our purposes.

```
Binned cts= 599869324.123505      from 599869324.172561
with f90 sum of bins= 5.9986931E+08
```

Some further investigation showed that replacing the **tthhb(n)** function (which originally calculated its own number) with a call to **tthln(n+1)** preserves the extra fraction of a count. The difference (I think) is due to not having exactly the same number for the top and bottom edges of the bin. Exactly why that difference arose remains unclear, but I've a feeling there was a rounding question somewhere.

```
temp.res file found and read in
-1.0000      < tth < 11.024      step= 0.10000E-02 npts= 12023
Binned cts= 599869324.172561      from 599869324.172561
with f90 sum of bins= 599869324.172561
```

Finally, a smallish program which is intended to take a single scan from within a SPEC file in BM16/ID31 variable step size and convert it to a commensurate stepsize, incorporating the detector offsets as it goes along.

Processing a series of lines of data from a SPEC file by repeatedly calling the **processline** routine is accomplished by a routine here called **processscan**. Given the number of a particular scan it will attempt to sum the counts from that scan into the array **ascan**.

A routine for processing an entire scan is included here. It is not part of any module, to prevent there being a dependency between the rebinning stuff and the specfiles stuff. This is the glue that holds them together. (so if we stopped using these beastly specfiles, the rebinning stuff will still be useful). It is tied to only processing turboscans for now, as processing the other types of scan with these methods would be a bit silly.

```
< processscan 33 > ≡
  < mma 36 >
  < logmotors 38 >
  < pointfilter 41 >
  subroutine processscan(n)
  use specfiles
  use rebin
  use outputfiles
  use pointfilter
  integer(kind=4),intent(in) :: n
  real(kind=8) :: tth1, tth2, tthf, x
  real(kind=8), allocatable :: chantot(:), a(:) ! ncolumns
  integer(kind=4) :: mon, ma0, ma8, itth, i
  integer(kind=4) :: igoodpoints, ibadpoints
  character(len=80) :: s
  call findscan(n) ! Positions the file on that scan
  call logmotors ! dumps all the starting motor positions to the log file
  if(ispecerr.ne.0) return
  write(*,'(i5,$)')iscan
  if(scantype.ne.'turboscan' .and. scantype.ne.'hookscan' &
  & .and. scantype.ne.'cscan' .and. scantype.ne.'zapline')then
    write(*,'(a)') 'is not a turboscan or a hookscan, ignoring it'
    ispecerr=-1
  return
  endif
```

---

<sup>1</sup>Have a look at "What every scientist should know about floating point arithmetic" from <http://docs.sun.com>.

```

tth1=getheadervalue('2_theta')
tthf=tth1
mon=-1;ma0=-1;itth=-1;itth=-1      ! FIXME
mon=whichcolumn(MONITORCOL)
if(mon.lt.1)goto 2
ma0=whichcolumn(FIRSTDET) ! FIXME deal with absent dets
if(ma0.lt.1)goto 2
ma8=whichcolumn(LASTDET) ! FIXME deal with absent dets.
itth=whichcolumn(TWOTTH)
!   write(*,*)"columns",mon,ma0,ma8,itth
if(mon.lt.1)goto 2
! Allocate a and chantot if necessary (cannot do this till after findscan)
if( allocated(a) ) deallocate(a) ! resets
if( allocated(chantot) ) deallocate(chantot)
allocate(a(ncolumns))
allocate(chantot(ncolumns))
chantot=0.0d0 ; igoodpoints=0 ; ibadpoints=0; a=0.0d0
! Set up any unit conversions which are requested
!   Send the current Q line from the specfile module to the rebin module
call convert_unit_setupQ( Q )
! Skip the first line of data
call getdata(a,ncolumns)
if(wincnt(1:6) .eq. 'Epoch') then
  write(*,*)'Scan Epoch starts at',a(wincol)
  winlow = winlowread+a(wincol)
  winhigh = winhighread+a(wincol)
  write(*,*)'Adjusting your Epoch window to',winlow,winhigh
endif
tth1=a(itth)
if(filterlogical)then ! initialise
  call filterinit()
  do i=1,3
    call getdata(a,ncolumns)
    if(ispecerr.ne.0)goto 2
    call pf(a,itth,ma0,ma8,mon)
    tth1=a(itth)
  enddo
endif
if( userandomstart ) then
  randomval = randomstart * rlcg() + tth1
  write(*,*)'rst using',randomval,'from',tth1
endif
if( userandomend ) then
  randomvalend = scanend - randomend * rlcg()
  write(*,*)'ren using',randomvalend,'from',scanend
endif
! Main loop
1  call getdata(a,ncolumns)
  if(filterlogical) call pf(a,itth,ma0,ma8,mon)
  do i=0,nchannel-1
    if(a(ma0+i).lt.0)then
      a(mon)=0.0d0 ! set monitor to zero if negative
! why not set all channels to zero for negative counts??
      write(*,*)'negative counts found for line'
      write(*,*)line(1:len_trim(line))
    endif
  enddo
  if(ispecerr.eq.0)then
    tth2=a(itth)
    if(a(mon).gt.minmon .and. window(a,ncolumns)) then

```

```

        igoodpoints=igoodpoints+1
        chantot=chantot+a ! sum on totals
        call processline(tth1,tth2,a(ma0:ma8),nchannel,a(mon))
        call mma(a,ncolumns,1)
    else
        ibadpoints=ibadpoints+1
    endif
    tth1=tth2
    goto 1          ! loops here
endif
2    continue
< reportsums 39 >
return
end subroutine processscan

```

◇

Fragment referenced in 40a, 57b, 89, 91.

Temperature processing has currently fallen into the domain of processing a scan. We are assuming that the main interest is to see the maximum, minimum and mean temperatures during a scan and so we only provide that information, to the **.log** file which is produced during the binning run. Since the temperature controllers have been renamed and are more widely available, we provide the minimum, maximum and average values of each column in the specfile. Some of them will be temperature.

```

⟨ mma 36 ⟩ ≡
  subroutine mma(arr,n,mode)
  use specfiles
  use outputfiles
  implicit none
  integer(kind=4),intent(in) :: mode,n
  real(kind=8),dimension(n),intent(in) :: arr
  real(kind=8),allocatable, dimension(:,:),save :: mmasum
  real(kind=8) :: rn
  integer(kind=4) :: i
  integer(kind=4),save :: nt, nwas
  character(len=80) :: s
  if(mode.eq.2 .and. nt.ne.0)then ! write to log file
    rn=real(nt,8)
    do i=1,nwas
      mmasum(i,1)=mmasum(i,1)/rn
    ! The absolute value is to avoid problems with rounding errors when there
    ! is zero variance in the data. eg: blower records a constant when the
    ! serial line is unplugged so we can end up trying to take a square root
    ! of a very tiny negative number, instead of a tiny positive number.
      mmasum(i,2)=sqrt(abs(rn*(mmasum(i,2)/rn-mmasum(i,1)**2)/(rn-1.0d0)))
      write(s,1000)iscan,columnlabels(i),mmasum(i,1),mmasum(i,2),nt
1000   format('Scan ',i5,' Ctr ',a10,' Avg =',F10.4,' +/- ',      &
    &   F10.4,' npts = ',i10)
      call wlogfile(s)
      write(s,1001)iscan,columnlabels(i),mmasum(i,4),mmasum(i,3)
1001   format('Scan ',i5,' Ctr ',a10,' T max = ',F10.4,' T min = ',F10.4)
      call wlogfile(s)
      enddo
      nt=0
      deallocate(mmasum)
    endif
    if(mode.eq.1)then
      if(.not. allocated(mmasum))then
        allocate(mmasum(ncolumns,4))
        mmasum=0.0d0
        nt=0
        nwas=ncolumns
      endif
      mmasum(:,1)=mmasum(:,1)+arr ! for avg
      mmasum(:,2)=mmasum(:,2)+arr*arr ! for std dev
      if(nt.eq.0)mmasum(:,3)=arr
      if(nt.eq.0)mmasum(:,4)=arr
      nt=nt+1 ! no of points
      if(nwas.eq.ncolumns)then
        do i=1,ncolumns
          if(mmasum(i,3).gt.arr(i))mmasum(i,3)=arr(i) ! min vals
          if(mmasum(i,4).lt.arr(i))mmasum(i,4)=arr(i) ! max vals
        enddo
      endif
    endif
    return; end subroutine mma

```

Fragment referenced in 33.

The formula used for the error on the temperature is:

$$\sigma = \sqrt{\sigma_{n-1}^2} = \sqrt{\frac{n}{n-1} \left[ \frac{\sum x^2}{n} - \langle x \rangle^2 \right]}$$

We are choosing to ignore that in practice the temperatures recorded in the SPEC file are not actually independent observations of the temperature. For most sample environments the temperature is measured much less frequently than the sampling time for the detectors and diffractometer position.

In the unlikely event(?) of temperature excursions during a scan a window function is created, which checks to see if the various counters are within an user defined window.

```

⟨ window 37 ⟩ ≡
    logical function window(ar,n)
    use specfiles !
    integer(kind=4),intent(in) :: n
    real(kind=8),dimension(n),intent(in) :: ar
    if(winlog)then
        wincol=whichcolumn(wincnt(1:len_trim(wincnt)))
        if(wincol.lt.0)then
            write(*,*)'PROBLEMS WITH YOUR WINDOW - cannot find column', &
& wincnt(1:len_trim(wincnt))
            write(*,*)'No further windowing will be attempted'
            winlog=.false.
            window=.true.
        else
            if( (ar(wincol).gt.winlow) .and. (ar(wincol).lt.winhigh)) then
                window=.true.
            else
                window=.false.
            endif
        endif
    else
        window=.true.
    endif
end function window

```

◇

Fragment referenced in 28.

Logging the motor positions (into a greppable format) might be handy enough for some users to make it a routine thing to do. Here is a routine which does just that.

```

< logmotors 38 > ≡
      subroutine logmotors
      use specfiles
      use outputfiles
      integer(kind=4) :: i, j
      character(len=80):: s
      do j=1,NLINES
        if(headerwords(j).gt.0)then
          do i=1,headerwords(j)
            write(s,1000)iscan,words(i,j),wordvalues(i,j)
1000      format('Scan ',i5,1X,a20,1X,f16.8)
            call wlogfile(s)
          enddo
        endif ! headerwords(j).gt.0)
      enddo
      write(s,1001)iscan, ncolumns
1001 format('Scan ',i5,' number of columns = ',i5)
      call wlogfile(s)
      do j=1,ncolumns
        write(s,1002)iscan,j,columnlabels(j)
1002 format('Scan ',i5,' column ',i3,1X,a30)
        call wlogfile(s)
      enddo
      return; end subroutine logmotors

```

Fragment referenced in 33.

The **processscan** subroutine will write out some summary information about each of the scans to the logfile. For clarity the code to do this is here - but is actually just a part of the **processscan** routine. This is included to replace the information which tended to flash past on the screen when the old programs were run.



```

⟨ reportsums 39 ⟩ ≡
    write(s,1000)iscan,chantot(mon)
1000 format('Scan ',i5, ' Total monitor ',F25.0)
    call wlogfile(s)
    write(s,1001)iscan,igoodpoints,igoodpoints+ibadpoints
1001 format('Scan ',i5, ' used ',i20,' of ',i20,' points' )
    call wlogfile(s)
    do i=ma0,ncolumns
        write(s,1002)iscan,columnlabels(i),chantot(i)
1002 format('Scan ',i5,1x,a20,' total ',F25.0)
        call wlogfile(s)
    enddo
    write(s,1003)iscan,chantot( whichcolumn('Seconds') )
1003 format('Scan ',i5,1x,'Total time          ',F25.0,' seconds')
    call wlogfile(s)
    ! Average step size = total range / total points
    x=abs((tthf-tth1)/real((igoodpoints+ibadpoints),8))
    write(s,1004)iscan,x
1004 format('Scan ',i5,1x,'Average stepsize',E25.15)
    call wlogfile(s)
    ! no of bins = step/(avg step)
    write(s,1005)iscan,step/x
1005 format('Scan ',i5,1x,'Average points per bin',E25.10)
    call wlogfile(s)
    write(s,1006)iscan,tthf,tth1
1006 format('Scan ',i5,1x,'low tth',E25.10,1x,'last tth',E25.10)
    call wlogfile(s)
    write(s,1007)iscan,igoodpoints,ibadpoints
1007 format('Scan ',i5,1x,'gd',i15,1x,'bad',i15)
    call wlogfile(s)
    call mma(a,ncolumns,2)

```

Fragment referenced in 33.

Now for the program:

```

"bindump.f90" 40a≡
  < specfiles 11a >
  < rebin 28 >
  < summation 56 >
  < outputfiles 72 >
  < processscan 33 >
  < useroptions 87a >
  < tidyup 98 >
  < helpmsg 87b >
  < bindumpmsg 40b >
    program bindump
    use specfiles
    use rebin
    use outputfiles ! pulls in rebin and summation
    character(len=256)::string
    integer(kind=4)::n
    call getarg(1,filnam)
! Get stepsize to use (defaults if not supplied)
    step=0.001 ; tthlow=-2.0 ; tthhigh=2.0
    call getarg(2,string)
    if(string(1:1).ne.' ')read(string,*)step
! Scan to treat (defaults if not supplied)
    iscan=1
    call getarg(3,string)
    if(string(1:1).ne.' ')read(string,*)n
    call getarg(4,string)
    if(string(1:1).ne.' ')read(string,*)tthlow
    call getarg(5,string)
    if(string(1:1).ne.' ')read(string,*)tthhigh
!
    call getfile
    call initialiserebin
    call processscan(n)
    call dumpscan(n)
    call tidyup
    end program bindump
    ◇

< bindumpmsg 40b > ≡
    1000 format('Usage would be bindump filename stepsize scan'/'
    & 'eg: bindump some_data.dat 0.003 23 ')
    return; end subroutine helpmsg
    ◇

```

Fragment referenced in 40a.

The program was tested with a couple of SPEC files and appears to function correctly. This is the first program which might actually be useful, for pulling out detector counts without any normalisation or summation. The routine **tidyup** is not defined until section 10, it will just run through and free any allocated memory, close any opened files etc.

## 5.8 Median filtering

When electronic noise appears in the detector channels a nasty workaround for the problem might be to take a 3 point median. The idea is that we read in lines of the data file 3 at a time and use the middle intensity value of the 3 for the middle data point, so if there is a big spike we will take one of the points on either side. The first and last points get their monitor set to zero to eliminate them.

```

⟨pointfilter 41⟩ ≡
  module pointfilter
    logical :: filterlogical=.false.
    real(kind=8), dimension(10,3) :: lt
    real(kind=8) :: tth, tthnew
    integer(kind=4) :: nc=0, np=1
    contains
    subroutine filterinit()
      lt(:,1)=0.0
      lt(:,2)=0.0
      lt(:,3)=0.0
      tth=-9999.0
      np=1
    return; end subroutine filterinit
    subroutine pf(a,itth,MA0,MA8,M)
      integer(kind=4), intent(in) :: itth, MA0, MA8, M
      real(kind=8), intent(inout), dimension(:) :: a
      integer :: i, j
      real(kind=8) :: tthnew
      tthnew=a(itth)
      a(itth)=tth
      tth=tthnew ! swap for last point
      lt(1:9,np)=a(MA0:MA8)
      lt(10,np)=a(M)
      do i = MA0, MA8
        j = i-MA0+1
        a(i)=lt(j,middle(lt(j,:)))
      enddo
      a(M)=lt(10,middle(lt(10,:)))
      np=np+1
      if(np.eq.4)np=1 ! Rolling buffer
    end subroutine pf
    integer function middle(x)
      real(kind=8), dimension(3), intent(in) :: x
      integer,dimension(1) :: i,j
      i=maxloc(x)
      j=minloc(x)
      middle=6-i(1)-j(1) ! 1+2+3 == 6
      if(middle.gt.3)middle=3
    end function middle
  end module pointfilter

```

◇

Fragment referenced in 33, 86b.

## 6 Summation and normalisation

Combining the data from multiple scans and or files is straightforward, provided they can be binned onto the same  $2\theta$  scale then the counts from the detectors and monitor are just added together. Merging scans from multiple detectors requires a detector efficiency correction. This section provides another module, which builds on the code in the binning module to carry out the merging of scans and detectors. It is intended that any multiple datasets will be summed together, keeping the detectors separate from each other, before counts from the various detectors are combined.

### 6.1 Combining data from multiple detectors

At the simplest level we just add the data together taking the detector efficiency into account. How precisely should the efficiency be taken into account? Imagine the situation where one of the efficiencies

is very low - the detector has been ignoring many of the photons it was supposed to be detecting <sup>2</sup>. We can think about this detector as effectively having been part of an experiment which had a lower incident flux, and so we would multiply the monitor by the efficiency and then sum the numbers in. In doing so we should be aware that error on the monitor spectrum will no longer be  $\sqrt{\text{counts}}$  and propagate the correct error. Provided the monitor counts are always far in excess of 10 counts this will be no great problem, the reason for applying the efficiency correction to the monitor is to avoid the  $\text{esd}=\sqrt{\text{counts}}$  when counts is very small problem. Nasty fudges like smoothing the monitor spectrum seem a bit dangerous as if there is any instability in the beam we would want to be able to divide this out, so we don't do anything like that. For these programs the monitor is going to have to be good.

So, we just multiply the monitor spectra by the efficiency and sum channels and monitor spectra to give a single, all channels and one monitor array (and monitor esd array)

Let's take things from the point where all the individual channel spectra are available (in **ascan**). We call a routine called **calibsum** which is eventually going to process that into an array called **hist** for us. This calls **effic** to get the efficiencies, **reporteffic** to write out a temp.res file and any other information. Then we call **sumthem** to combine the data from the channels together in a **sumdata** array and **normerr** to figure out the scale factor and make the final histogram. Optionally we can write out the ascan array via **bcmfile** at this point for later re-interpretation with more clever statistical methods.

---

<sup>2</sup>Sometimes they really do that!

```

< calibsum 43a > ≡
    subroutine calibsum
    use specfiles
    use rebin
    integer(kind=4) :: m,n,ierr
    real(kind=8),dimension(NCHAN) :: tmult,tmulterr
    call effic(n,m,tmult,tmulterr)
    call reporteffic(n,m,tmult,tmulterr)
    ierr=0
    if(.not.allocated(sumdata))allocate(sumdata(3,npts),stat=ierr)
! 3 ! cts, mon, e(mon)
    if(ierr.ne.0)then
        write(*,'(a)') 'Memory allocation error in calibsum'
        stop
    endif
    if(.not.allocated(hist))allocate(hist(2,npts),stat=ierr)
! final signal and esd
    if(ierr.ne.0)then
        write(*,'(a)') 'Memory allocation error in calibsum'
        stop
    endif
    call sumthem      ! combine detectors in sumdata array
    if(medianofchannels)call medianchannels
    call normerr      ! determine error bars and fill in hist
    call checkdets    ! check ascan matches hist
1  if(zapping)then
    write(*,*)
    write(*,'(A,F8.5,A)') "After zapping a sigma level ",zap," ... "
    call sumthem      ! combine detectors in sumdata array
    call normerr      ! determine error bars and fill in hist
    call checkdets    ! check ascan matches hist
    endif
    if(nzap.gt.0)goto 1
2  if(superzap .and. nsuperzap.gt.0)then
    write(*,*)
    write(*,'(A,F8.5,A)') "After superzapping at level ", superzaplevel,"..."
    call sumthem      ! combine detectors in sumdata array
    call normerr      ! determine error bars and fill in hist
    call checkdets    ! check ascan matches hist
    endif
    if(nsuperzap.gt.0)goto 2
    return
end subroutine calibsum

```

Fragment referenced in 56.

A routine for working out the detector efficiencies is provided by the **effic** subroutine. It goes through the **ascan** array from the **rebin** module and works out how many points there are where all of the channels are overlapping. These points are then all summed together (effectively into a single bin) and the efficiency is calculated in order to have all channels giving the same signal for this megabin. (FIXME) alternative for when there is no overlap of nine channels (FIXME) Summation must also be in some way aware of excluded channels to determine efficiencies if one of the channels is unplugged. - test this, I think it is there bar the esds.

```

< effic 43b > ≡
    subroutine effic(n,m,tmult,tmulterr)
    use rebin ! mult, multerr, nchan, ascan, tempres, npts
    integer(kind=4),intent(inout)::n,m
    real(kind=8),intent(inout),dimension(NCHAN)::tmult,tmulterr

```

```

integer(kind=4) :: i,j,k,nex
real(kind=8), dimension(4,NCHAN) :: signal
real(kind=8) :: sumsig,sumsige
signal=0.0d0; n=0
write(*,'(a)') 'Determining detector efficiencies'
nex=0; do j=1,nchannel; if(logexdet(j).eq.1)nex=nex+1; enddo
do i=1,npts
  k=0
  do j=1,nchannel      ! must have more than 1 monitor in bin to use
    if(ascan(2*j,i).gt.1.0d0)then
      if(ascan(2*j-1,i).gt.minrenormsig)k=k+1
    endif
  enddo
  if(k.eq.(nchannel-nex))then
    n=n+1
    do j=1,nchannel    ! signal is a big bin for all overlapping points
      if(logexdet(j).eq.1)cycle
! If a channel is excluded for the whole tth range in a file then this will fail?
      signal(1,j)=signal(1,j)+ascan(2*j-1,i) ! counts
      signal(2,j)=signal(2,j)+ascan(2*j,i) ! monitors
    enddo
  endif
enddo
if(n.gt.0)then
  do j=1,nchannel      ! Normalise Signal and get the error bar on it
    if(logexdet(j).eq.1)cycle
    signal(3,j)=signal(1,j)/signal(2,j)
    signal(4,j)=signal(3,j) *                                     &
    & sqrt(1.0d0/signal(1,j)+1.0d0/signal(2,j))
  enddo
  sumsig=sum(signal(3,:)) ! sum of all signals in overlap region
  sumsige=0.0d0
  do i=1,nchannel      ! Error in sum of all signals
    sumsige=sumsige+signal(4,i)*signal(4,i)
  enddo
  sumsige=sqrt(sumsige)
  do j=1,nchannel      ! Finally get the channel efficiencies
    if(logexdet(j).eq.1)then
      tmult(j)=1.0d0; tmulterr(j)=0.0d0
    else
      tmult(j)=(real((nchannel-nex),8)*signal(3,j)/sumsig)
      tmulterr(j)=tmult(j)*sqrt((signal(4,j)/signal(3,j))**2 +      &
      & (sumsige/sumsig)**2)
    endif
  enddo
  if(.not.tempres)then  ! Copy these to rebin module if no temp.res file
    mult=tmult
    multerr=tmulterr
    m=1
  else ! tempres
    m=2                ! Flag the need to print the temporary (unused vals)
  endif ! tempres
else ! n.gt.0
  m=1                  ! no overlap so only one to print
  if(.not.tempres)then
    mult=1.0d0        ! Make sure defaults are always 1.0d0
    multerr=0.0d0
  endif
endif ! if n.gt.0
return

```

end subroutine effic

◇

Fragment referenced in 56.

Reporting the results of the calibration is relatively involved as we handle a variety of possible situations, depending whether channel overlap was found and whether a temp.res file was found. The intention is to print numbers for comparison if a temp.res file is used (in case it doesn't agree with the current dataset). Also to indicate where the efficiencies are coming from and create a new temp.res file when necessary. The arguments **n** and **m** indicate the number of points used to find the efficiency and whether or not some efficiencies were already available.

⟨reporteffic 45⟩ ≡

```

subroutine reporteffic(n,m,tmult,tmulterr)
use rebin ! mult, multerr, NCHAN
integer(kind=4),intent(in) :: n,m
real(kind=8),intent(in),dimension(nchan)::tmult,tmulterr
integer(kind=4) :: i
logical :: trex
if(n.gt.0)then
  write(*,'(a,i9,a)')'Channel efficiencies found from ',n,      &
  & ' points where all detectors overlap'
  if(m.eq.1)then
    write(*,'(a)')'Det      Offset      Effic      <Effic>'
    do i=1,NCHANNEL
! i-1 instead of i to start at channel zero
      write(*,'(i3,3(1X,F10.7))')i-1,offset(i),mult(i),multerr(i)
    enddo
  endif ! m.eq.1
  if(m.eq.2)then
    write(*,'(a)')'Efficiencies from temp.res file,'//      &
    & ' the values found now are compared'
    if(.not.renorm)then ; write(*,'(a)')      &
    & 'Det      Offset      Effic      <Effic> current unused values'      &
    else ; write(*,'(a)')      &
    & 'Det      Offset      Old Eff      Old <E>      New Eff      New <E>'      &
    endif
    do i=1,NCHANNEL
! i-1 instead of i to start at channel zero
      write(*,'(i3,5(1X,F10.7))')i-1,offset(i),mult(i),multerr(i),      &
      & tmult(i),tmulterr(i)
    enddo
    if(renorm)then ; mult=tmult ; multerr=tmulterr ; endif
  endif ! m.eq.2
  inquire(file='temp.res',exist=trex)
  if(.not.trex .or. renorm)call tempreswrite
  if(renorm) renorm=.false. ! for sumall - can only renorm once
else ! n.gt.0
  write(*,'(a,$)')'No detector overlap found, efficiencies'
  if(tempres)then
    write(*,'(a)')' from temp.res file'
  else
    write(*,'(a)')' are probably wrong'
    call tempreswrite
  endif
  write(*,'(a)')'Det      Offset      Effic      <Effic>'
  do i=1,NCHANNEL
! i-1 instead of i to start at channel zero
    write(*,'(i3,3(1X,F10.7))')i-1,offset(i),mult(i),multerr(i)

```

```

        enddo
    endif
    return
end subroutine reporteffic

```

◇

Fragment referenced in 56.

A short subroutine to write out the temp.res file....

```

⟨ tempreswrite 46 ⟩ ≡
    subroutine tempreswrite
        use rebin ! offset, nchan, mult
        integer(kind=4) :: ierr, i
        open(unit=16,file='temp.res',status='UNKNOWN',access='SEQUENTIAL',&
        & form='FORMATTED',iostat=ierr)
        if(ierr.ne.0)then
            write(*,'(a)')'Couldn't open temp.res to write!!!'
            return !! bugs out
        endif
        write(*,'(a)')'Created temp.res file'
        do i=1,nchannel
            write(16,1000)offset(i),mult(i),multerr(i),0.0d0
        enddo
1000    format(4(f11.8,' ',''))
        close(16)
        tempres=.true. ! Should exist and be readable now
        return; end subroutine tempreswrite

```

◇

Fragment referenced in 56.

Once the detector efficiencies have been determined the data in the **ascan** array can be combined into the **sumdata** array using the efficiency numbers stored in **mult**. We will assume that the error on the efficiency is actually zero for determining the final error bars on the histogram. Strictly this is incorrect, but unless someone can figure out the details we'll just ignore the problem. The efficiency errors are correlated, in some way, with the errors we are after, so simple formulae will not suffice. **sumdata** column 1 will hold the summed counts, column 2 holds the summed monitor×efficiency and column 3 holds the error on the monitor column. The corrected monitor is given by

$$m = ce$$

where  $c$  is the number of monitor counts,  $m$  is the corrected monitor signal and  $e$  is the efficiency for the channel. *Note that this is a product, the original version screwed things up here. If the detector has a low efficiency, that means the monitor was effectively less.* The error for the monitor of one channel (after efficiency correction) is:

$$\langle m \rangle = e\sqrt{c}$$

if the efficiency is assumed to have zero error. However, if the error in the efficiency is propagated through as well then we get

$$\begin{aligned}
 \left(\frac{\langle m \rangle}{m}\right)^2 &= \left(\frac{\sqrt{c}}{c}\right)^2 + \left(\frac{\langle e \rangle}{e}\right)^2 \\
 \langle m \rangle &= ce\sqrt{\frac{1}{c} + \left(\frac{\langle e \rangle}{e}\right)^2} \\
 \langle m \rangle &= \sqrt{ce^2 + (c\langle e \rangle)^2} \\
 \langle m \rangle &= \sqrt{c(e^2 + c\langle e \rangle^2)}
 \end{aligned}$$



This clearly goes towards  $\sqrt{c}$  as the efficiency tend towards one, with something added if the efficiency is not known precisely. So the error on the sum of all of these monitors is:

$$\begin{aligned} \langle m_{tot} \rangle^2 &= \sum_{i=1, nchan} \langle m_i \rangle^2 \\ \langle m_{tot} \rangle^2 &= \sum_{i=1, nchan} c (e^2 + c \langle e \rangle^2) \end{aligned}$$

The **sumthem** routine is supposed to perform that calculation, putting the results into the **sumdata** array.

```

< sumthem 47 > ≡
  subroutine sumthem
  use rebin ! for ascan and mult
  integer(kind=4) :: i,j
  sumdata=0.0d0 ! the big array where the sum of channels goes
  do i=1,npts
    do j=1,nchannel
      ! counts
      sumdata(1,i)=sumdata(1,i)+ascan(2*j-1,i)
      ! normalised mon = m1*e1 + m2*e2 + ...
      sumdata(2,i)=sumdata(2,i)+ascan(2*j,i)*mult(j)
      ! emon**2 =
      sumdata(3,i)=sumdata(3,i)+
      !      c      (      e**2      +      c      *      <e>**2      )      &
      &      ascan(2*j,i) * ( mult(j)**2 + ascan(2*j,i)*multerr(j)**2)
      enddo
      ! emon=sqrt(emon**2)
      sumdata(3,i)=dsqrt(sumdata(3,i))
      enddo
      return
    end subroutine sumthem
  ◇

```

Fragment referenced in 56.

The final data is going to be signal and esd on signal where signal is (summed counts)/(summed normalised monitor counts). We derive the esd on the final signal here. During summation we calculated the summed counts for the detectors and monitor and propagated an error bar on the summed monitor counts. The final signal is

$$signal = counts/monitor$$

(FIXME - move this up a bit) We use the following formulae for determining errorbars on uncorrelated random variables.

- If  $z = x + y$  or  $z = x - y$  then  $e_z^2 = e_x^2 + e_y^2$
- If  $z = x/y$  or  $z = xy$  then  $e_z^2/z^2 = (e_x/x)^2 + (e_y/y)^2$
- If  $z = ax$  where  $a$  is a constant  $e_z = ae_x$

The error bar on the counts is estimated as  $\sqrt{counts + \alpha}$ , with  $\alpha$  normally being 0.5 but this is to be a modifiable parameter For 0 counts the errorbar on the counts would be zero, which is appears to be silly. Hence the addition of a parameter,  $\alpha$ , tries to avoid this problem. There are some papers and Bayesian arguments which need to be referenced here (FIXME) For counts of more than about 10 the difference is unimportant, but in cases where there are wide expanses of very low background (and small stepsize) the user can run into problems with  $\chi^2$  significantly less than one, as the  $\sqrt{counts}$  stuff breaks down. A warning message could be issued if a lot of the channels have this problem.

Taking the counts from the **sumarray** and putting them into **hist** means just dividing the counts by the monitor counts to get the hist array. We set unobserved points to be negative and let the writing out

routines worry about what that means (or not write them out!) If we use  $y$  and  $\langle y \rangle$  to be the signal and error in the **hist** array, and  $c$ ,  $m$  and  $\langle m \rangle$  to be the contents of the sumdata array then we have:

$$y = c/m$$

and for the error:

$$\left(\frac{\langle y \rangle}{y}\right)^2 = \left(\frac{\langle c \rangle}{c}\right)^2 + \left(\frac{\langle m \rangle}{m}\right)^2$$

Since:

$$\langle c \rangle = \sqrt{c + \alpha}$$

we can substitute for  $\langle c \rangle$  and  $y$  giving:

$$\begin{aligned} \langle y \rangle^2 &= \frac{c^2}{m^2} \left( \frac{c + \alpha}{c^2} + \frac{\langle m \rangle^2}{m^2} \right) \\ \langle y \rangle^2 &= \frac{c + \alpha}{m^2} + \left( \frac{c \langle m \rangle}{m^2} \right)^2 \end{aligned}$$

The constant  $\alpha$  prevents the error ever becoming zero. Note that if no counts are observed at all then the error is actually undefined as  $\langle y \rangle / y$  would involve a divide by zero.

A quick example would be  $10^4$  counts in the bin, with  $10^6$  monitor counts and the efficiencies having all been 1, so that the error on the counts is  $10^2$  and on the monitor is  $10^3$  (ignoring alpha). In this case the signal is 0.01, the first term of the error is  $10^4/10^{12}$  and the second term is  $10^{20}/10^{24}$  giving  $10^{-8} + 10^{-4} \sim 10^{-4}$ , which is roughly what we expect from the monitor being a constant scale factor. We place a certain amount of faith in the equations above and floating point arithmetic here.

```

< normerr 48 > ≡
  subroutine normerr
    ! alp, hist & sumdata available as member of summation
    use rebin ! npts
    real(kind=8):: msq, s
    integer(kind=4) :: i, n, n3
    n=0; s=0.0d0; n3=0
    do i=1,npts
      if(sumdata(2,i).gt.minmon)then ! needs at least 1. mon count
        hist(1,i)=sumdata(1,i)/sumdata(2,i) ! correct
        msq=(sumdata(2,i)*sumdata(2,i)) ! always gt zero
        hist(2,i)=sqrt(
& (sumdata(1,i)+alp)/msq + (sumdata(1,i)*sumdata(3,i)/msq)**2)
        n=n+1
        s=s+hist(1,i)**2/hist(2,i)**2
        if(hist(1,i)/hist(2,i).lt.3.0d0)n3=n3+1
      else
        hist(1,i)= 0.0d0 ! unobserved regions are filled with nonsense
        hist(2,i)=-1.0d0
      endif
    enddo
    write(*,1000)100.0d0*sqrt(real(n,8)/s),n3,n
1000 format('R_exp = ',f7.3,' with ',i7,
& ' pts having I/<I> less than 3, from ',i7,' pts obs')
    end subroutine normerr

```

Fragment referenced in 56.

After the **hist** array has been filled in we can determine whether or not each of the individual channels is in agreement with the summed up scan. We just need to compute the final histogram on the basis of only using one channel (and it's error bar) and then take the difference between this single channel scan and the sum total, then see if they agree to within  $3\sigma$ . We can write out at the end the number of points

which differ by more than some number of sigmas (for each channel), or alternatively decide on some criterion for a problem and only print anything then. If a problem is found we should create a diagnostic plot which can aid the user for examining the details.

Checking that the various detectors all agree with the final summed dataset is a useful diagnostic. Sometimes on BM16, when using the cryostat at hard energies, some of the channels contained a significant background compared to others. The programs should automatically detect this kind of problem and inform the user that things are going wrong. Ideally the instrument should never produce such poor quality data, nevertheless, if it ever does we should be ready to catch it and work around the problem until the hardware is fixed. The **checkdets** routine will perform this check by forming the signal from each detector channel separately and comparing it to the contents of the **hist** array. Whether or not the two agree depends on the difference between them, scaled to the esd on the difference. For nine channels, the signal  $y$  as estimated by a single channel with counts  $c$ , monitor  $m$  and efficiency  $e$  is given by:

$$y = \frac{c}{em}$$

so just a scale factor of  $1/e$ , the esd on  $c/m$  is

$$\left( \frac{\langle c/m \rangle}{c/m} \right)^2 = \left( \frac{\sqrt{c+\alpha}}{c} \right)^2 + \left( \frac{\sqrt{m}}{m} \right)^2$$

$$\langle c/m \rangle = \sqrt{\frac{c+\alpha}{m^2} + \frac{c^2}{m^3}}$$

so that

$$\langle y \rangle = \frac{1}{e} \sqrt{\frac{c+\alpha}{m^2} + \frac{c^2}{m^3}}$$

Those equations are encapsulated in the following fortran which is used both here and later in the diagnostic plot output routine. The error in the efficiency is neglected in calculating the individual channels error bars. It has already been incorporated into the overall errorbar on the **hist** array, but clearly it's more important here. However if a signal was calculated from only one channel then it's error in efficiency is logically zero - as that correlates with the output units. That is as unclear to me as it is to you, anyway, we forget about the error in efficiency for this calculation.

```

< getwd2 49 > ≡
      y=ascan(2*j-1,i)/ascan(2*j,i)/mult(j)    ! OK
      ! assumes zero error in the efficiency
      ey2= ((ascan(2*j-1,i)+alp)/ascan(2*j,i)**2
            &      + ascan(2*j-1,i)**2/ascan(2*j,i)**3)/mult(j)**2
      wd2=(y-hist(1,i))**2/(ey2+hist(2,i)**2)

```

Fragment referenced in 50, 64b, 66.

```

< checkdets 50> ≡
  < ctchan 51>
    subroutine checkdets
      use rebin ! for ascan and mult arrays
      integer(kind=4) :: i, j, ipt3s, ipt6s, iapts, myctchan
      real(kind=8) :: c,y,ey2,wd2 ! chi2,sig,error^2 and wtd diff^2
      real(kind=8),dimension(nchan)::sz ! superzap array
      ipt3s=0 ; ipt6s=0; c=0.0d0; iapts=0; nzap=0; nsuperzap=0;
      do i = 1, npts
        myctchan=ctchan(i)
        if(myctchan.ge.2)then
          if(superzap)sz=-1.0
          do j = 1, nchannel
            if(ascan(2*j,i).gt.1.0d0)then ! need one monitor count to bother
1000  < getwd2 49>
              if(superzap)sz(j)=y !
              c=c+wd2
              iapts=iapts+1
              if(wd2.gt.9.0d0) ipt3s=ipt3s+1 ! 3s^2=9.0
              if(wd2.gt.36.0d0)ipt6s=ipt6s+1 ! 6s^2=36.0
              if(wd2.gt.zap*zap .and. zapping)then
                ascan(2*j,i)=0.0 ! set mon and det to zero for zapped points
                ascan(2*j-1,i)=0.0
                nzap=nzap+1
              endif
            endif
          enddo
          ! now check with superzap
          if(superzap .and. myctchan.ge.3)call superzapem(sz,ascan(:,i),nchannel)
          endif
        enddo
        if(iapts.gt.0) then
          c=c/real(iapts,8)
          write(*,1000)c,iapts
1000  format('Reduced chi**2 for channel merge =',F9.4,' from ',
            & i10,' pairs of pts')
          write(*,1001)100.0d0*real(ipt3s)/real(iapts),
            & 100.0d0*real(ipt6s)/real(iapts)
1001  format(f5.2,'% differ by >3 sigma, ',f7.4,
            & '% by >6 sigma (ideally 0.04% and 0.0000%)' )
          i6s=ipt6s ! module variable
          if(zapping)then
            write(*,'(A,I8,A)')'Zapping ',nzap,' points'
          endif
          if(superzap)write(*,'(A,I8,A)')'Superzapping ',nsuperzap,' points'
          else
            write(*,*)'No channel overlap, whatsoever, was found'
            i6s=0
          endif
        return
      end subroutine checkdets

```

Fragment referenced in 56.

We introduce the concept of a statistic to measure the quality of the data merging here. It is defined as follows... for each point in each separate channel, take the difference between that point and the summed data at that point. Clearly this will break down rather badly when only one detector is contributing to a particular point, so we ignore points where that is the case. The **ctchan** function just gets the number of active channels for us.

```

⟨ ctchan 51 ⟩ ≡
    integer(kind=4) function ctchan(i)
    use rebin
    integer(kind=4),intent(in) :: i
    integer(kind=4) :: j
    ctchan=0
    do j=1,nchannel
        if(ascan(2*j,i).gt.minmon)ctchan=ctchan+1
    enddo; return; end function ctchan

```

Fragment referenced in 50. ◇

If this is greater than 2 we proceed to calculate the difference between the signal estimated from a signal channel and the signal from the summed data. The error on this is also calculated and we sum up the difference squared over the error squared to get a  $\chi^2$  type statistic. This then reduced by the number of datapoints which went into it.

## 6.2 Taking the median of the channels

The code in **sumthem** averages the different channels by summing up the counts in the different channels and computes the corresponding error bars. With outliers appearing in the data sometimes it could be good to try to take the median of the contributing channels rather than the mean. The error bars seem only to be something that can be fudged by taking the old values?

This means that the **sumdata** array will hold the median value multiplied by the number of channels and weighted by the monitor... we had `sumdata(0)=sum(ascan(0))` and replace this with the median column...

$\langle \text{medianchannels } 52 \rangle \equiv$

```

! netlib code
 $\langle \text{dsort } 136b \rangle$ 
! netlib code
subroutine medianchannels
use rebin
integer(kind=4) :: i,j,k
real(kind=8), dimension(nchan) :: signal ! local copy
integer(kind=4), dimension(nchan) :: iactive ! local copy
sumdata=0.0d0 ! the big array where the sum of channels goes
do i=1,npts
  k=0
  do j=1,nchannel
    if(ascan(2*j,i).gt.0)then
      k=k+1
      signal(k)=ascan(2*j-1,i)/(ascan(2*j,i)*mult(j))
      iactive(k)=j ! which channel was active
    endif ! signal is computed
  enddo
  if(k.eq.0)then ! nothing recorded at all
    sumdata(1,i)=0.
    sumdata(2,i)=0.
    sumdata(3,i)=0.
  endif
  if(k.eq.1)then ! only one channel active
    j=iaactive(1)
    sumdata(1,i)=ascan(2*j-1,i)
    sumdata(2,i)=ascan(2*j,i)*mult(j)
    sumdata(3,i)= &
&    ascan(2*j,i) * ( mult(j)**2 + ascan(2*j,i)*multerr(j)**2)
    sumdata(3,i)=dsqrt(sumdata(3,i))
  endif
  if(k.gt.1)then ! more channels active
    CALL DISORT(SIGNAL,IACTIVE,K,2) ! 2 means increasing order
    j=iaactive((k+1)/2) ! median channel
    sumdata(1,i)=ascan(2*j-1,i)*k
    sumdata(2,i)=ascan(2*j,i)*mult(j)*k
    sumdata(3,i)= k* &
&    ascan(2*j,i) * ( mult(j)**2 + ascan(2*j,i)*multerr(j)**2)
    sumdata(3,i)=dsqrt(sumdata(3,i))
  endif
enddo
return
end subroutine medianchannels

```

◇

Fragment referenced in 56.

### 6.3 Saving the ascan array

Due to the problems which sometimes occur with noise or crosstalk between channels it is sometimes desirable to massage the data when combining them together. Rather than implement all possible algorithms for carrying out these procedures, we dump the problem onto someone else by offering to write out a "bcm" file, which is the binned counts and monitor.

The file will contain a header line describing the data which are to be written out and then lines with the data themselves. We only need write as many columns as are used and we should label them in a generic way.

```

⟨bcmfile 53⟩ ≡
  subroutine bcmfile(n)
    use rebin ! for ascan and mult
    use summation
    use specfiles
    implicit none
    integer(kind=4) :: n,ma0
    character(len=WORDLENGTH) :: c
    integer(kind=4) :: ilow,ihigh,i,j
    call filext('.bcm',n)
    open(unit=ioutunit,status='UNKNOWN',file=outfile)
    ilow=getfirstpoint(hist,2,npts,2) ! range of real points
    ihigh=getlastpoint(hist,2,npts,2)
    c=columnlabels(1)
    ma0=whichcolumn(FIRSTDET) ! FIXME deal with absent dets
    write(ioutunit,'(a,2x,a,$)')'#',c(1:len_trim(c))
    do i=0,nchannel-1
      if(logexdet(i+1).eq.1)cycle ! skip if excluded completely
      c=columnlabels(ma0+i)
      write(ioutunit,'(2X,a,2x,a,$)')c(1:len_trim(c)),
& c(1:len_trim(c))//'_mon'
    enddo
    write(ioutunit,*) ! end of line
    do i=ilow,ihigh
      write(ioutunit,'(G14.8,2X,$)')bincen(i)
      do j=1,nchannel
        if(logexdet(j).eq.1)cycle ! skip if excluded completely
        write(ioutunit,'(2(G14.8,2X),$)')ascan(2*j-1,i),ascan(2*j,i)
      enddo
      write(ioutunit,*) ! end of line
    enddo
    close(ioutunit)
    write(*,'(a)')'Wrote '//outfile(1:len_trim(outfile))
    return
  end subroutine bcmfile

```

Fragment referenced in 72.

## 6.4 Superzapem

Yet another approach to zapping. The idea is that when there are at least 3 channels overlapping, we can calculate the mean and variance of the 3 independent signals and compare it to the error bar derived from the counting statistics. If the channels "agree" then all will be fine, but if the channels are markedly different we have a problem. We eliminate the highest value (it is always too much that is a problem) and recalculate. There must always be 2 channels left. The "improvement" after deleting a channel is measured in terms of the percentage reduction in the n channel esd versus the n-1 channel esd, where the esds are normalised to the "real esd" in each case.

```

⟨ superzapem 54 ⟩ ≡
  subroutine superzapem(sz,as,n)
    integer(kind=4), intent(in) :: n
    real(kind=8),dimension(n),intent(in) :: sz    ! signal
    real(kind=8),dimension(2*n),intent(inout) :: as ! ascan array
    integer(kind=4) :: i, ib, nc
    real(kind=8) :: s, ss, mean1, mean2, esd1, esd2, pc
    ! make mean and esd from sz
    nc=0;s=0.0;ss=0.0
    ib=1; pc=0.0
    do i = 1,n
      if(sz(i).ge.0.0)then
        s=s+sz(i)
        ss=ss+sz(i)*sz(i)
        nc=nc+1
        if(sz(i).gt.sz(ib))ib=i
      endif
    enddo
    mean1=s/nc
    esd1=(ss-mean1*mean1)/nc
    ! now without ib
    s=0.0;ss=0.0;nc=0
    do i = 1,n
      if(sz(i).ge.0.0 .and. i.ne.ib)then
        s=s+sz(i)
        ss=ss+sz(i)*sz(i)
        nc=nc+1
      endif
    enddo
    mean2=s/nc
    esd2=(ss-mean2*mean2)/nc
    if(esd1.gt.0.0) pc=(esd1-esd2)/esd1 ! should use mean info too.
    if(pc.gt.superzaplevel)then ! zap the ib channel
      as(2*ib)=0.0
      as(2*ib-1)=0.0
      nsuperzap=nsuperzap+1
    endif
    return
  end subroutine superzapem

```

Fragment referenced in 56.

## 6.5 Combining data from different scans

We should provide some statistics to determine whether or not the various scans which are being combined are in agreement with each other. The `id31sum` program will just blindly add the scans together into the **ascan** array for now. A second pass through the SPEC file will allow this to be accomplished, or we could insist that the scans are summed separately (binem style) and then those scans are compared to the final summed total.

## 6.6 Pattern scaling

Deciding on a final multiplicative factor to give the units of the final dataset remains fairly vague. Some suggestions for deciding what the final scale factor should be are:

- The total number of counts in the final histogram should equal the total number of counts detected.
- The height of the highest peak in the final histogram should have the correct number of counts.



- The units should be counts per monitor count, with the data being dimensionless and just looking smoother as you count longer (so all scans are comparable)
- The counts should not be scaled at all, instead of writing esd's out you should supply a scaling column, s, such that the signal is counts/s and the esd is  $\sqrt{\text{counts}}/\text{s}$

Each of these suggestions has its own advantages and disadvantages. For looking at data during online collection we should implement a method which puts out the raw counts, so people can see how long they need to spend on different regions of the dataset. For final datasets with binem type output must be dimensionless (counts per monitor count) so that many datasets can be directly compared. For final summed refineable datasets we think that the first option is attractive, although the second clearly is appealing as well. Whatever happens, if you count for much longer at high angles and then plot normalised data you will not appreciate the effort you have made on the high angle data. Routines to implement each of the methods suggested above are given here. They take the data from **sumdata** and **hist**, calculate the appropriate scale factor and applying it, if necessary. The default format in **hist** is counts per monitor count, so in that case nothing needs to be done.

For scaling the total counts to be the same as the original total number of counts the **scaltot** routine uses the information in the **sumdata** array to rescale the **hist** array.

```

< scaltot 55a > ≡
      subroutine scaltot
      ! uses hist and sumdata arrays in summation module
      real(kind=8) :: factor, sum1, sum2
      sum1=sum(sumdata(1,:)) ! sum up real counts
      sum2=sum(hist(1,:))   ! sum up normed counts
      factor=sum1/sum2      ! determine multiplier
      hist=hist*factor      ! rescale hist array
      return
      end subroutine scaltot

```

Fragment referenced in 56.

To scale the height of the highest peak as being correct the **scalpk** routine just finds the highest peak in the normalised data and then gets the scale factor by comparing this to the original counts.

```

< scalpk 55b > ≡
      subroutine scalpk
      ! uses hist and sumdata arrays in summation module
      real(kind=8) :: factor
      integer(kind=4), dimension(3) :: i
      i=maxloc(sumdata,dim=2) ! get index of highest number of counts
      factor=sumdata(1,i(1))/hist(1,i(1)) ! get scaling factor
      hist=hist*factor
      return
      end subroutine scalpk

```

Fragment referenced in 56.

To write out counts and a scaling column the current scaling of the **hist** array will be irrelevant. So long as signal and esd on signal are both available the numbers can be determined when we come to write the **vct** file.

## 6.7 A module for the summation stuff

Drawing together the routines for summing and normalising the contents of the **ascan** array in the **rebin** module to reach final refineable data sets. A space for the summed array is also needed, it needs at three columns - for the summed data and summed monitor spectrum and an error on the summed monitor

(for propagating the detector efficiency into the error bars). To create this array we will need to use the detector efficiencies. `eff` and `effsqrd` hold information for when we do the detector normalisation, `mult` holds the multiplier numbers which are needed at the final summing stage. **hist** will hold the final scaled counts and error bars.

```

< summation 56 > ≡
    module summation
    integer(kind=4):: isc=0, i6s=0 ! which scale factor, 6sig pts
    integer(kind=4):: nzap=0, nsuperzap=0 ! for zapping
    real(kind=8),allocatable :: sumdata(:,:),hist(:,:)
    real(kind=8) :: alp=0.5d0 ! Bayesian zero counts fudge
    real(kind=8) :: scalinp=1.0d5 ! Scale factor for .inp files
    real(kind=8) :: zap=6.0d0 ! level for outlier elimination (median filter?)
    real(kind=8) :: superzaplevel=1.0d0 ! level for zinger elimination
    logical :: renorm=.false. ! to update temp.res file efficiencies
    logical :: zapping=.false.
    logical :: superzap=.false.
    logical :: medianofchannels=.false.
    contains
    < calibsum 43a >
    < effc 43b >
    < sumthem 47 >
    < normerr 48 >
    < checkdets 50 >
    < reporteffc 45 >
    < tempreswrite 46 >
    < scaltot 55a >
    < scalp 55b >
    < superzapem 54 >
    < medianchannels 52 >
    end module summation

```

◇

Fragment referenced in 40a, 57b, 86b, 89, 91.

## 7 Detector offset calibration

We'll take a rebinned scan (or scans) and assume the rebinning introduces a negligible extra broadening in the peaks (small bins!). Then compute the overlap integral as a function of the offset of the scans with respect to each other. Figure 5 shows graphically the function to be computed.

Since the different channels might not be correctly scaled with respect to each other, the overlap will be normalised to the sum of the two individual channels over the given range. The simplest way to do this seems to be to compute a correlation integral as a function of the detector offset. Assuming we have filled in the `ascan` array (from module `rebin`) with a set of counts and monitor counts we will just need a function to compute this integral for a given offset. The two theta bin centers and high and low limits will be given by the hard coded offsets in the program. If we subsequently decide these offsets are wrong then the scan will still need to be rebinned from the original data. This is a calculation which should only be done once in a while (the offsets are physically fixed on the instrument anyway).




Figure 5: Illustration of the overlap integral to be computed as a function of detector offset. The grey area represents the minimum of the two patterns.

```

< crrfun 57a > ≡
    real(kind=8) function crrfun(j)
    ! Computes correlation function offset of n2 versus n1
    use rebin ! contains ascan array and npts, ibin, tthhb etc
    implicit none
    integer(kind=4), intent(in) :: j
    integer(kind=4) :: n1, n2
    common/cfcom/ n1, n2
    real(kind=8) :: s1, s2
    integer(kind=4) :: n, i
    crrfun=0.0d0; n=0
    do i=1,npts
        if((i+j).gt. 0 .and. (i+j).lt.npts)then
            if(ascan(2*n1,i).gt.0.0d0 .and. ascan(2*n2,i+j).gt.0.0d0) then
                n=n+1
    ! signal in n1 and n2, n22
                s1= ascan(2*n1-1,i) /ascan(2*n1,i)
                s2= ascan(2*n2-1,i+j) /ascan(2*n2,i+j)
    ! crrfun is product of the two signals
                crrfun=crrfun + s1*s2
            endif; endif
        enddo
        crrfun=crrfun/real(n,8)
    return
    end function crrfun

```

Fragment referenced in 57b.

As a test we should read a SPEC file in, rebinning it along the way and then compute the correlation function for each detector versus channels four. Use a sharp dataset with good statistics and small binsize for this to be effective.

```

"id31offsets.f90" 57b≡
    < specfiles 11a >
    < rebin 28 >
    < summation 56 >
    < outputfiles 72 >
    < processscan 33 >

```

```

<useroptions 87a>
<tidyup 98>
<crrfun 57a>
<helpmsg 87b>
<id31offsetsmsg 59>
  program id31offsets
    ! read and bin the scan(s)
    ! calculate crrfun for each channel on a small tth range
    use specfiles
    use rebin
    !   use summation ! for user options - FIXME - module deps
    use useroptions
    use outputfiles
    integer(kind=4) :: n1, n2
    common/cfcom/ n1, n2
    integer(kind=4)::i, j, ioffmin
    real(kind=8) :: off, crr, crrmin, d1, d2, d21
    real(kind=8) :: offm1, offp1, crrm1, crrp1
    real start_time, end_time
    real(kind=8) crrfun
    external crrfun
    call cpu_time(start_time)
    call getcmdline                      ! Get users options
    ! fill out names of monitor, tth, first and last columns
    if(snbl)then
      MONITORCOL="Mon";FIRSTDET="Det1";LASTDET="Det6";TWOTTH="TwoTheta"
      NCHANNEL=6
    else
      FIRSTDET="MA0";LASTDET="MA8";TWOTTH="2_theta"
      NCHANNEL=9
    endif
    call initialiserebin                ! Allocate space and check parameters
    call getfile                        ! Opens the spec file
    call openlogfile                   ! for counter details to be logged when binning
    write(*,'(a,$)')'Binning scan '
1   i=nextscan()                      ! Gets next scan to be processed
    if(i.gt.0) then
      call processscan(i)              ! sums into ascan array
      if(ispecerr.eq.-1)then
        ispecerr=0
        goto 1                        ! next only if current was OK
      endif
    endif
    write(*,*) ! after non advancing io list of scans binned
    ! Now for the actual channel calibration - data in ascan array of rebin
    open(unit=20,file='offsets.out',status='UNKNOWN')
    do i=1,9 ! all nine channels
      if(logexdet(i).eq.1)then
        write(*,1000)0.0d0,offset(i) ! dont compute offset for channel
        write(20,1000)0.0d0,offset(i)
        cycle
      endif
      if(i.eq.5) then
        write(*,1000)0.0d0,offset(i) ! compute the new offset for channel
        write(20,1000)0.0d0,offset(i)
        cycle
      endif
      n1=i; n2=5; crrmin=0.0d0
      do j=-5,5,1
        crr=crrfun(j)

```

```

        if(crr.gt.crrmin)then; ioffmin=j; crrmin=crr; endif
    enddo
! Got the maximum of the function, now take the points on either side of it
    offm1=real(ioffmin-1,8)*step; offp1=real(ioffmin+1,8)*step
    off=real(ioffmin,8)*step; crr=crrfun(ioffmin)
    crrm1=crrfun(ioffmin-1); crrp1=crrfun(ioffmin+1)
! data three data points
! (x0,y0), (x1,y1), (x2,y2) with xi != xj for i,j in {0,1,2}
! set
! d1=(y1-y0)/(x1-x0)
! d2=(y2-y1)/(x2-x1)
! d21=(d2-d1)/(x2-x0)
! p2(x) = y0 + (x-x0)*(d1+(x-x1)*d21)
!       = x*d1 + x*(x-x1)*d21 - x0*x*d21 ! terms in x only
! dp2/dx = 0 at maximum (point of interest for us)
!       = d1 + 2 x d21 - x1 d21 - x0 d21 = 0
!       x(at max) = ((x1 + x0)*d21 - d1 )/2*d21
    d1=(crr-crrm1)/step; d2=(crrp1-crr)/step
    d21=(d2-d1)/(2.0d0*step)
    off = (d21*(offm1+off) - d1) / (2.0d0*d21)
    offset(i)=offset(i)-off
    write(*,1000)off,offset(i) ! compute the new offset for channel
    write(20,1000)off,offset(i)
1000  format(10g20.10)
    enddo
    open(file='temp.new',status='UNKNOWN',form='FORMATTED',unit=18)
    write(*,'(a)')'temp.new file to be created (copy over temp.res)'
    if(snb1)then
        do i=1,nchannel
            write(*,1001)offset(i)-offset(1),mult(i),multerr(i),0.0d0
            write(18,1001)offset(i)-offset(1),mult(i),multerr(i),0.0d0
        enddo
    else
        do i=1,nchannel
            write(*,1001)offset(i),mult(i),multerr(i),0.0d0
            write(18,1001)offset(i),mult(i),multerr(i),0.0d0
        enddo
    endif
1001 format(4(f11.8,' ',''))
    close(18)
    call tidyup ! Frees allocated memory
    call cpu_time(end_time)
    write(*,'(a,f6.2,a)')'Time taken was ',end_time-start_time,'/s'
    end program id31offsets

```

◇

⟨ id31offsetsmsg 59 ⟩ ≡

```

1000  format('example: ',a,' file.dat 0.01 1 10  [snbl] ',
&
&/'will process the scans 1 to 10 from file.dat with binsize 0.01'&
&/'to determine optimal detector offsets'&
&/' snbl if your data are in the Swiss Norwegian format'&
&/'New values will be written to a temp.new file. You need a very'&
&/'good dataset to attempt this. A small stepsize is recommended.'/&
&/'Use the optional argument step=x.xxxx to force a small stepsize'&
&/'We recommend checking the results from this program.')
    stop; end subroutine helpmsg

```

◇

Fragment referenced in 57b.

## 8 Writing output files

Some bits and pieces for outputting the results of our efforts. We derive filenames from the variable **filnam** in the **specfiles** module. If it ends in a .dat extension then the extension is replaced with something appropriate for the current output file, otherwise the extension is appended. Whether or not to write each kind of file will be determined by the **useroptions** module. A title, where it is needed can either come from the user options module, or better something like the #F line in the original SPEC file. Summary temperature information could also be included where allowed.

Variables to hold relevant parameters of interest are

```

< outputfilesvars 60a > ≡
      character(len=1024) :: outfile
      integer(kind=4) :: ioutunit=11, ilow, ihigh, itot, logfile=21
      logical :: gsas=.false. , spf=.false. , pds=.false., diag=.true.
      logical :: epf=.false. , xye=.true.
      character(len=1024) :: title

```

Fragment referenced in 72.

**ilow** and **ihigh** are the first and last bins in the summed array which are going to be output.

### 8.1 Various file formats

Most of this code is cut and pasted from Andy's **epfout** program, with adjustments to save reading in from an x,y,esd file.

```

< scantitle 60b > ≡
      if(is.gt.0)then
        write(str,'(i10)')is; str=adjustl(str)
        title=title(1:len_trim(title))//' Scan '//str(1:len_trim(str))
      endif

```

Fragment referenced in 61, 63, 64a.

Firstly we attack the GSAS format. We will need the number of points being written, the minimum angle and stepsize, both in centidegrees. For now we will not allow data from negative angles or zero to be written to the file (bitter experience suggests some old GSAS version didn't appreciate these points).

```

< gsasout 61 > ≡
  subroutine gsasout(is)
  use rebin ! ibin, bincen
  use specfiles
  use summation
  integer(kind=4),intent(in)::is
  integer(kind=4) :: low, ipts, nrec, i, k, n
  real(kind=8), dimension(2,5) :: vals
  real(kind=8) :: top ! for scaling to format
  character(len=10) :: str
  if(units.ne.'T')then
    write(*,*)'Can only write GSAS format for two theta data'
    stop
  endif
  call filext('.gsa',is)
  open(file=outfile,unit=ioutunit,status='UNKNOWN')
  write(title,1002) specsfilename(1:len_trim(specsfilename)), &
& specsdate(1:len_trim(specsdate))
1002 format(a,1X,a)
< scantitle 60b >
  write(ioutunit,'(a80)')title
  write(ioutunit,1015)
1015 format("Instrument parameter      id31.prm          ", &
& " ")
  ilow=getfirstpoint(hist,2,npts,2)
  ihigh=getlastpoint(hist,2,npts,2)
  low=ibin(0.0d0)
  if(low.gt.ilow) then ! the data goes below zero !
    low=low+1          ! one point after zero
  else
    low=ilow           ! value from module
  endif
  ipts=ihigh-low      ! how many points for GSAS
  nrec=(ipts-mod(ipts,5))/5 ! how many records
  if(mod(ipts,5).ne.0)nrec=nrec+1
  write(ioutunit,1000)ipts,nrec,bincen(low)*100.0d0,step*100.0d0
!      1234567 +16      123456      20      123456789012 = 61(+19=80)
1000 format('BANK 1 ',2(I7,1X),'CONST ',2(F9.3,1X),' 0.0 0.0 ESD' &
!      1234567890123456789
& ', ' ')
! scale factor to fit into F8.1 format? Max must be less than 999999
do i = low,ihigh
  if(hist(1,i).gt.top)top=hist(1,i)
enddo
if (top.gt.999999.9d0)then
  top=999999.9d0/top
  write(*,'(a)')'Rescaled your data for GSAS to avoid overflows'
  write(*,*)' The data were multiplied by',top
else
  top=1.0d0
endif
n = ipts-mod(ipts,5)
do i = low,n+low-1,5
  do k=1,5
    vals(1,k)=hist(1,k+i-1)*top
    vals(2,k)=hist(2,k+i-1)*top
    if (vals(2,k).lt.0.)then
      vals(1,k)=0.0
      vals(2,k)=999999.9
    endif
  enddo
  write(ioutunit,'(10F8.1)')(vals(1,k),vals(2,k),k=1,5)
enddo
if(mod(ipts,5).ne.0) then

```

Now for something with the extension spf.



```

< spfout 63 > ≡
subroutine spfout(is)
  use rebin
  use summation
  use specfiles
  integer(kind=4),intent(in)::is
  integer(kind=4)::low,ipts
  character(len=18) :: time
  character(len=10) :: str
  integer(kind=4) :: n, i, k
  character(len=3) :: months(12)
  data months /'JAN','FEB','MAR','APR','MAY','JUN',
&              'JUL','AUG','SEP','OCT','NOV','DEC'/
  if(units.ne.'T')then
    write(*,*)'Can only write SPF format for two theta data'
    stop
  endif
  call filext('.spf',is)
  open(file=outfile,unit=ioutunit,status='UNKNOWN')
  ilow=getfirstpoint(hist,2,npts,2)
  ihigh=getlastpoint(hist,2,npts,2)
  low=ibin(0.0d0)
  if(low.gt.ilow) then ! the data goes below zero !
    low=low+1          ! one point after zero
  else
    low=ilow           ! value from module
  endif
  ipts=ihigh-low       ! how many points
  n = ipts-mod(ipts,10)
  write(title,1002)specsdate(1:len_trim(specsdate)),
&  specsfilename(1:len_trim(specsfilename))
1002 format(a,1X,a)
< scantitle 60b >
  write(ioutunit,'(a80)')title
  write(ioutunit,2)ipts,bincen(low),bincen(ihigh),step
2  format("      1      1      1"/1i8,"      1",
& 3f10.3," 0.000000 0.000000")
!      CCYYMMDD      HHMMSS.SSS
  call date_and_time(date=time(1:8),time=time(9:18))
! 123456789012345678
! hh:mm:ss dd-mm-yy
  read(time(5:6),'(i2)')i
  time=time(9:10)///': '///time(11:12)///': '///time(13:14)///' '///
& time(7:8)///'- '///months(i)///'- '///time(3:4)
  write(ioutunit,3)time(1:18)
3  format(" 100000.0      0.0      ",1a18,
& "      0      0 SYNCHROTRON ID31 ESRF"/8("      0.000"))
  do i = low,n+low-1,10
    write(ioutunit,4)bincen(i),(nint(hist(1,k)),k=i,i+9)
    write(ioutunit,5)(hist(2,k),k=i,i+9)
4  format(1f8.3,10i7)
5  format(8x,10f7.2)
  enddo
  if(mod(ipts,10).ne.0)then
    write(ioutunit,4)bincen(n+low),(nint(hist(1,k)),k=n+low,ipts)
    write(ioutunit,5)(hist(2,k),k=n+low,ipts)
  endif
  write(ioutunit,6)
6  format('      0 HKL reflection(s)'/
& '      0 excluded region(s)')
  close(ioutunit)
  write(*,'(a)')'Wrote '//outfile(1:len_trim(outfile))
end subroutine spfout

```

◇

Winmprof likes to have files in something called pds format. Here's something to produce that format.

```

< pdsout 64a > ≡
    subroutine pdsout(is)
    use rebin
    use summation
    use specfiles
    integer(kind=4),intent(in)::is
    integer (kind=4) :: low, ipts, i, k, n
    character(len=10):: str
    if(units.ne.'T')then
        write(*,*)'Can only write PDS format for two theta data'
        stop
    endif
    call filext('.pds',is)
    open(file=outfile,unit=ioutunit,status='UNKNOWN')
    write(title,1002)specsdate(1:len_trim(specsdate)),
    & specsfilename(1:len_trim(specsfilename))
1002 format(a,1X,a)
< scantitle 60b >
    write(ioutunit,'(1a80)')title
    write(ioutunit,7)
    ilow=getfirstpoint(hist,2,npts,2)
    ihigh=getlastpoint(hist,2,npts,2)
    low=ibin(0.0d0)
    if(low.gt.ilow) then ! the data goes below zero !
        low=low+1      ! one point after zero
    else
        low=ilow      ! value from module
    endif
    ipts=ihigh-low
7   format("   Start      Step      End      Monitor")
    write(ioutunit,8)bincen(low),step,bincen(ihigh)
8   format(3f9.3,'      10000')
    n = ipts-mod(ipts,10)
    do 800 i = low,n+low-1,10
        write(ioutunit,9)(hist(2,k),k=i,i+9)
        write(ioutunit,10)(nint(hist(1,k)),k=i,i+9)
9   format(10f8.4)
10  format(10i8)
800 continue
    if(mod(ipts,10).ne.0)then
        write(ioutunit,9)(hist(2,k),k=n+1,ipts)
        write(ioutunit,10)(nint(hist(1,k)),k=n+1,ipts)
    endif
    write(ioutunit,11)
11  format('   -1000'/'   -10000')
    close(ioutunit)
    write(*, '(a)')'Wrote '//outfile(1:len_trim(outfile))
    return; end subroutine pdsout

```

Fragment referenced in 72.

## 8.2 Diagnostic plots

When there is some kind of problem with the binning we will need to write out a diagnostic file showing what each of the separate channels has recorded.

```

< outputdiagnostic 64b > ≡

```

```

subroutine outputdiagnostic(wd)
use rebin ! ascan
use summation
real(kind=8),intent(in)::wd
integer(kind=4) :: i, ihigh, ilow, j
real(kind=8) :: y, ey2, wd2
! write a file in plotmtv format with the 9 channels and sum
open(unit=ioutunit,file='diag.mtv',status='UNKNOWN',
& form='FORMATTED', access='SEQUENTIAL',iostat=i)
if(i.ne.0)then
write(*,'(a)')'error opening diagnostic file'
endif
write(ioutunit,'(a)')'$ DATA = CURVE2D'
write(ioutunit,'(a)')'% xlabel = "Two Theta"'
write(ioutunit,'(a)')'% ylabel = "Cts/Monitor"'
write(ioutunit,'(a)')'% toplabel= "Diagnostic plot"'
do j=1,nchannel
! j-1 to fix the zeroth channel being first
write(ioutunit,'(a,i1,a)')'% linelabel = " MA',j-1,' "'
write(ioutunit,'(a,i2)')'% linecolor = ',j
write(ioutunit,'(a)')'% linetype=1 markertype=0'
ilow=getfirstpoint(ascan,2*nchannel,npts,2*j)
ihigh=getlastpoint(ascan,2*nchannel,npts,2*j)
do i=ilow,ihigh
if(ascan(2*j,i).gt.1.0d0)
&
write(ioutunit,'(2F15.8)')bincen(i),
&
& ascan(2*j-1,i)/ascan(2*j,i)/mult(j)
enddo
write(ioutunit,*)
enddo
write(ioutunit,'(a,i2,a)')'% linelabel = "Total"'
write(ioutunit,'(a,i2)')'% linecolor = ',nchannel+1
write(ioutunit,'(a)')'% linetype=1 markertype=0'
ilow=getfirstpoint(hist,2,npts,2)
ihigh=getlastpoint(hist,2,npts,2)
do i=ilow,ihigh
write(ioutunit,'(2F15.8)')bincen(i),hist(1,i)
enddo
write(ioutunit,*)
write(ioutunit,'(a,i2,a)')'% linelabel = "outliers"'
write(ioutunit,'(a,i2)')'% markercolor = ',nchannel+2
write(ioutunit,'(a)')'% linetype=0 markertype=2'
do i=ilow,ihigh
if(ctchan(i).ge.2)then
do j=1,nchannel
if(ascan(2*j,i).gt.1.0d0)then
< getwd2 49>
if(wd2.gt.wd)write(ioutunit,'(2F15.8)')bincen(i),y
endif
enddo
endif
enddo
write(ioutunit,*)
write(ioutunit,'(a)')'$ END'
write(*,'(a,G12.5)')'Wrote diag.mtv file, outliers at ',sqrt(wd)
close(ioutunit)
return
end subroutine outputdiagnostic

```

◇

The w32 version of this program writes out the diagnostic plot file in a format suitable for the presto plotting program (a windows port of xmgr/grace).

```

< outputdiagnosticw32 66 > ≡
    subroutine outputdiagnosticw32(wd)
    use rebin ! ascan
    use summation
    real(kind=8),intent(in)::wd
    integer(kind=4) :: i, j, k
    real(kind=8) :: y, ey2, wd2
    ! write a file in plotmtv format with the 9 channels and sum
    open(unit=ioutunit,file='diag.mtv',status='UNKNOWN',
    & form='FORMATTED', access='SEQUENTIAL',iostat=i)
    if(i.ne.0)then
        write(*,'(a)') 'error opening diagnostic file'
        stop
    endif
    write(ioutunit,'(a)')'#@set[0].point.style none'
    write(ioutunit,'(a)')'#@set[1].point.style none'
    write(ioutunit,'(a)')'#@set[2].point.style none'
    write(ioutunit,'(a)')'#@set[3].point.style none'
    write(ioutunit,'(a)')'#@set[4].point.style none'
    write(ioutunit,'(a)')'#@set[5].point.style none'
    write(ioutunit,'(a)')'#@set[6].point.style none'
    write(ioutunit,'(a)')'#@set[7].point.style none'
    write(ioutunit,'(a)')'#@set[8].point.style none'
    write(ioutunit,'(a)')'#@set[9].point.style none'
    write(ioutunit,'(a)')'#@set[10].point.style circle'
    write(ioutunit,'(a)')'#@set[0].line.style solid'
    write(ioutunit,'(a)')'#@set[1].line.style solid'
    write(ioutunit,'(a)')'#@set[2].line.style solid'
    write(ioutunit,'(a)')'#@set[3].line.style solid'
    write(ioutunit,'(a)')'#@set[4].line.style solid'
    write(ioutunit,'(a)')'#@set[5].line.style solid'
    write(ioutunit,'(a)')'#@set[6].line.style solid'
    write(ioutunit,'(a)')'#@set[7].line.style solid'
    write(ioutunit,'(a)')'#@set[8].line.style solid'
    write(ioutunit,'(a)')'#@set[9].line.style solid'
    write(ioutunit,'(a)')'#@set[10].line.style none'
    write(ioutunit,'(a)')'#@set[0].line.color custom 255 0 0'
    write(ioutunit,'(a)')'#@set[1].line.color custom 215 30 0'
    write(ioutunit,'(a)')'#@set[2].line.color custom 180 60 0'
    write(ioutunit,'(a)')'#@set[3].line.color custom 150 90 0'
    write(ioutunit,'(a)')'#@set[4].line.color custom 120 120 0'
    write(ioutunit,'(a)')'#@set[5].line.color custom 90 150 0'
    write(ioutunit,'(a)')'#@set[6].line.color custom 60 180 0'
    write(ioutunit,'(a)')'#@set[7].line.color custom 30 215 0'
    write(ioutunit,'(a)')'#@set[8].line.color custom 0 25 0'
    write(ioutunit,'(a)')'#@set[9].line.color custom 0 0 0'
    write(ioutunit,'(a)')'#@set[10].point.color custom 0 0 0'
    do i=1,npts
    if(ctchan(i).lt.1)cycle
    write(ioutunit,'(11F15.8)',advance='no')bincen(i),
    & (ascan(2*j-1,i)/(1.0d0+ascan(2*j,i)*mult(j)), j=1,nchannel),
    & hist(1,i)
    k=0
    if(ctchan(i).ge.2)then
        do j=1,nchannel
            if(ascan(2*j,i).gt.1.0d0)then
                < getwd2 49 >
                if(wd2.gt.wd .and. k.eq.0)then

```

```

        write(ioutunit,'(F15.8)')
&        ascan(2*j-1,i)/(ascan(2*j,i)*mult(j))
        k=1
    endif
endif
enddo
endif
if(k.eq.0)write(ioutunit,'(F15.8)')-0.001d0
enddo
write(ioutunit,*)
write(*,'(a,G12.5)')'Wrote diag.mtv file, outliers at ',sqrt(wd)
close(ioutunit)
return; end subroutine outputdiagnosticw32
◇

```

Fragment referenced in 72.

```

⟨ dumpscan 67 ⟩ ≡
subroutine dumpscan(is)
use rebin
integer(kind=4) :: i, j
integer(kind=4),intent(in)::is
real(kind=8),dimension(9) :: sig
call filext('.dum',is)
open(unit=ioutunit,file=outfile,status='UNKNOWN',
& form='FORMATTED',access='SEQUENTIAL',iostat=i)
if(i.ne.0) stop 'error opening output file'
do i=1,npts
    do j=1,9
        if(ascan(2*j,i).gt.0)then
            sig(j)=ascan(2*j-1,i)/ascan(2*j,i)
        else
            sig(j)=0.0
        endif
    enddo
    write(ioutunit,'(10(F12.8,1X))')bincen(i),(sig(j),j=1,9)
enddo
close(ioutunit)
return
end subroutine dumpscan
◇

```

Fragment referenced in 72.

We should have an "undumpscan" routine which can take up where dumpscan has left off. It just needs to include some additional information in the dumped out format. This was only ever included for debugging and on line viewing of detector channels - so perhaps a bit pointless. FIXME - range of scan to dump.

### 8.3 epf and inp formats

```

< outputepf 68a > ≡
    subroutine epfout(n)
    use rebin
    use summation
    integer(kind=4) :: i, ilow, ihigh, n
    call filext('.epf',n) ! epfs never refer to a scan, always a sum
    open(unit=ioutunit,status='UNKNOWN',file=outfile)
    ilow=getfirstpoint(hist,2,npts,2)
    ihigh=getlastpoint(hist,2,npts,2)
    if(ilow.eq.-1 .or. ihigh.eq.-1)then
        ! No points with valid data found !
        return
    endif
    do i=ilow,ihigh
        write(ioutunit,*)bincen(i),hist(1,i),hist(2,i)
    enddo
    close(ioutunit)
    write(*,'(a)')'Wrote '//outfile(1:len_trim(outfile))
    return
end subroutine epfout

```

◇

Fragment referenced in 72.

```

< outputxye 68b > ≡
    subroutine xyeout(n)
    use rebin
    use summation
    integer(kind=4) :: i, ilow, ihigh, n
    character(len=4) :: extn
    if(units .eq. 'T') extn = '.xye'
    if(units .eq. 'Q') extn = '.qye'
    if(units .eq. 'R') extn = '.q2'
    call filext(extn,n) ! epfs never refer to a scan, always a sum
    open(unit=ioutunit,status='UNKNOWN',file=outfile)
    ilow=getfirstpoint(hist,2,npts,2)
    ihigh=getlastpoint(hist,2,npts,2)
    if(ilow.eq.-1 .or. ihigh.eq.-1)then
        ! No points with valid data found !
        return
    endif
    do i=ilow,ihigh
        write(ioutunit,'(F12.6,2(1X,G14.8))')bincen(i),hist(1,i),hist(2,i)
    enddo
    close(ioutunit)
    write(*,'(a)')'Wrote '//outfile(1:len_trim(outfile))
    return
end subroutine xyeout

```

◇

Fragment referenced in 72.

```

< outputinp 69a > ≡
      subroutine outputinp(n)
      use rebin
      use summation
      integer(kind=4),intent(in)::n
      integer(kind=4) :: i, ilow, ihigh
      character(len=15):: string
      write(string,'(i10)')n
      string=adjustl(string)
      if(units.eq.'T') string=string(1:len_trim(string))//'.inp'
      if(units.eq.'Q') string=string(1:len_trim(string))//'.inq'
      if(units.eq.'R') string=string(1:len_trim(string))//'.inq2'
!      write(*,*)'String was ',string
      open(unit=ioutunit,status='UNKNOWN',file=string)
      ilow=getfirstpoint(hist,2,npts,2)
      ihigh=getlastpoint(hist,2,npts,2)
      if(ilow.eq.-1 .or. ihigh.eq.-1)then
        ! No points with valid data found !
        return
      endif
      do i=ilow,ihigh
        write(ioutunit,*)bincen(i),hist(1,i),hist(2,i)
      enddo
      close(ioutunit)
      write(*,'(a)')'Wrote '//string(1:len_trim(string))
      return
    end subroutine outputinp

```

◇

Fragment referenced in 72.

## 8.4 Logfiles

When running through the binning program there is a fair bit of information which we might want to preserve - the kind of thing which belongs in a lab notebook. Much of this will flash past on the screen before anyone has a chance to see it, so the plan is to write things like motor positions, dates, temperatures and so on to a logfile for the particular binning run.

We will need (at least) three routines - one to start a log, another to write to it and finally a close at the end. The closing will be done by the general **tidyup** routine, so in fact we only need to provide the first two.

```

< openlogfile 69b > ≡
      subroutine openlogfile
      integer(kind=4) :: i
      call filext('.log',0) ! logfiles for all scans, not one each
      open(unit=logfile,file=outfile,status='UNKNOWN',
& form='FORMATTED',access='SEQUENTIAL',iostat=i)
      if(i.ne.0)stop 'could not open logfile'
      return; end subroutine openlogfile

```

◇

Fragment referenced in 72.

Since the outputfiles module uses the various other modules, this routine cannot be available to any routine which is contained in another module. Only main programs or the processscan stuff can access the log file. This is a sort of intentional design decision... the subroutines which are in modules aren't meant to be doing any user type io, even if they currently are.

```

< wlogfile 70a > ≡
  subroutine wlogfile(s)
    character(len=*) :: s
    write(logfile,'(a)')s
    return; end subroutine wlogfile

```

◇

Fragment referenced in 72.

## 8.5 Utils

```

< filext 70b > ≡
  subroutine filext(extn,is)
    use specfiles
    character(len=4),intent(in) :: extn
    integer(kind=4),intent(in) :: is
    integer(kind=4) :: n, i
    character(len=15):: string
    if(is.gt.0)then
      write(string,'(i10)')is ! scan number into string
      string=adjustl(string) ! shift to left end of string
    endif
    n=len_trim(filnam)
    if(filnam(n-3:n).eq.'.dat')n=n-4 ! overwrite .dat if exists
    if(is.gt.0) then
      write(outfile,'(a)')filnam(1:n)//'_ '//string(1:len_trim(string)) &
      & //extn(1:4)
    else
      write(outfile,'(a)')filnam(1:n)//extn(1:4)
    endif
    ! strip any / or \ from the start of the filename so that all output
    ! is in the working directory.
    n=len_trim(outfile)
    do i=n,1,-1
      if(outfile(i:i).eq.'\' .or. outfile(i:i).eq.'/')then
        outfile=outfile(i+1:n); exit; endif; enddo
    end subroutine filext

```

◇

Fragment referenced in 72.

```

< getfirstpoint 70c > ≡
  integer(kind=4) function getfirstpoint(array,m,n,l)
    integer(kind=4), intent(in) :: l, m, n
    real(kind=8), intent(in), dimension(m,n) :: array
    integer(kind=4) :: i
    getfirstpoint=1
    do i=1,n
      if(array(1,i).gt.0.0d0)then
        getfirstpoint=i
        return
      endif
    enddo
    return
  end function getfirstpoint

```

◇

Fragment referenced in 72.



```

< getlastpoint 71a > ≡
    integer(kind=4) function getlastpoint(array,m,n,l)
    integer(kind=4), intent(in) :: m,n,l
    real(kind=8), intent(in), dimension(m,n) :: array
    integer(kind=4) :: i
    getlastpoint=n
    do i=n,1,-1
        if(array(1,i).gt.0.0d0)then
            getlastpoint=i
            return
        endif
    enddo
    return
end function getlastpoint

```

Fragment referenced in 72.

In general the ID31sum program will produce an epf file and the ID31sumall program will produce a series of inp files. In case we want to output GSAS/PDS/SPF etc files from the binning program, instead of running something like epfout, then we need a general routine to decide which ones to output and actually write them. The default will be to only write the inp/epf files, but if a command line flag is specified then we will write out the other files as well.

```

< outputformats 71b > ≡
    subroutine outputformats(is)
    integer(kind=4),intent(in)::is
    if(epf) call epfout(is)
    if(xye) call xyeout(is)
    if(gsas) call gsasout(is)
    if(spf)call spfout(is)
    if(pds)call pdsout(is)
    return; end subroutine outputformats

```

Fragment referenced in 72.

## 8.6 Output files module

A module collecting together all the various output subroutines.

```

< outputfiles 72 > ≡
    module outputfiles
    !      use summation ! to get the final dataset
    < outputfilesvars 60a >
        logical :: bcm=.false.
        contains
        < outputformats 71b >
        < gsasout 61 >
        < spfout 63 >
        < pdsout 64a >
        < outpute pf 68a >
        < outputxye 68b >
        < outputinp 69a >
        < outputdiagnostic 64b >
        < outputdiagnosticw32 66 >
        < getfirstpoint 70c >
        < getlastpoint 71a >
        < dumpscan 67 >
        < filext 70b >
        < openlogfile 69b >
        < wlogfile 70a >
        < bcmfile 53 >
        < rstset 76 >
        < renset 77 >
    end module outputfiles

```

Fragment referenced in 40a, 57b, 86b, 89, 91.

## 9 Driver programs

These are to be the final user interface to the binning routines. Something to replace the scripts `binit` and `binem` is required and perhaps some new programs which do genuinely new things. All will initially need a routine for interpreting the command line. The usage is intended to be similar to the previous situation, so that the first argument is a filename, the second is the step size and the third and fourth are the first and last scans to process. The `glic` and `zinger` eliminations together with any other esoteric options will need flags to label them.

### 9.1 Specifying user options

These program might one day be used behind a graphical interface, or read input files to determine the user options. However, as an initial method they will be set up to be driven from the command line. The whole bundle of information needed will be supplied by a `useroptions` module, which can be replaced by something else implementing the same functions in the future.

Interpreting the command line should be relatively straightforward, we just use the `iargc` and `getarg` functions to determine the first four arguments and notify if any are missing while supplying a sensible defaults if we at least have a filename. Some extra arguments to replace the `resum` command are also needed. These will be `ed=n1,n2,n3` where `n1,n2,n3` are detectors to exclude and `es=m1,m2,m3` where `m1,m2,m3` are scans to exclude. When `glic` and `zinger` elimination are tackled they will also need some options.

This should perhaps be updated to use Lawson Wakefield's `f2kcli` module, as the command line is supposed to be getting standardised in the next few years.

```

⟨ getcmdline 73 ⟩ ≡
    subroutine getcmdline
    use rebin
    use specfiles
    integer(kind=4) :: iarg, i
    integer(kind=4),external :: iargc
    character(len=256) :: string ! massive string in case of silly user
    step=0.001d0; ifirstscan=0; ilastscan=0
    iarg=iargc()
    if(iarg.eq.0)then ; call helpmsg ; endif
    call getarg(1,filnam) ! get's filename
    if(iarg.eq.1) then
        write(*,'(a)')'No stepsize, assuming 0.001, and all scans'
        goto 100 ! return
    endif
    if(iarg.ge.2)then
        call getarg(2,string) ! get the stepsize
        read(string,*,err=10,end=10)step
        user_step = step
    endif
    if(iarg.ge.3)then
        call getarg(3,string) ! get the first scan to bin
        read(string,*,err=10,end=10)ifirstscan
        if(ifirstscan.eq.0)write(*,*)
    endif
    &'Zero for first scan is going to cause problems... sorry' ! FIXME?
    if(iarg.ge.4)then
        call getarg(4,string) ! get last scan to bin
        read(string,*,err=10,end=10)ilastscan
    endif
    if(iarg.ge.5)then
        do i=5,iarg
            call getarg(i,string); call option(string)
        enddo
    endif
    goto 100
10  write(*,*)'Problems interpreting your command line'
    call helpmsg
100  return
    end subroutine getcmdline

```

Fragment referenced in 87a.

Hopefully that will catch errors like "binit step file silly option".

For adding as many bells and whistles as we like after the basic options there is an option routine which takes a string and calls a routine to deal with that particular option. Not all options are going to make sense for all programs, nevertheless we'll process them anyway. Currently we have ed and es to specify excluded detectors and scans.

```

⟨ option 74 ⟩ ≡
  subroutine option(string)
  use outputfiles ! logicals for which files to write
  character(len=*) , intent(in) :: string
  if(string(1:1).ne.' ') then
    select case (string(1:3))
      case('mon') ; call setmonitorcol(string)
      case('wvl') ; call setwavelength(string)
      case('uni') ; call setunits(string)
      case('ed=') ; call exdet(string)
      case('es=') ; call exscan(string)
      case('is=') ; call incscan(string)
      case('ef=') ; call exfile(string)
      case('low') ; call lowtth(string)
      case('hig') ; call hightth(string)
      case('sca') ; call scale(string)
      case('ste') ; call minstepset(string)
      case('wd=') ; call wdset(string)
      case('alp') ; call alpset(string)
      case('ren') ; call renormset(string)
      case('mm=') ; call minmonset(string)
      case('mr=') ; call minrenormset(string)
      case('win') ; call windowset(string)
      case('zap') ; call zapset(string)
      case('sup') ; call superzapset(string)
      case('med') ; call medianchannelset(string)
      case('3pf') ; call filterset(string)
      case('bcm') ; call bcmset(string)
      case('snb') ; call snblset(string)
      case('rst') ; call rstset(string)
      case('rnd') ; call renset(string)
      case('nod') ; if(string(1:6).eq.'nodiag')then
        diag=.false. ; else ;
        write(*,*)'Sorry, I did not understand the command line'
        write(*,*)string(1:len_trim(string)) ; endif
      case('gsa') ; gsas=.true.
      case('spf') ; spf=.true.
      case('pds') ; pds=.true.
      case('epf') ; epf=.true.
      case default
        write(*,*)'Sorry, I did not understand the command line'
        write(*,*)string(1:len_trim(string))
        call helpmsg
      end select
    endif
  return; end subroutine option

```

◇

Fragment referenced in 87a.

To read the list of excluded detectors as a string **ed=1,2,3** would mean that detectors 1, 2 and 3 are to be excluded. **exdetlist** will be array provided by this module.

```

⟨ exdet 75 ⟩ ≡
      subroutine exdet(string)
      ! Read a comma separated list of detectors to ignore for the
      ! final sum
      use rebin
      character(len=*), intent(in) :: string
      integer(kind=4) :: i, j
      i=ncommas(string)+1
      allocate(exdetlist(i))
      read(string(4:len(string)),*,err=10)exdetlist
      do j=1,i
      ! Added a plus one here to go from channel zero
      logexdet(exdetlist(j)+1)=1
      write(*,'(a,i2)')'Intending to exclude channel ',exdetlist(j)
      enddo
      goto 100
10      STOP 'Could not understand your list of excluded detectors'
100     return; end subroutine exdet

```

Fragment referenced in 87a.

Add a jitter offset at the start of each scan to hide steps in the background when different channels have differing background. Philosophically this is not a good solution. Should actually measure the background for each channel independently and use that...

```

⟨rstset 76⟩ ≡
  subroutine rstset(string)
  use rebin
  character(len=*), intent(in) :: string
  integer itok, iprev, itot, ichan, i
  itot = len_trim(string)
  itok = index( string(5:itot), ",")
  if (itok.eq.0) then
    read(string(5:itot),*,err=10) randomstart
    do i = 1, nchan
      rstchan(i) = .true.
    enddo
    goto 100
  else
    read(string(5:5+itok-2),*,err=10) randomstart
    iprev = 5+itok
    do
      itok = index( string(iprev:itot), ",")
      if (itok.eq.0) then
        if(iprev.le.itot)then
          read( string(iprev:itot), *, err=10, end=10) ichan
          call rstadd(ichan)
        endif
        exit ! loop
      endif
      read( string(iprev:iprev+itok-1), *, err=10, end=10) ichan
      call rstadd(ichan)
      iprev = iprev+itok
    enddo
  endif

  goto 100
10  write(*,*)itok,iprev,itot,string
  STOP 'Could not understand your rst=x.xxx request'
100 userandomstart = .true.
  write(*,'(A,1X,f7.5,1X,A,$)') 'Applying random start of',randomstart, &
  & 'to channels:'
  do i = 1,nchan
    if (rstchan(i)) write(*,'(1X,I2,$)') i-1
  enddo
  write(*,*)
  return; end subroutine rstset

  subroutine rstadd(ichan)
  use rebin
  integer ichan
  if ((ichan .lt. 0).or. (ichan .ge. NCHAN)) then
    write(*,*) 'rst channel out of range',ichan
    stop
  endif
  ! the +1 makes it go from zero
  rstchan(ichan+1) = .true.
  return
  end subroutine rstadd

```

◇

Fragment referenced in 72.

Here is the same code copied and pasted for a random ending point on a scan to complement the random starting point. We will take the end point from the #S declaration at the beginning of the scan.

```

⟨ reset 77 ⟩ ≡
  subroutine reset(string)
  use rebin
  character(len=*), intent(in) :: string
  integer itok, iprev, itot, ichan, i
  itot = len_trim(string)
  itok = index( string(5:itot), ",")
  if (itok.eq.0) then
    read(string(5:itot),*,err=10) randomend
    do i = 1, nchan
      renchan(i) = .true.
    enddo
    goto 100
  else
    read(string(5:5+itok-2),*,err=10) randomend
    iprev = 5+itok
    do
      itok = index( string(iprev:itot), ",")
      if (itok.eq.0) then
        if(iprev.le.itot)then
          read( string(iprev:itot), *, err=10, end=10) ichan
          call renadd(ichan)
        endif
        exit ! loop
      endif
      read( string(iprev:iprev+itok-1), *, err=10, end=10) ichan
      call renadd(ichan)
      iprev = iprev+itok
    enddo
  endif

  goto 100
10  write(*,*)itok,iprev,itot,string
  STOP 'Could not understand your ren=x.xxx request'
100 userandomend = .true.
  write(*, '(A,1X,f7.5,1X,A,$)') 'Applying random end of',randomend, &
  & 'to channels:'
  do i = 1,nchan
    if (renchan(i)) write(*, '(1X,I2,$)') i-1
  enddo
  write(*,*)
  return; end subroutine reset

  subroutine renadd(ichan)
  use rebin
  integer ichan
  if ((ichan .lt. 0).or. (ichan .ge. NCHAN)) then
    write(*,*) 'ren channel out of range',ichan
    stop
  endif
  ! the +1 makes it go from zero
  renchan(ichan+1) = .true.
  return
  end subroutine renadd

```

◇

Fragment referenced in 72.

Unfortunately the instrument occasionally gives spurious bits of background only in some channels for small twotheta ranges. Since these have not yet been fixed there is a need to just throw away channels

for small regions. A subroutine to deal with that will be needed.

```

< exfile 78 > ≡
  subroutine exfile(string)
    ! Given a filename in string read in the file and set up for
    ! excluding regions of channels
    use rebin
    character(len=*), intent(in) :: string
    integer(kind=4)i,ier,ns
    real(kind=8)x1,xh
    open(unit=29,file=string(4:len_trim(string)),
    & status='OLD',iostat=ier)
    if(.not.(ier.eq.0))then
      write(*,'(a,i5)')'Error opening your file '//
    & string(4:len_trim(string))//' iostat=',ier
      write(*,'(a)')
    &
    &'want exclusion file, lines must have channel lowtth hightth scan'
    write(*,'(a)')'Make scan number less than zero for all'
      stop
    endif
    write(*,'(a)')
    &
    &'Reading exclusion file, lines must have channel lowtth hightth scan'
    write(*,'(a)')'Make scan number less than zero for all'
1    read(29,*,end=2,err=2)i,x1,xh,ns
    if((i.ge.0).and.(i.le.NCHAN)) then
      if(logexdet(i+1).eq.0) then
        logexdet(i+1)=2 ! set that this has excluded regions
      else
        write(*,'(a,i3,a)')'Channel',i,' already excluded'
      endif
    else
      write(*,'(a,i3,a)')'Error in your exfile, channel',i,
    &
    & 'not allowed'
      stop ! Operator error - give up
    endif
    iexrc=iexrc+1 ! increment number of excluded regions
    goto 1
2    allocate(exarray(iexrc,3)) ! holds low tth, high tth pairs
    allocate(iexarray(iexrc)) ! holds channel number
    rewind(29)
    do i=1,iexrc
      read(29,*)iexarray(i),exarray(i,1:2),exarray(i,3)
      write(*,'(a,i2,a,f10.6,a,f10.6)')'Excluding channel ',
    &
    & iexarray(i),' from ', exarray(i,1),' to ',exarray(i,2)
      if(exarray(i,3).lt.1) then
        write(*,*)"in all scans"
      else
        write(*,'(a,i4)')'in scan',int(exarray(i,3))
      endif
    enddo
    close(29)
    return
  end subroutine exfile

```

◇

Fragment referenced in 87a.

As for excluding detectors the same format is to be used for excluding scans. **exscanlist** comes from this module.



```

< exscan 79a > ≡
  subroutine exscan(string)
    ! Read a comma separated list of scans to skip over when summing
    ! Place these in an array which only nextscan knows about
    character(len=*),intent(in) :: string
    nexclد=ncommas(string)+1
    allocate(exscanlist(nexclد))
    read(string(4:len(string)),*)exscanlist
    return; end subroutine exscan

```

Fragment referenced in 87a.

Opposite of exscan - if you want to only include a certain list of scans.

```

< incscan 79b > ≡
  subroutine incscan(string)
    ! Read a comma separated list of scans to skip over when summing
    ! Place these in an array which only nextscan knows about
    character(len=*),intent(in) :: string
    integer(kind=4),allocatable :: temp(:)
    integer(kind=4) :: i, j, ninclد
    ninclد=ncommas(string)+1
    allocate(temp(ninclد))
    read(string(4:len(string)),*)temp
    ! total = ilastscan - ifirstscan + 1 ... eg 2 -> 10 is 2,3,4,5,6,7,8,9,10
    !                                     1 2 3 4 5 6 7 8 9
    nexclد=ilastscan - ifirstscan + 1 - ninclد
    allocate(exscanlist(nexclد))
    i=1

    do j=ifirstscan,ilastscan
      if (.not. inlist(j,temp,ninclد)) then
        exscanlist(i)=j
        i=i+1
      endif
    enddo
    if (i-1 .ne. nexclد) then
      write(*,*)' problem - debug is= please'
    endif
    deallocate(temp)
    return; end subroutine incscan

```

Fragment referenced in 87a.

Allowing the user to specify high and low two theta limits for the binning program means we can save on memory and override the default values if it is ever necessary. Ideally the program should figure these limits out for itself from the SPEC file, but I still don't see any way to do this without reading the entire file before doing any binning, so it is down to the user.

```

⟨ limits 80a ⟩ ≡
    subroutine lowtth(s)
    use rebin
    character(len=*),intent(in) :: s
    if(s(1:7).eq.'lowtth=')read(s(8:len_trim(s)),*)tthlow
    user_tthlow = tthlow
    return ; end subroutine lowtth
    subroutine hightth(s)
    use rebin
    character(len=*),intent(in) :: s
    if(s(1:8).eq.'hightth=')read(s(9:len_trim(s)),*)tthhigh
    user_tthhigh = tthhigh
    return; end subroutine hightth

```

Fragment referenced in 87a.

For the diagnostic plots, outliers are placed at  $6\sigma$  by default. If the user prefers a different criterion, the cutoff can be adjusted by specifying **wd=xx.xx** where **xx.xx** is the number of esd's a data point must differ by to be considered an outlier.

```

⟨ wdset 80b ⟩ ≡
    subroutine wdset(s)
    character(len=*),intent(in)::s
    if(s(1:3).eq.'wd=')read(s(4:len_trim(s)),*)wd
    wd=wd*wd
    return ; end subroutine wdset

```

Fragment referenced in 87a.

The minimum stepsize was hard wired into the program for a default for the expected resolution for the ID31 beamline. In the unlikely event that someone genuinely needs a finer step size they can force the program to accept the small number by specifying **step=xxx** on the command line, where step is their stepsize. This comes in handy for determining detector offsets, where very fine binning might be required with a reduced angular range.

```

⟨ minstepset 80c ⟩ ≡
    subroutine minstepset(s)
    use rebin
    character(len=*),intent(in):: s
    if(s(1:5).eq.'step=')read(s(6:len_trim(s)),*)aminstep
    user_step = aminstep
    return; end subroutine minstepset

```

Fragment referenced in 87a.

For deciding whether or not a particular  $2\theta$  value has been observed or not we insist on there having been a certain threshold number of counts arriving on the monitor spectrum at that angle. The variable **minmon** in the **rebin** module decides that threshold value, although the user may modify it by supplying an argument **mm=xx** where **xx** is the threshold value.

```

⟨ minmonset 80d ⟩ ≡
    subroutine minmonset(s)
    use rebin
    character(len=*),intent(in) :: s
    if(s(1:3).eq.'mm=')read(s(4:len_trim(s)),*)minmon
    return; end subroutine minmonset

```

Fragment referenced in 87a.

```

< minrenormset 81a > ≡
    subroutine minrenormset(s)
    use rebin
    character(len=*),intent(in) :: s
    if(s(1:3).eq.'mr=')read(s(4:len_trim(s)),*)minrenormsig
    write(*,1)"Only using where all channels have more than ",minrenormsig, &
    & " counts for renorm"
1    format(a,f8.2,a)
    return; end subroutine minrenormset

```

Fragment referenced in 87a.

```

< windowset 81b > ≡
    subroutine windowset(s)
    use rebin
    character(len=*),intent(in) :: s
    character(len=256) :: mys
    integer i
    if(s(1:7).eq.'window=')then
        do i=8,len_trim(s)
            if(s(i:i).eq.',')then
                mys(i:i)=' '
            else
                mys(i:i)=s(i:i)
            endif
        enddo
        read(mys(8:len_trim(s)),*)wincnt,winlow,winhigh
    endif
    write(*,*)'windowing counter ',wincnt(1:len_trim(wincnt)), ' from ', &
    & winlow, ' to ', winhigh
    winhighread=winhigh
    winlowread=winlow
    winlog=.true.
    return
end subroutine windowset

```

Fragment referenced in 87a.

## 9.2 Get monitor col

Normally the monitor column is called “Monitor”, but you might want something else.

```

< setmonitorcol 81c > ≡

    subroutine setmonitorcol(s)
    use specfiles
    character(len=*), intent(in) :: s
    if(s(1:4).eq."mon=")then
        MONITORCOL=s(5:len_trim(s))
        write(*,*)"Using monitor counts from column ",MONITORCOL
    else
        write(*,*)"I do not understand your argument",s(1:len_trim(s))
        write(*,*)"Perhaps you mean mon=Monitor ?"
    endif
end subroutine setmonitorcol

```

◇

Fragment referenced in 87a.

The default determination of the error bars includes a parameter to avoid problems with zero counts. Essentially the problem is that the square root of zero counts is a very poor approximation to the uncertainty on the number zero, so we take the square root of counts plus some parameter, called  $\alpha$ . By default the parameter is 0.5, but the user can choose some other number if they prefer, by adding **alp=xx.xx** to their command line.

```
< alpset 82a > ≡
  subroutine alpset(s)
  use summation
  character(len=*), intent(in) :: s
  if(s(1:4).eq.'alp=')read(s(5:len_trim(s)),*)alp
  return; end subroutine alpset
```

Fragment referenced in 87a.

Zap any outlying data point with the zapping option!

```
< zapping 82b > ≡
  subroutine zapset(s)
  use summation
  character(len=*), intent(in) :: s
  if(s(1:4).eq.'zap=')then
    read(s(5:len_trim(s)),*)zap
    zapping=.true.
    write(*,*)"Aiming to zap points which off by ",zap," sigma"
  else
    write(*,*)"did not understand your option ",s
  endif
  return; end subroutine zapset
```

Fragment referenced in 87a.

Try to do a better job on the zinger elimination by "cheating" when we read in lines of data from the datafile. We will take the median of three consecutive points, discarding the first and last points in the file, in the hope that this eliminates electrical spikes.

```
< filter 82c > ≡
  subroutine filterset(s)
  use pointfilter
  character(len=*), intent(in) :: s
  if(s(1:3).eq.'3pf')filterlogical=.true.
  write(*,*)"Aiming to apply a 3 point median filter"
  return; end subroutine filterset
```

Fragment referenced in 87a.

Yet another approach to zapping. The idea is that when there are at least 3 channels overlapping, we can calculate the mean and variance of the 3 independent signals and compare it to the error bar derived from the counting statistics. If the channels "agree" then all will be fine, but if the channels are markedly different we have a problem. We eliminate the highest value (it is always too much that is a problem) and recalculate. There must always be 2 channels left. The "improvement" after deleting a channel is measured in terms of the percentage reduction in the n channel esd versus the n-1 channel esd, where the esds are normalised to the "real esd" in each case.

```

⟨ superzapset 83a ⟩ ≡
  subroutine superzapset(s)
  use summation
  character(len=*),intent(in) :: s
  if(s(1:9).eq.'superzap=')then
    read(s(10:len_trim(s)),*,err=1)superzaplevel
    superzap=.true.
    write(*,'(A,F8.5)') "Superzapping at level ",superzaplevel
  else
    write(*,*) 'Do you mean to say "superzap=0.9" ? '
    write(*,*) 'So say it then, how else can I read the number!!'
    goto 1
  endif
  return
1 write(*,*) 'Could not figure out what you meant by'
  write(*,*) s(1:len_trim(s))
end subroutine superzapset

```

Fragment referenced in 87a.

```

⟨ medianchannelset 83b ⟩ ≡
  subroutine medianchannelset(s)
  use summation
  character(len=*),intent(in) :: s
  if(s(1:len_trim(s)).eq.'medianofchannels')then
    medianofchannels=.TRUE.
    write(*,'(A)') "Signal and statistics reflect the median", &
    & " channel * number_of_active_channels"
  else
    goto 1
  endif
  return
1 write(*,*) 'Could not figure out what you meant by'
  write(*,*) s(1:len_trim(s))
  write(*,*) 'Nearest guess is medianofchannels ??'
  write(*,*) 'IGNORING:',s(1:len_trim(s))
end subroutine medianchannelset

```

Fragment referenced in 87a.

By default the program takes detector efficiencies from the temp.res files, however it makes sense to be able to update the temp.res file if you decide you now have better data (particularly if the offsets are changing frequently). Adding the flag **renorm** to the command line should make the programs reset the detector efficiencies to the new set of numbers (when it has such a set of numbers to use).

```

⟨ renormset 83c ⟩ ≡
  subroutine renormset(s)
  use summation
  character(len=*)::s
  if(s(1:6).eq.'renorm') renorm=.true.
  return; end subroutine renormset

```

Fragment referenced in 87a.

A flag to indicate that we want a binned counts and monitor counts file.

```

< bcmset 84a > ≡
  subroutine bcmset(s)
  use outputfiles
  character(len=*)::s
  if(s(1:6).eq.'bcm') bcm=.true.
  return; end subroutine bcmset

```

Fragment referenced in 87a.

Flag the swiss norwegian beamline data format (nchan=6, with different names).

```

< snblset 84b > ≡
  subroutine snblset(s)
  use rebin
  character(len=*)::s
  if(s(1:4).eq.'snbl') snbl=.true.
  nchannel=6
  logexdet(7)=1
  logexdet(8)=1
  logexdet(9)=1
  write(*,'(a,i2)') "Exclude last 3 channels, as they don't exist"
  return; end subroutine snblset

```

Fragment referenced in 87a.

Determining what the user wants for the overall scale factor allows for some subjectivity. We will provide a routine to allow them to specify whether they want the highest peak to be the right height or the number of counts in the pattern to be conserved.

```

< scale 84c > ≡
  subroutine scale(string)
  use summation
  character(len=*),intent(in) :: string
  if(string(1:6).eq.'scalpk' )then ; isc=1 ; return ;endif
  if(string(1:7).eq.'scaltot')then ; isc=2 ; return ;endif
  ! Replace scalinp string with scalmon string
  if(string(1:8).eq.'scalmon=')then
    read(string(9:len_trim(string)),*)scalinp;isc=0
    write(*,'(a,E12.6,a)') 'Units will be counts per ',scalinp,
    & ' monitor ct' ; return ; endif
  if(string(1:7).eq.'scalmon')then ; isc=0 ; return ;endif
  write(*,'(a)') 'Didn't understand your optional argument:'
  write(*,'(a)') string(1:len_trim(string))
  return
  end subroutine scale

```

Fragment referenced in 87a.

Determining which scan the user wants to process next (dealing with excluded scans) is handled by a routine which uses a pair of variables to hold onto the first and last scan and returns with the next scan to process or -1 for all done.

```

< nextscan 85a > ≡
      integer(kind=4) function nextscan()
      ! Returns to next scan to process or -1 for all done
1      if(ifirstscan.eq.0 .and. ilastscan.eq.0) then           ! all scans
          icurrscan=icurrscan+1
          nextscan=icurrscan
          return ! exit here !
      elseif (ilastscan.eq.0) then                             ! one scan only
          if(icurrscan.eq.ifirstscan)then
              icurrscan=-1 ! stop
          else
              icurrscan=ifirstscan
              nextscan=ifirstscan
              return ! exit here !
          endif
      else
          if(icurrscan.eq.0)then                                ! needs to deal with ifirstscan=0
              icurrscan=ifirstscan
          else
              icurrscan=icurrscan+1
          endif
      endif
      if(allocated(exscanlist))then ! check if this is an excluded scan
          if(inlist(icurrscan,exscanlist,nexcld)) goto 1
      endif
      if(icurrscan.le.ilastscan)then
          nextscan=icurrscan
      else
          nextscan=-1
      endif
      return
      end function nextscan

```

Fragment referenced in 87a.

Is the integer *i* in the array *iarr* which has length *n*? A simple routine which might have some from a standard library if we'd found it.

```

< inlist 85b > ≡
      logical function inlist(i,iarr,n)
      integer(kind=4),intent(in) :: i, n
      integer(kind=4),dimension(n),intent(in)::iarr
      integer(kind=4) :: j
      inlist=.FALSE.
      do j=1,n
          if(iarr(j).eq.i) inlist=.TRUE.
      enddo
      return; end function inlist

```

Fragment referenced in 87a.

Counts the number of commas in a string - should perhaps generalise this for any character. Used for deciphering flags like **ed=2,5,17**

```

< ncommas 86a > ≡
    integer(kind=4) function ncommas(string)
    character(len=*),intent(in) :: string
    integer(kind=4) :: i
    ncommas=0
    do i=1,len(string)
        if(string(i:i).eq.',')ncommas=ncommas+1
    enddo
    return; end function ncommas

```

◇

Fragment referenced in 87a.

A quick test program to check that nextscan does what we intend:

```

"testnextscan.f90" 86b≡
    < specfiles 11a >
    < rebin 28 >
    < summation 56 >
    < outputfiles 72 >
    < pointfilter 41 >
    < useroptions 87a >
    < helpmsg 87b >
1000 format(' example: ',a,' testnextscan file.dat 0.01 1 10 es=5,8 ')
    STOP
    end subroutine helpmsg
    program testnextscan
    use specfiles
    use rebin
    use useroptions
    use summation
    integer(kind=4)::i ; i=0
    call getcmdline
    if(ifirstscan.eq.0 .and. ilastscan.eq.0)then
        call helpmsg
    endif
1    i=nextscan() ;
    if(i.gt.0)then
        write(*,'(1x,i5,$)')i
        goto 1
    endif
    end program testnextscan

```

◇

The options routines are pulled together in the module user options. When information needs to be passed from here into the individual subroutines it should be done by useroptions using that module and setting variables in it, not the other way around (so none of the other modules should depend on the useroptions stuff).



```

< useroptions 87a> ≡
    module useroptions
    integer(kind=4) :: ifirstscan, ilastscan, icurrscan=0,nexcl
    integer(kind=4),allocatable :: exscanlist(:),exdetlist(:)
    real(kind=8)::wd=36.0d0
    logical :: snbl=.false.
    contains
    < getcmdline 73>
    < option 74>
    < exdet 75>
    < exscan 79a>
    < incscan 79b>
    < exfile 78>
    < ncommas 86a>
    < limits 80a>
    < renormset 83c>
    < wdset 80b>
    < minstepset 80c>
    < minmonset 80d>
    < minrenormset 81a>
    < windowset 81b>
    < scale 84c>
    < nextscan 85a>
    < inlist 85b>
    < alpset 82a>
    < zapping 82b>
    < superzapset 83a>
    < medianchannelset 83b>
    < filter 82c>
    < bcmset 84a>
    < snblset 84b>
    < setmonitorcol 81c>
    < unitoptions 29>
    end module useroptions

```

◇

Fragment referenced in 40a, 57b, 86b, 89, 91.

The **isc** variable holds a number which tells us which scale factor to apply to any data which is written out. Binit style programs will want 1 or 2 and binem style will want zero.

### 9.3 helpmsg

A simple routine to write out an error message and then stop the program. Intended to be used for helping out the user. Needs to know the name of the calling program to give an example of usage, in the unlikely event that someone decides to rename their executable. Should give version information and should be writing something specific for a particular program. The idea is to make this a generic thing for copying and pasting around as required, with each program putting it's own version in. This is not included in the useroptions module to ake it possible to have a specific message for each program.

```

< helpmsg 87b> ≡
    subroutine helpmsg
    ! Writes a helpful message to stdout
    character(len=80)::name
    call getarg(0,name)
    write(*,'(a)')name(1:len_trim(name))// ' version 10-05-2011'
    write(*,*)
    write(*,1000)name(1:len_trim(name))

```

◇

Fragment referenced in 40a, 57b, 86b, 89, 91.

## 9.4 Something like binit

A program to replace the binit script. We ought to think carefully about the naming of these programs so as not to clash with any existing executables and provide something short and logical to type. `id31sum` appeals for now. The program starts out much as the `bindump` program from the binning section but proceeds to include the routines from the summing section as well. I'd like this program to just be a series of simple routine calls, with everything else wrapped up inside previously defined code.

```

" id31sum.f90" 89≡
  < specfiles 11a >
  < rebin 28 >
  < summation 56 >
  < outputfiles 72 >
  < processscan 33 >
  < useroptions 87a >
  < tidyup 98 >
  < helpmsg 87b >
  < id31summsg 90 >
  program id31sum
  use specfiles
  use rebin
  use useroptions
  use summation
  use outputfiles
  integer(kind=4)::i
  real start_time, end_time
  call cpu_time(start_time)
  isc=1 ! scalpk
  call getcmdline ! Get users options
! fill out names of monitor, tth, first and last columns
  if(snb1)then
    MONITORCOL="Mon";FIRSTDET="Det1";LASTDET="Det6";TWOTTH="TwoTheta"
  else
    FIRSTDET="MA0";LASTDET="MA8";TWOTTH="2_theta"
  endif
  call initialiserebin ! Allocate space and check parameters
  call getfile ! Opens the spec file
  call openlogfile ! Get the file for recording T etc open
  write(*,'(a,$)')'Processing scan '
1  i=nextscan() ! Gets next scan to be processed
  if(i.gt.0) then
    call processscan(i) ! sums into ascan array
    if(ispecerr.eq.-1)then
      ispecerr=0
      goto 1 ! next only if current was OK
    endif
  endif
  write(*,*) ! after non advancing io list of scans binned
  call calibsum
  if(bcm)call bcmfile(0)
  !if(diag) call outputdiagnosticw32(wd) ! for windows and prestoplot
  if(diag)call outputdiagnostic(wd)!If problems in combining chans
  hist=hist*scalinp ! apply scale factor after diagnostic plot
  select case(isc) ! rescale if necessary
  case(0) ; continue ! do nothing
  case(1) ; call scalpk ! scale to peak height
  case(2) ; call scaltot ! scale to total counts
  case default ; continue ! do nothing
  end select
! call outputxye ! Writes out sumdata array
  call outputformats(0) ! Any other (GSAS etc formats)
  call tidyup ! Frees allocated memory
  call cpu_time(end_time)
  write(*,'(a,f6.2,a)')'Time taken was ',end_time-start_time,'s'
  end program id31sum

```

◇

```

< id31summsg 90) ≡
1000 format('example: ',a,' file.dat 0.01 1 10' /                                &
&'will process the scans 1 to 10 from file.dat with binsize 0.01' &
&'Optional arguments:' &
&' ed=n1,n2 to exclude detectors n1 and n2 (default=none)' &
&' es=m1,m2 to exclude scans m1 and m2 (default=none)' &
&' is=m1,m2 to include only scans m1 and m2 (default=none)' &
&' lowtth=xx.xx to set min two theta to use (default=-30.0)' &
&' hightth=xx.xx to set max two theta to use (default=160.0)' &
&' step=x.xxxx to force a stepsize (if very small)' &
&' wd=xx to set a limit for outliers on diagnostic file', &
&' , diag.mtv' &
&' wvln=xx.xx to set the wavelength used for unit conversion', &
&' units=Q bins into  $Q = 4\pi\sin(\theta)/\text{wavelength}$ ', &
&' units=Q2 bins into  $Q^2 = 16\pi^2\sin^2(\theta)/\text{wavelength}^2$ ', &
&' alp=xx for  $\text{esd}=\sqrt{\text{cts}+\text{alp}}$ , default is 0.5' &
&' mm=xx sets the minimum monitor counts threshold (default=5)' &
&' mr=xx for minimum counts needed in all channels to use ' &
&' points for determining detector efficiencies (default=1)' &
&' zap=xx for esd level in filtering operation ' &
&' superzap=xx for zinger elimination,  $0.0 < \text{xx} < 1.0$  ' &
&' where xx is the fractional reduction in esd for ' &
&' removing the highest point' &
&' medianofchannels to get median_channel*n_active_channels ' &
&' 3pf for a 3 point median filter (you need to be desperate) ' &
&' rst=x.xx to randomly offset scan start points by rnd*x.xx ' &
&' rst=x.xx,1,2,3 to do that for channels 1,2,3 only ' &
&' rnd=x.xx to randomly offset scan end points by rnd*x.xx ' &
&' rnd=x.xx,1,2,3 to do that for channels 1,2,3 only ' &
&' renorm to use current effics instead of values from temp.res' &
&' nodiag to prevent creation of diagnostic file, diag.mtv' &
&' gsa to output a .gsa file for gsa's' &
&' spf to output a .spf file for profil' &
&' pds to output a .pds file for a pds file' &
&' epf to output a .epf file for a epf file' &
&' bcm to output a .bcm file (binned counts,monitor)' &
&' window=counter,low,high to reject bad points' &
&' eg: window=lake,4,6 for temperatures between 4 and 6 K' &
&' snbl if your data are in the Swiss Norwegian format' &
&' *** PATTERN SCALING ****' &
&' scalpk to scale highest peak to correct number of obs counts' &
&' scaltot to make total number of counts in scan correct ' &
&' scalmon=xxx for counts per xxx monitor counts (default=100000)' &
&' ****' &
stop
end subroutine helpmsg ◇

```

Fragment referenced in 89.

The program in the next section ought to be pretty much identical except it will sumchannels and outputscan within the processscan loop. Should also reinitialise things to zero.

## 9.5 Something like binem

Have a look at rocketeer <sup>3</sup> for plotting 2D datasets.

Essentially an identical program to id31sum, it just writes out and resets the summing arrays to zero between scans.

<sup>3</sup>[http://www.csar.uiuc.edu/F\\_software/rocketeer/Rocketeer\\_Users\\_Guide.htm](http://www.csar.uiuc.edu/F_software/rocketeer/Rocketeer_Users_Guide.htm)

```

"i31sumall.f90" 91≡
  < specfiles 11a >
  < rebin 28 >
  < summation 56 >
  < outputfiles 72 >
  < processscan 33 >
  < useroptions 87a >
  < tidyup 98 >
  < zerodata 92b >
  < helpmsg 87b >
  < i31sumallmsg 92a >
    program i31sumall
    use specfiles
    use rebin
    use useroptions
    use summation
    use outputfiles
    integer(kind=4)::i
    real start_time, end_time
    call cpu_time(start_time)
    call getcmdline                                ! Get users options
! fill out names of monitor, tth, first and last columns
    if(snb1)then
      MONITORCOL="Mon";FIRSTDET="Det1";LASTDET="Det6";TWOTTH="TwoTheta"
    else
      FIRSTDET="MA0";LASTDET="MA8";TWOTTH="2_theta"
    endif
    call initialiserebin                          ! Allocate space and check parameters
    call getfile                                  ! Opens the spec file
    call openlogfile
    write(*, '(a,$)') 'Processing scan '
1    i=nextscan()                                ! Gets next scan to be processed
    if(i.gt.0) then
      call processscan(i)                        ! sums into ascan array
      if(ispecerr.eq.-1)then
        ispecerr=0
        write(*,*)                               ! after non advancing io list of scans binned
        call calibsum                             ! Fills hist array from ascan array
        if(bcm)call bcmfile(i)
        hist=hist*scalinp                         ! scales to scalinp
        call outputinp(i)                         ! Writes out hist array
        call outputformats(i)
        call zerodata                             ! zero the data
        goto 1                                    ! next only if current was OK
      endif
    endif
    call tidyup                                  ! Frees allocated memory
    call cpu_time(end_time)
    write(*, '(a,f6.2,a)') 'Time taken was ',end_time-start_time,'/s'
    end program i31sumall

```

◇

```

< id31sumallmsg 92a > ≡
1000 format('example: ',a,' file.dat 0.01 1 10' /                                &
&'will process the scans 1 to 10 from file.dat with binsize 0.01' &
&'Optional arguments:' &
&' ed=n1,n2 to exclude detectors n1 and n2 (default=none)' &
&' es=m1,m2 to exclude scans m1 and m2 (default=none)' &
&' is=m1,m2 to include only scans m1 and m2 (default=none)' &
&' lowtth=xx.xx to set min two theta to use (default=-30.0)' &
&' hightth=xx.xx to set max two theta to use (default=160.0)' &
&' step=x.xxxx to force a stepsize (if very small)' &
&' wvln=xx.xx to set the wavelength used for unit conversion', &
&' units=Q bins into Q = 4*pi*sin(theta)/wavelength', &
&' units=Q2 bins into Q^2 = 16*pi^2*sin^2(theta)/wavelength^2', &
&' alp=xx for esd=sqrt(cts+alp), default is 0.5' &
&' zap=xx for esd level in filtering operation ' &
&' superzap=xx for zinger elimination, 0.0 < xx < 1.0 ' &
&' where xx is the fractional reduction is esd for ' &
&' removing the highest point' &
&' medianofchannels to get median_channel*n_active_channels ' &
&' 3pf for a 3 point median filter (you need to be desperate) ' &
&' rst=x.xx to randomly offset scan start points by rnd*x.xx ' &
&' window=counter,low,high to reject bad points' &
&' eg: window=lake,4,6 for temperatures between 4 and 6 K' &
&' mm=xx for to set minimum monitor counts threshold', &
&' (default=5)' &
&' mr=xx for minimum counts needed in all channels to use ' &
&' points for determining detector efficiencies (default=1)' &
&' snbl if your data are in the Swiss Norwegian format' &
&' scalmon=xxx to rescale .inp files (default=10000)' &
&' gsa to output a .gsa files' &
&' spf to output a .spf files' &
&' pds to output a .pds files' &
&' epf to output a .epf files' &
&' bcm to output a .bcm file (binned counts,monitor)')
stop; end subroutine helpmsg ◇

```

Fragment referenced in 91.

The zerodata subroutine just sets the various arrays which are filled in during binning to zero.

```

< zerodata 92b > ≡
subroutine zerodata
use rebin
use summation
ascan=0.0d0
sumdata=0.0d0
hist=0.0d0
! sums of various counts should also be zeroed, if they are used.
return
end subroutine zerodata ◇

```

Fragment referenced in 91.

## 9.6 Checking scans for consistency

Might have been better to do this within the binning programs, but a separate shell script is to be used. First it runs id31sum with scalinp=1.0 so we have counts per monitor count, then id31sumall with the same scalinp. Then it uses the little program id31check, which compares a series of .inp files to a .epf file.

The shell script is here:

```
"id31sumcheck" 93a≡
  echo EXECUTING: id31sum $* scalinp=1.0
  id31sum $* scalinp=1.0
  echo EXECUTING: id31sumall $* scalinp=1.0
  id31sumall $* scalinp=1.0
  echo "id31check " > id31cmd
  ls -rt *.epf | tail -1 >> id31cmd
  testnextscan $* | grep -v "Intending" >> id31cmd
  perl -e 'while(<>){chomp;print}' < id31cmd > id31cmd2
  echo >> id31cmd2
  echo EXECUTING:
  cat id31cmd2
  source id31cmd2
  /bin/rm id31cmd id31cmd2
  ◇
```

The fortran for the checking program is straightforward. Here it is:

```
"id31check.f90" 93b≡
  program id31check
  implicit none
  character(len=256)::line
  integer :: ier, iarg, ioutunit, jjj
  integer, external :: iargc
  integer(kind=4) :: i, ndata, iscan, n, j, k, n3s, n6s
  real(kind=8) :: x,y,esd, diff, c2, wd, tthlow, tthhigh, step
  real(kind=8),allocatable,dimension(:,:) :: data
  integer(kind=4),external :: jbin
  ! Read epf file in
  i=0; ier=0
  call getarg(1,line) ! name of epf file
  open(unit=20,file=line(1:len_trim(line)),status='OLD',iostat=ier)
  if(ier.ne.0)then
    write(*,'(a)')'Error opening '//line(1:len_trim(line))
    stop
  endif
  open(unit=21,status='SCRATCH',form='UNFORMATTED')
1  read(20,*,end=10,err=11)x,y,esd
  write(21)x,y,esd
  i=i+1
  goto 1
11 stop 'Data format problem'
10 allocate(data(3,i),stat=ier)
  if(ier.ne.0)stop 'Memory allocation problem'
  ndata=i
  rewind(21)
  do i=1,ndata
    read(21)data(1,i),data(2,i),data(3,i)
  enddo
  close(21)
  close(20)
  jjj=1; ioutunit=19
  ! write a file in plotmtv format with the 9 channels and sum
  open(unit=ioutunit,file='scans.mtv',status='UNKNOWN',
  & form='FORMATTED', access='SEQUENTIAL',iostat=ier)
  if(ier.ne.0)stop 'error opening diagnostic file'
  write(ioutunit,'(a)')'$ DATA = CURVE2D'
  write(ioutunit,'(a)')'% xlabel = "Two Theta"'
  write(ioutunit,'(a)')'% ylabel = "Cts/Monitor"'
```

```

write(ioutunit,'(a)')'% toplevel= "Scans and total plot"'
write(ioutunit,'(a,i2,a)')'% linelabel = "Total" '
write(ioutunit,'(a,i2)')'% linecolor = ',jjj
write(ioutunit,'(a)')'% linetype=1 markertype=0'
do j=1,ndata
  write(ioutunit,'(2F15.8)')data(1,j),data(2,j)
enddo
write(ioutunit,*)
! Epf file is now in array data
step=(data(1,ndata)-data(1,1))/real(ndata-1,8)
tthlow=data(1,1)-step/2.0d0
tthhigh=data(1,ndata)+step/2.0d0
write(*,*)'Opened summed file ',line(1:len_trim(line)), ' npts ', &
& ndata
write(*,1000)'Step ',step,' tthlow ',tthlow,' tthhigh ',tthhigh
1000 format(3(a,F15.8))
! Now read in the series of scans
iarg=iargc() ! total number of cmdline arguments
if(iarg.le.1)then
  stop 'need to supply a list of scans to check'
endif
!
1234567890123456789012345678901234567890
write(*,999)
999 format('      Scan      npts  3-sigma 6-sigma  chi**2')
do i=2,iarg
  call getarg(i,line)
  line=adjustl(line)
  n=len_trim(line)
  write(line,'(a)')line(1:n)//'.inp'
  open(unit=20,file=line(1:len_trim(line)),status='OLD',iostat=ier)
  if(ier.ne.0)then
    write(*,'(a)') 'Error opening '//line(1:len_trim(line))
    stop
  endif
!
  write(*,*)'Opened '//line(1:len_trim(line))
  write(ioutunit,*)
  write(ioutunit,'(a,a,a)')'% linelabel = ', &
& line(1:len_trim(line)),''
  jjj=jjj+1
  write(ioutunit,'(a,i3)')'% linecolor = ',jjj
  write(ioutunit,'(a)')'% linetype=1 markertype=0'
  c2=0.0d0; n3s=0; n6s=0; n=0
! read and check the .inp against the .epf now
2 read(20,*,end=21,err=22)x,y,esd
  write(ioutunit,'(2F15.8)')x,y
  k=jbin(x,tthlow,tthhigh,step)
! test binning ! write(*,*)x,data(1,k)
  if(k.gt.0)then
    n=n+1
    diff=y-data(2,k)
    wd=diff*diff/(esd*esd+data(3,k)*data(3,k))
    c2=c2+wd
    if(wd.gt.9.0d0)n3s=n3s+1
    if(wd.gt.27.0d0)n6s=n6s+1
  endif
  goto 2
22 write(*,'(a)')'Data format problem in '//line(1:len_trim(line))
  stop
21 close(20)
  write(*,1001)line(1:len_trim(line)),n,n3s,n6s,c2/real(n,8)

```



```

1001  format(a10,3i10,F16.8)
      enddo
      write(ioutunit,'(a)')'$ END'
      close(ioutunit)
      end program id31check

      integer(kind=4) function jbin(tth,tthlow,tthhigh,step)
      implicit none
      real(kind=8),intent(in)::tth, tthlow, tthhigh, step
      real(kind=8)::x
      if(tth.ge.tthlow .and. tth.le.tthhigh)then
        x=(tth-tthlow)/step+0.5d0
        jbin=nint(x)
      else
        jbin=-1
      endif
      return
      end function jbin

```

◇

## 9.7 Merging inp files

When zapping outliers there is a problem that when adding many scans together the effect of the zinger is reduced. In an individual scan they are easier to remove. The idea is to run id31sumall generating a series of .inp files and then to merge those together in a little program called id31inpsum. This is very similar to the program for checking that channels are equivalent.

```

"id31inpsum.f90" 95≡
      program id31inpsum
      implicit none
      character(len=256)::line
      integer :: ier, iarg, ioutunit, jjj
      integer, external :: iargc
      integer(kind=4) :: i, ndata, iscan, n, j, k, n3s, n6s
      real(kind=8) :: x,y,esd, diff, c2, wd, tthlow, tthhigh, step
      real(kind=8),allocatable,dimension(:,,:) :: data
      integer(kind=4),external :: jbin
! Read epf file in
      i=0; ier=0
      call getarg(1,line) ! name of first file
      open(unit=20,file=line(1:len_trim(line)),status='OLD',iostat=ier)
      if(ier.ne.0)then
        write(*,'(a)')'Error opening '//line(1:len_trim(line))
        write(*,'(a)')'Usage: id31inpsum file1 file2 file3 ...'
        stop
      endif
      write(*,'(a)')'Reading'//line(1:len_trim(line))
      read(20,*,end=10,err=11)x0,y0,esd0
      read(20,*,end=10,err=11)x,y,esd
      step=x-x0
      nbins=360.0d0/step
10    allocate(data(3,nbins),stat=ier)
      if(ier.ne.0)stop 'Memory allocation problem'
      i=i+1
      goto 1
11    stop 'Data format problem'
!

```

```

        ndata=i
        rewind(21)
        do i=1,ndata
            read(21)data(1,i),data(2,i),data(3,i)
        enddo
        close(21)
        close(20)
        jjj=1; ioutunit=19
! write a file in plotmtv format with the 9 channels and sum
        open(unit=ioutunit,file='scans.mtv',status='UNKNOWN',
        & form='FORMATTED', access='SEQUENTIAL',iostat=ier)
        if(ier.ne.0)stop 'error opening diagnostic file'
        write(ioutunit,'(a)')'$ DATA = CURVE2D'
        write(ioutunit,'(a)')'% xlabel = "Two Theta"'
        write(ioutunit,'(a)')'% ylabel = "Cts/Monitor"'
        write(ioutunit,'(a)')'% toplabel= "Scans and total plot"'
        write(ioutunit,'(a,i2,a)')'% linelabel = "Total" '
        write(ioutunit,'(a,i2)')'% linecolor = ',jjj
        write(ioutunit,'(a)')'% linetype=1 markertype=0'
        do j=1,ndata
            write(ioutunit,'(2F15.8)')data(1,j),data(2,j)
        enddo
        write(ioutunit,*)
! Epf file is now in array data
        step=(data(1,ndata)-data(1,1))/real(ndata-1,8)
        tthlow=data(1,1)-step/2.0d0
        tthhigh=data(1,ndata)+step/2.0d0
        write(*,*)'Opened summed file ',line(1:len_trim(line)), ' npts ',
        & ndata
        write(*,1000)'Step ',step,' tthlow ',tthlow,' tthhigh ',tthhigh
1000 format(3(a,F15.8))
! Now read in the series of scans
        iarg=iargc() ! total number of cmdline arguments
        if(iarg.le.1)then
            stop 'need to supply a list of scans to check'
        endif
!
        1234567890123456789012345678901234567890
        write(*,999)
999 format('      Scan      npts   3-sigma  6-sigma    chi**2')
        do i=2,iarg
            call getarg(i,line)
            line=adjustl(line)
            n=len_trim(line)
            write(line,'(a)')line(1:n)//'.inp'
            open(unit=20,file=line(1:len_trim(line)),status='OLD',iostat=ier)
            if(ier.ne.0)then
                write(*,'(a)') 'Error opening '//line(1:len_trim(line))
                stop
            endif
!
            write(*,*)'Opened '//line(1:len_trim(line))
            write(ioutunit,*)
            write(ioutunit,'(a,a,a)')'% linelabel = "',
            & line(1:len_trim(line)),'"'
            jjj=jjj+1
            write(ioutunit,'(a,i3)')'% linecolor = ',jjj
            write(ioutunit,'(a)')'% linetype=1 markertype=0'
            c2=0.0d0; n3s=0; n6s=0; n=0
! read and check the .inp against the .epf now
2      read(20,*,end=21,err=22)x,y,esd
            write(ioutunit,'(2F15.8)')x,y

```

```

        k=jbin(x,tthlow,tthhigh,step)
! test binning !      write(*,*)x,data(1,k)
        if(k.gt.0)then
            n=n+1
            diff=y-data(2,k)
            wd=diff*diff/(esd*esd+data(3,k)*data(3,k))
            c2=c2+wd
            if(wd.gt.9.0d0)n3s=n3s+1
            if(wd.gt.27.0d0)n6s=n6s+1
        endif
        goto 2
22  write(*,'(a)')'Data format problem in'//line(1:len_trim(line))
    stop
21  close(20)
    write(*,1001)line(1:len_trim(line)),n,n3s,n6s,c2/real(n,8)
1001 format(a10,3i10,F16.8)
    enddo
    write(ioutunit,'(a)')'$ END'
    close(ioutunit)
end program id31check

integer(kind=4) function jbin(tth,tthlow,tthhigh,step)
implicit none
real(kind=8),intent(in)::tth, tthlow, tthhigh, step
real(kind=8)::x
if(tth.ge.tthlow .and. tth.le.tthhigh)then
    x=(tth-tthlow)/step+0.5d0
    jbin=nint(x)
else
    jbin=-1
endif
return
end function jbin

```

◇

## 9.8 AOB

Are there any other main programs which we need? binit, binem and resum are replaced so far. We still need flat plate versions or flags in user options. Do we still need a peak fitting thing for the strain scanners?

## 10 Tidying up

This routine should close all opened files and free all allocated memory. Use grep on the main source file (binit.w) to be sure to catch all possible allocations and opened files.

```

< tidyup 98 > ≡
      subroutine tidyup
      ! Free all allocated memory and close all opened files
      use specfiles      ! iunit
      use rebin          ! ascan
      use summation      ! hist, sumdata
      use outputfiles    ! ioutunit
      use useroptions    ! exdetlist, exscanlist
      if(allocated(ascan  ))deallocate(ascan)
      if(allocated(hist   ))deallocate(hist )
      if(allocated(sumdata ))deallocate(sumdata)
      if(allocated(exdetlist ))deallocate(exdetlist)
      if(allocated(exscanlist))deallocate(exscanlist)
      ! do an inquire on these to see if they are actually open
      close(iunit)      ! specfile
      close(ioutunit) ! output data
      close(logfile)    ! logfile
      close(16) ! temp.res unit
      return
      end subroutine tidyup

```

Fragment referenced in 40a, 57b, 89, 91.

## 11 Peak Fitting

A fortran 90 module for fitting peaks to data. The intention is to generate a routine which can be called with an array of data as an argument and it will fit a peak to that data

### 11.1 Introduction

Fitting a single peak to a dataset means defining a peak function and a background function, estimating initial values for any parameters and carrying out a non-linear least squares fit of the parameters to the data.

### 11.2 Peak function

Initially we will just use Larry Finger's low angle asymmetry correcting peak shape. This acts as a conventional pseudo-voigt when the asymmetry parameters are set to zero. A general interface to a peak shape function is required. We need to supply some number of parameters and the  $2\theta$  value at which the peak function is required. The function (or subroutine) should return the function value and it's derivatives with respect to each of the parameters. A generalised block is as follows:

```

⟨ peakfunction 99a ⟩ ≡
  subroutine peakfunction(x,phi,pars,n,dpars)
  real, intent(in) :: x      ! two theta value
  real, intent(out):: phi    ! function value
  integer, intent(in):: n    ! Number of parameters
  real, intent(in), dimension(n) :: pars ! parameters
  real, intent(out), dimension(n) :: dpars
  real(kind=4):: Eta , Gamma , S_L , D_L , TwoTH
  real(kind=4):: TwoTH0 , dPRdT, dPRdG, dPRdE , dPRdS , dPRdD
! also uses asy module variable
  logical Use_asy
  real(kind=4)::Profval
  external Profval
  TwoTH0=pars(1);Gamma=pars(2);Eta=pars(3)
  S_L=pars(4);D_L=pars(5); TwoTH=x
  if(asy.and.(S_L.gt.0.0 .or. D_L.gt.0.0))Use_Asym=.TRUE.
! Optimise to cut at 50 fwhm tails
  if(abs(TwoTH0-TwoTH).lt.50.0*Gamma)then
    phi = Profval(Eta,Gamma,S_L,D_L,TwoTH,TwoTH0 ,
    &          dPRdT, dPRdG, dPRdE , dPRdS , dPRdD , Use_Asym)
    dpars(1)=dPRdT;dpars(2)=dPRdG;dpars(3)=dPRdE
    dpars(4)=dPRdS;dpars(5)=dPRdD
  else
    phi=0.; dpars(1:5)=0.
  endif
  return
end subroutine peakfunction

```

Fragment referenced in 101b.

This needs to link to the **profval.f** file containing the function.

### 11.3 Background function

A background function has the same interface as a peak function, just the code and the parameters will be different. It is a polynomial defined by:

$$y_{back} = \sum_{i=0}^n a_i x^i$$

$$\frac{dy_{back}}{da_i} = x^i$$

Other functions can clearly be added later

```

⟨ backfunction 99b ⟩ ≡
  subroutine backfunction(x,phi,pars,n,dpars)
  real, intent(in) :: x      ! two theta value
  real, intent(out):: phi    ! function value
  integer, intent(in):: n    ! Number of parameters
  real, intent(in), dimension(n) :: pars ! parameters
  real, intent(out), dimension(n) :: dpars
  integer :: i                ! derivs of phi wrt parameters
  phi=0.0; dpars=0.0 ! set out args to zero
  do i=1,n
    phi=phi+pars(i)*x**(i-1)      ! polynomial
    dpars(i)=x**(i-1)            ! derivatives
  enddo
  return
end subroutine backfunction

```

Fragment referenced in 101b.

## 11.4 Defining the problem

Forming the observed, calculated and difference patterns, and the derivatives of the pattern needs to be repeated for each cycle and is placed in a subroutine called **calpat** which takes care of setting up the least squares matrix. Formally the problem is to minimise

$$\chi^2 = \sum_i \frac{(y_{obs} - y_{calc})^2}{\sigma_i}$$

where

$$y_{calc} = y_{back} + s\phi(2\theta, \Gamma, \eta, S/L, D/L)$$

and  $s$  is a scale factor and  $\phi$  is the peak function. At a minimum the derivative of  $\chi^2$  with respect to any model parameter will be zero and we follow a standard least squares procedure to get there. [ elaborate sometime - essentially form matrix of  $d\chi^2/da_i * d\chi^2/da_j$ . This is an approximation to the second derivatives (products of first derivatives, the true second derivative terms are multiplied by  $y_{obs} - y_{calc}$  and so they theoretically will drop out if the model works. A vector of first derivatives is also formed as  $d\chi^2/da_i$  and we solve the equation  $\mathbf{Ax} = \mathbf{b}$ , where  $\mathbf{A}$  is the matrix and  $\mathbf{b}$  is the vector of first derivatives times differences between observed and calculated. Inverting the matrix gives an inverse and so we multiply both sides by that to get the shift vector,  $\mathbf{x}$ . Esd's are given by  $\sqrt{\chi^2 A_{ii}}$ . ]

FIXME - Tidy plus add  $\chi^2$  reasoning about altering the weights to suit the original assumptions.

The **calpat** takes care of calculating the matrix and right hand side vector, as well as the  $y_{calc}$  ready to get a set of shifts.

```
< calpat 100 > ≡
subroutine calpat
integer :: i, j, k
real :: wt
do i=1,ndata
  d=0.0 ! set derivs to zero
  call backfunction(data(1,i),yback,pars(bk1:bk2),nbkpars,d(bk1:bk2))
  call peakfunction(data(1,i),ypeak,pars(pos:d_1),npkpars,d(pos:d_1))
  ycalc = yback + pars(scal) * ypeak ! calc and diff
  d(pos:d_1)=d(pos:d_1) * pars(scal) ! correct peaks derivs for scale
  d(scal) = ypeak ! deriv w.r.t to scal
  ydiff = data(2,i) - ycalc ! obs - calc
  fit(i)=ycalc ! store calc
  if(gau)d(eta)=0.
  if(posonly)d(pos:d_1)=0.
  if(poswd)d(s_1:d_1)=0.
  wt=1./data(3,i)**2
! make wt a function of difference like WIFD ??
  do j=1,npars ! sum in lsq contributions
    do k=1,npars
      amat(j,k)=amat(j,k)+d(j)*d(k)*wt
    enddo
    bvec(j)=bvec(j)+ydiff*d(j)*wt
  enddo
  chi2=chi2+ydiff*ydiff*wt
enddo
chi2=chi2/real(ndata-npars) ! make reduced chi2
return; end subroutine calpat
```

◇

Fragment referenced in 101b.

## 11.5 Marquardt Damping

```

<marq 101a> ≡
! Apply Marquardt damping (effectively observations of no change)
subroutine marq
! fac, amult, icyc, chi2, chiold, npars, amat
integer::i
if(icyc.eq.0)fac=0           ! icyc=0 is last cycle flag
if(icyc.gt.2) then
  if(chimin.gt.chi2)fac=fac/amult
  if(chimin.le.chi2 .and.icyc.gt.3)then
    fac=fac*amult**2        ! go back two lambda steps and make smaller
    amult=1.+amult/2.      ! reducing steps on lambda >= 1
  endif
endif
if(icyc.gt.1)then
  do i=1,npars; amat(i,i)=amat(i,i)*(1.0+fac) ; enddo
endif
end subroutine marq      ◇

```

Fragment referenced in 101b.

## 11.6 Modules

```

<pkfit 101b> ≡
module pkfit
! fac, amult, icyc, chi2, chiold, npars, amat
integer, parameter :: npars=8 , nbkparms=2, npkparms=5
character(len=10), dimension(npars) :: names
data names /'Intensity ', 'Const Bk ', 'Slope Bk ', 'Position ', &
&'Width ', 'Eta ', 'S/L ', 'D/L ' /
integer, parameter :: scal=1, bk1=2, bk2=3, pos=4, wid=5, eta=6, &
& s_1=7, d_1=8
real, dimension(npars) :: bvec, d, s, e, pars, minpars
real, dimension(npars,npars) :: amat
real, dimension(:),allocatable :: fit
real,allocatable :: data(:, :)
real, parameter :: facst=10., ast=10.0
real :: ycalc, ypeak, yback, ydiff, wydiff
real :: fac, amult, chi2, chiold, chimin
logical :: asy=.false. ! whether or not to use asymetry
logical :: gau=.false. ! Gaussian only - eta is fixed
logical posonly, poswd
integer :: icyc, ndata ! number of data points
contains
<backfunction 99b>
<peakfunction 99a>
<calpat 100>
<marq 101a>
<matscal 104>
<lsqmagic 102b>
<estimatepars 105>
<invert 106>
<updatepars 102a>
end module pkfit      ◇

```

Fragment referenced in 107.

```

< updatepars 102a > ≡
      subroutine updatepars
      real :: step
      pars=pars+bvec
      if(pars(s_1).lt.0.)pars(s_1)=1.0e-6
      if(pars(d_1).lt.0.)pars(d_1)=1.0e-6
      step=abs(data(1,1)-data(1,2))
      if(pars(wid).lt.step)pars(wid)=step
      if(pars(eta).lt.-0.5)pars(eta)=-0.5
      if(pars(eta).gt.2.0)pars(eta)=2.0 ! hard limits
      if(chi2.gt.chimin)then
        pars=(minpars+pars)/2.
      endif
      end subroutine updatepars

```

Fragment referenced in 101b.

## 11.7 Forming the least squares matrix

We will assume that a function will be called with a data array as an argument. This will then fill in a least squares matrix by forming the calculated function at each data point, and using the difference and derivatives.

```

< lsqmagic 102b > ≡
      subroutine lsqmagic()
      integer :: i, j
      allocate(fit(ndata))
      call estimatepars(pars(bk1:bk2),pars(pos:d_1),pars(scal))
      fac=facst; amult=ast; chiold=-1.; icyc=0 ; poswd=.true.; posonly=.true.
      write(*,*)fac
1      icyc=icyc+1
      amat=0.0; bvec=0.0; chi2=0.
      call calpat
      call marq
      call matscal
      call invert      ! overwrites amat with inverse
      s=0.
      do i=1,npars
        do j=1,npars; s(i)=s(i)+amat(i,j)*bvec(j) ; enddo
        if(d(i).gt.0)then ; s(i)=s(i)/d(i) ; else ; s(i)=0. ;endif
        if(d(i).lt.0.)s(i)=0.
      enddo
      do i=1,npars
        if(d(i).gt.0)then
! Scale back to original matrix
          amat(:,i)=amat(:,i)/d(i); amat(i,:)=amat(i,:)/d(i)
        else
          amat(:,i)=0.; amat(i,:)=0.
        endif
      enddo
      bvec=s ; e=0. ; s=0.
      do i=1,npars
        if(d(i).gt.0.) then
          e(i)=sqrt(chi2*amat(i,i))
          s(i)=bvec(i)/e(i)
        endif
      enddo
      write(*,1000)icyc,chi2,fac,amult,maxval(s)
1000 format('Cycle ',i5,' Chi2 ',F15.6,3E15.6)

```



```

if(icyc.gt.0) then ! not last cycle so update pars
  if(icyc.eq.1.and.posonly)then
    chimin=chi2 ; minpars=pars
  endif
  chiold=chi2
  if(chi2.lt.chimin)then
    chimin=chi2
    minpars=pars ! save best parameters found
  endif
  call updatepars
  if(maxval(s).lt.0.01 .and. icyc.gt.5)then ! converged
    if(posonly)then
      posonly=.false. ; fac=facst ; amult =ast ; icyc=1
      write(*,*)'width now varying'
    else if(poswd) then
      poswd=.false. ; fac=facst ; amult = ast ; icyc=1
      write(*,*)'all parameters now varying'
    else
      write(*,*)'Going to minpars'
      pars=minpars ! use best set found
      icyc=-1
    endif
  endif
  endif
  if(icyc.gt.100) icyc=-1 ! run out of patience
  goto 1
else
  write(*,1001)chi2,icyc
1001  format('Perhaps converged somewhere with Chi2 = ',F15.8,1x,i5)
  write(21,*)
  write(21,'(a)')'% linecolor=3'
  do i=1,ndata
    write(21,*)data(1,i),fit(i)
  enddo
  write(21,*)
  write(21,'(a)')'% linecolor=5'
  do i=1,ndata
    write(21,*)data(1,i),data(2,i)-fit(i)
  enddo
  write(21,'(a)')'$ END'
  close(21)
  do i=1,npars
    j=i
    write(*,1002)names(i),pars(i),e(i),s(i)
  enddo
1002  format(a10,G15.8,' +/- ',G15.8,' with sh/esd =',G15.8)
  endif
  return
end subroutine lsqmagic

```

◇

Fragment referenced in 101b.

```

⟨ matscal 104 ⟩ ≡
  ! Scale matrix and take out non variable parameters
  subroutine matscal
    integer::i
    do i=1,npars
      if(amat(i,i).gt.0.) then
        d(i)=sqrt(amat(i,i))
        amat(:,i)=amat(:,i)/d(i); amat(i,:)=amat(i,:)/d(i)
        bvec(i)=bvec(i)/d(i)
      else
        d(i)=-1.0; amat(:,i)=0. ; amat(i,:)=0. ; amat(i,i)=1.0
      endif
    enddo
  end subroutine matscal
  ◇

```

Fragment referenced in 101b.

## 11.8 Parameter estimation

We could try to form some quite elegant approximations to guessing the peak before fitting. Initially something rather awful is done.

```

< estimatepars 105 > ≡
      subroutine estimatepars(bkpars, pkpars, scal)
      integer, parameter :: nbkpars=2, npkpars=5, npars=8 ! 1+5+2
      real, dimension(nbkpars), intent(out) :: bkpars
      real, dimension(npkpars), intent(out) :: pkpars
      real, intent(out) :: scal
      integer :: i, mxdata(1)
      real :: y, y1, y2
      ! background is linear a+bx+... so make a=0 and b=minimum value
      bkpars(1)=minval(data(2,1:ndata))
      bkpars(2:nbkpars)=0.0
      ! peak pars are tth0, width, eta, s/l and h/l
      ! tth0 - weighted average of data
      s=0.; sy=0.; syt=0.; syt2=0.
      mxdata=maxloc(data(2,:))
      pkpars(1)=data(1,mxdata(1))      ! centre is max of scan
      ! Width from (mxdata back(1))/2 and walk array up and down to find 1/2
      y=(data(2,mxdata(1))+bkpars(1))/2.0
      ! walk array up
      do i=mxdata(1),ndata
         y1=data(1,i)
         if(data(2,i).lt.y) exit ! next loop
      enddo
      do i=mxdata(1),1,-1
         y2=data(1,i)
         if(data(2,i).lt.y) exit
      enddo
      pkpars(2)=abs(y1-y2)
      write(*,*) 'tth guess= ', pkpars(1), pkpars(2)
      pkpars(3)=0.9 ! guess a default eta ?
      if(gau) pkpars(3)=0.0
      if(asy) then
         pkpars(4)=0.001
         pkpars(5)=0.002
      else
         pkpars(4)=0.0
         pkpars(5)=0.0
      endif
      scal=1.
      return
      end subroutine estimatepars

```

Fragment referenced in 101b.

## 11.9 Matrix inversion

We invert the least squares matrix using a singular value decomposition from the IMSL library (available for both windows and linux, hopefully). It's actually a generalised inverse to save me book keeping.

The **invert** routine takes a matrix and vector as arguments and is expected to return the inverse matrix and the solution vector, both overwriting their arguments.

```

< invert 106 > ≡
      subroutine invert
      !use imsl ! to get lsgr
      real, dimension(npars,npars) :: ainv
      integer :: i
      ainv=0. ; i=1
      do i=1,5;call eraset(i,1,0);enddo ! don't give up!
      call lsgr(npars,npars,amat,npars,1.0e-5,i,ainv,npars)
!      print *, "Returned from LSGRR"
      amat=ainv
      return
      end subroutine invert

```

◇

Fragment referenced in 101b.

### 11.10 Driver program

This program just needs to obtain an x,y,esd array of data from somewhere and call the least squares magic routine to fit that data.

```

"fitit.f90" 107≡
  <pkfit 101b>
    program fitit
      use pkfit
      real x,y,esd,start_time,end_time, lowtth, hightth
      integer i,j
      character(len=256) :: line
! read from stdin to a scratch file
      call cpu_time(start_time)
      i=0
      asy=.false.; lowtth=-360.0 ; hightth=360.0
      call getarg(1,line)
      read(line,*,end=20,err=20)lowtth
      call getarg(2,line)
      read(line,*,end=20,err=20)hightth
      call getarg(3,line)
      if(line(1:2).eq.'+a') asy=.true.
      if(line(1:2).eq.'+g') gau=.true.
20    open(unit=20,status='SCRATCH',form='UNFORMATTED')
      open(unit=21,status='UNKNOWN',form='FORMATTED',
& ACCESS='SEQUENTIAL',FILE='temp.mtv')
!      open(unit=19,file='pk.dat')
!      do j=1,5
!        read(*,'(a256)')line
!        write(21,'(a)')line(1:len_trim(line))
!      enddo
1    read(*,*,end=10,err=100)x,y,esd
      if(x.gt.lowtth .and. x.lt.hightth) then
        write(20)x,y,esd
        write(21,*)x,y
        i=i+1
      endif
      goto 1
! Copy from scratch to allocatable data array
10    allocate(data(3,i),stat=j)
      ndata=i
      if(j.ne.0)stop 'Memory allocation error'
      rewind(20)
      do j=1,i
        read(20)data(1,j),data(2,j),data(3,j)
      enddo
      close(20)
      write(*,*)'No of datapoints=',ndata
! Data is read in, now do the fitting
      call lsqmagic()
      deallocate(data)
      call cpu_time(end_time)
      write(*,1000)end_time-start_time
1000  format('Time taken was ',F15.3,'/s')
      goto 101
100  STOP 'Data format problem'
101  end program fitit

```

## 12 Sifit

A short program to fit the wavelength and zero shift given a set of measured peak positions for silicon. We will hard wire in the hkl values and lattice parameter. The program will receive data in the following

format:

|           |          |          |          |         |         |         |         |
|-----------|----------|----------|----------|---------|---------|---------|---------|
| 13.999281 | 0.000060 | 0.004433 | 0.000105 | 1.04809 | 0.03724 | 1.00000 | 0.00000 |
| 22.957993 | 0.000058 | 0.004236 | 0.000100 | 1.27979 | 0.04473 | 8.00000 | 0.00000 |
| 26.989207 | 0.000060 | 0.005161 | 0.000106 | 1.06040 | 0.03314 | 1.00000 | 0.00000 |
| 32.691785 | 0.000104 | 0.006308 | 0.000193 | 0.99243 | 0.05133 | 5.00000 | 0.00000 |
| 35.719265 | 0.000097 | 0.007127 | 0.000160 | 0.81956 | 0.03477 | 5.00000 | 0.00000 |
| 40.324606 | 0.000090 | 0.008111 | 0.000152 | 0.76555 | 0.02770 | 5.00000 | 0.00000 |
| 42.887492 | 0.000127 | 0.008951 | 0.000210 | 0.73656 | 0.03638 | 5.00000 | 0.00000 |
| 46.906190 | 0.000173 | 0.009835 | 0.000263 | 0.65374 | 0.04489 | 5.00000 | 0.00000 |
| 49.194373 | 0.000137 | 0.010920 | 0.000205 | 0.58492 | 0.02972 | 5.00000 | 0.00000 |

The first two columns are the observed two theta and esd.

The literate programming business got shelved here as the author was in a hurry - essentially it just calculates a trial wavelength off of the last peak and then does least squares. The  $2 \times 2$  matrix inversion is done manually. Interesting things to write up in here will be the equation for the derivative of peak position in  $2\theta$  with respect to wavelength.

"sifit.f90" 108≡

```

program sifit
implicit none
integer, parameter :: npks=20
integer, dimension(3,npks) :: hkls
integer :: n, i, j, ncyc, nobs
data hkls / 1,1,1, 2,2,0, 3,1,1, 4,0,0, 3,3,1,      &
&          4,2,2, 5,1,1, 4,4,0, 5,3,1, 6,2,0,      &
&          5,3,3, 4,4,4, 5,5,1, 6,4,2, 5,5,3,      &
&          8,0,0, 7,3,3, 8,2,2, 5,5,5, 8,4,0      /
real(kind=8) :: alambda, dalam, zero, dzero ! lambda, deriv
real(kind=8), dimension(2,npks) :: peaks ! obs, esd
real(kind=8), dimension(2,npks) :: cpks ! calc peaks, d, tth
real(kind=8), dimension(npks) :: d ! differences obs-calc
real(kind=8), dimension(2,2) :: a, ai ! LSQ matrix & inverse
real(kind=8), dimension(2) :: b, x ! vectors in A.x=b
real(kind=8), parameter :: alatt=5.4311946d0 ! si lattice par
real(kind=8) :: PI, RAD ! convert deg to rad by multiplying
real(kind=8) :: t, chi2 ! temporary

!
peaks=0.0d0 ! initialise peaks array to zero
cpks=0.0d0
PI=2.0d0*datan2(1.0d0,0.0d0)
RAD=PI/180.0d0
! write(*,*)'RAD = ',RAD, ' PI= ',PI
i=1
! write(*, '(a)') ' h k l obs esd d'
1 read(*,*,end=2,err=2)peaks(1,i),peaks(2,i) ! fill in obs and esd
n=hkls(1,i)*hkls(1,i)+hkls(2,i)*hkls(2,i)+hkls(3,i)*hkls(3,i)
cpks(1,i)=1.0d0/dsqrt(real(n,8)/alatt/alatt) ! in d-spacing
! write(*,1000)(hkls(j,i),j=1,3),(peaks(j,i),j=1,2),cpks(1,i)
1000 format(3(i3,1x),3F16.8)
i=i+1
if(i.le.npks)then
goto 1
else
write(*,*)'Only the first ',npks,' peaks can be used'
endif
2 nobs=i-1
! Estimate lambda from highest hkl peak with zero as 0.0
zero=0.0d0; ncyc=0

```

```

        alambda=2.0d0*cpks(1,nobs)*dsin(RAD*peaks(1,nobs)/2.0d0)
!      write(*,*)'Initial guess lambda = ',alambda,' zero ', zero
3      a=0.0d0 ; x=0.0d0; b=0.0d0; chi2=0.0d0
! Fill in calc positions with this lambda and zero
      do i=1,nobs
        t=alambda/2.0d0/cpks(1,i)
        cpks(2,i)=2.0d0*dasin(t)/RAD - zero
        dalam=(1.0d0/cpks(1,i))*(1.0d0/dsqrt(1.0d0-t*t))/RAD
        dzero=-1.0d0
        dalam=dalam/peaks(2,i)
        dzero=dzero/peaks(2,i)
        a(1,1)=dalam*dalam+a(1,1)
        a(2,2)=dzero*dzero+a(2,2)
        a(2,1)=dzero*dalam+a(2,1)
        a(1,2)=dzero*dalam+a(1,2)
        d(i)=(peaks(1,i)-cpks(2,i))/(peaks(2,i)) ! wtd difference
        b(1)=b(1)+d(i)*dalam
        b(2)=b(2)+d(i)*dzero
        chi2=chi2+d(i)*d(i)
      enddo
! invert matrix, and apply shifts
      t=1.0d0/(a(1,1)*a(2,2)-a(1,2)*a(2,1))
        ! determinant - should catch singularity here
!      write(*,*)'Determinant = ',t
      ai(1,1)=a(2,2)*t
      ai(2,1)=-a(1,2)*t
      ai(1,2)=-a(2,1)*t
      ai(2,2)=a(1,1)*t
! fill in shifts in x
      x(1)=ai(1,1)*b(1)+ai(1,2)*b(2)
      x(2)=ai(2,1)*b(1)+ai(2,2)*b(2)
      alambda=alambda+x(1)
      zero=zero+x(2)
! more cycles? say 5 for now - seems to converge in 1 (is it linear?)
!      write(*,*)alambda,zero
!      write(*,*)'Cycle ',ncyc,' Chi2 ',chi2
      ncyc=ncyc+1
      if(ncyc.lt.6)goto 3
      write(*,*)' h k l obs esd calc '// &
      & 'diff/esd diff'
      do i=1,nobs
        write(*,2000)(hkls(j,i),j=1,3),peaks(1,i),peaks(2,i),cpks(2,i),d(i), &
          & peaks(1,i)-cpks(2,i)
2000 format(3(i3,1x),3(F12.7,1x),1F8.2,1x,F12.7)
      enddo
      write(*,'(a)')'Lambda = 2 * d * sin[(two_theta + zero)/2] '
      write(*,2001)'a(Si) ',alatt
2001 format(a7,F16.8,a,F16.8)
      write(*,2001)'Chi2 ',chi2 ', reduced Chi2', chi2/real(nobs-2,8)
!
      write(*,'(a)')'+-----+
      write(*,2002)'Lambda',alambda,' Zero',zero
      write(*,2002)' +/-',dsqrt(chi2*ai(1,1)), ' +/-', dsqrt(chi2*ai(2,2))
2002 format('|',a7,F16.8,8x,a7,F16.8,' |')
      write(*,'(a)')'+-----+
      write(*,2003)ai(1,2)/dsqrt(ai(1,1)*ai(2,2))
2003 format('Correlation coefficient of wavelength and zero = ',F16.8)
      end program sifit

```

◇

## 13 Code building

Building all of this code and documentation should be possible from a single monolithic file called binit.w, which contains mainly latex source with fortran interspersed. Running nuweb<sup>4</sup> produces a tex file and a series of fortran files, which all need to be compiled and checked. Some files for doing that checking are included here.

Firstly an ms-dos batch file for windows systems. We assume the file test.dat exists and is a valid SPEC file.

```
"doscheck.bat" 110a≡
echo off
set OPTS=/warn:all /assume:accuracy_sensitive /fast
nuweb binit.w
latex binit
nuweb binit.w
latex binit
dvi2pdf -p a4 binit
dvips binit
df %OPTS% /exe:id31sum.exe      id31sum.f90
df %OPTS% /exe:id31sumall.exe  id31sumall.f90
df %OPTS% /exe:id31offsets.exe id31offsets.f90
df %OPTS% /exe:sifit.exe       sifit.f90      ◇
```

```
"dos95check.bat" 110b≡
echo off
nuweb binit.w
echo id31sum
g95 id31sum.f90 -o test/id31sum.exe
echo id31offsets
g95 id31offsets.f90 -o test/id31offsets.exe
echo id31sumall
g95 id31sumall.f90 -o test/id31sumall.exe
rem df %OPTS% /exe:sifit.exe      sifit.f90      ◇
```

The unix equivalent of for building the programs is as follows. (-en means enforce standard -m0 means all warnings and comments, -lU77 means link to a library for getarg and iargc if necessary). Use it with the command **source absoftlinuxcheck** . The **-X -xstatic** is an option for the linker to get it to link the runtime libraries in statically. This appears to circumvent problems of different linux computers in the amber cluster having different libraries installed. It appears that amber-a and amber-b have problems compiling the code (to be investigated) but that amber-c is happy to do it. Probably this is due to different compiler versions... but provided the routines are linked statically (no libc dependencies) then they appear to work on all the amber linux machines. The -M200 suppresses the \$ edit descriptor warning.

```
"bld" 110c≡
#more /opt/intel/README.esrf
source /opt/intel/icsxe/2013.0.028/bin/ictvars.sh
# absoft f90 options were:
OPTS="-en -m0 -M200 -lU77 -O2 -X -static"
# ifort options are
export OPTS="-static -traceback -implicitnone"
echo "Compiling " $1 "with options" $OPTS
ifort -o $1 $OPTS $1.f90
strip $1      ◇
```

---

<sup>4</sup>A system for literate programming available from <http://nuweb.sourceforge.net>



```
"absoftlinuxcheck" 111a≡
    nuweb binit
    latex binit
    nuweb binit
    latex binit
    dvipdf binit
    dvips binit
    chmod a+x ./bld
    ./bld plotit
    ./bld columns
    ./bld c2xye
    ./bld plotmesh
    ./bld testbins1
    ./bld testbins2
    ./bld bindump
    ./bld testnextscan
    ./bld id31sum
    ./bld id31sumall
    ./bld id31offsets
    ./bld sifit
    ./bld id31check
    chmod a+x id31sumcheck
    /opt/intel/fc/current/bin/ifort
    f90 -c profval.f -o profval.o -O2
    f90 -o fitit fitit.f90 -en -m0 -lU77 -O2 -X -static profval.o -limsl
    strip fitit
```

◇

Some example commands to test the programs are as follows (either platform).

```
"testcommands.bat" 111b≡
    testspec test.dat > testspec.out
    testbins1      > testbins1.out
    testbins2      > testbins2.out
    echo bindumpedfile.out > in
    bindump test.dat 0.01 1 < in > bindump.out
    id31sum test.dat 0.01 1 1
    id31sumall test.dat 0.01 1 1
```

◇

At the time of writing (there might have been a full moon) all of the programs appeared to be working correctly on both platforms.

## 14 Indices

Source code files created from this document

```
"absoftlinuxcheck" Defined by 111a.
"bindump.f90" Defined by 40a.
"bld" Defined by 110c.
"c2xye.f90" Defined by 15.
"columns.f90" Defined by 14b.
"dos95check.bat" Defined by 110b.
"doscheck.bat" Defined by 110a.
"fitit.f90" Defined by 107.
"id31check.f90" Defined by 93b.
"id31inpsum.f90" Defined by 95.
```

"id3loffsets.f90" Defined by 57b.  
 "id3lsum.f90" Defined by 89.  
 "id3lsumall.f90" Defined by 91.  
 "id3lsumcheck" Defined by 93a.  
 "plotit.f90" Defined by 13.  
 "plotmesh.f90" Defined by 16.  
 "profval.f" Defined by 114.  
 "sifit.f90" Defined by 108.  
 "specfiles.f90" Defined by 11b.  
 "testbins1.f90" Defined by 31a.  
 "testbins2.f90" Defined by 32.  
 "testcommands.bat" Defined by 111b.  
 "testnextscan.f90" Defined by 86b.  
 "testspec.f90" Defined by 12.

#### Macros defined in this document

<alpset 82a> Referenced in 87a.  
 <backfunction 99b> Referenced in 101b.  
 <bcmfile 53> Referenced in 72.  
 <bcmset 84a> Referenced in 87a.  
 <bin 25> Referenced in 28.  
 <bincen 21b> Referenced in 28.  
 <bindumpmsg 40b> Referenced in 40a.  
 <calibsum 43a> Referenced in 56.  
 <calpat 100> Referenced in 101b.  
 <checkdets 50> Referenced in 56.  
 <checkrebinpars 22> Referenced in 28.  
 <cmdline 14a> Referenced in 13, 15, 16.  
 <convertunitfunction 30> Referenced in 28.  
 <crrfun 57a> Referenced in 57b.  
 <ctchan 51> Referenced in 50.  
 <dsort 136b> Referenced in 52.  
 <dumpscan 67> Referenced in 72.  
 <effic 43b> Referenced in 56.  
 <estimatepars 105> Referenced in 101b.  
 <exdet 75> Referenced in 87a.  
 <exfile 78> Referenced in 87a.  
 <exscan 79a> Referenced in 87a.  
 <filext 70b> Referenced in 72.  
 <filter 82c> Referenced in 87a.  
 <findscan 5> Referenced in 11a.  
 <getcmdline 73> Referenced in 87a.  
 <getdata 10b> Referenced in 11a.  
 <getfile 4> Referenced in 11a.  
 <getfirstpoint 70c> Referenced in 72.  
 <getheadervalue 9b> Referenced in 11a.  
 <getlastpoint 71a> Referenced in 72.  
 <getwd2 49> Referenced in 50, 64b, 66.  
 <gsasout 61> Referenced in 72.  
 <helpmsg 87b> Referenced in 40a, 57b, 86b, 89, 91.  
 <ibin 20a> Referenced in 28.  
 <id3loffsetmsg 59> Referenced in 57b.  
 <id3lsumallmsg 92a> Referenced in 91.  
 <id3lsummsg 90> Referenced in 89.  
 <incscan 79b> Referenced in 87a.  
 <initialiserebin 23> Referenced in 28.  
 <inlist 85b> Referenced in 87a.  
 <invert 106> Referenced in 101b.  
 <junk 9a> Not referenced.  
 <junk1 8> Not referenced.

<led 27> Referenced in 28.  
 <limits 80a> Referenced in 87a.  
 <logmotors 38> Referenced in 33.  
 <lsqmagic 102b> Referenced in 101b.  
 <marq 101a> Referenced in 101b.  
 <matscal 104> Referenced in 101b.  
 <medianchannels 52> Referenced in 56.  
 <medianchannelset 83b> Referenced in 87a.  
 <minmonset 80d> Referenced in 87a.  
 <minrenormset 81a> Referenced in 87a.  
 <minstepset 80c> Referenced in 87a.  
 <mma 36> Referenced in 33.  
 <ncommas 86a> Referenced in 87a.  
 <nextscan 85a> Referenced in 87a.  
 <normerr 48> Referenced in 56.  
 <OFFSET 17> Referenced in 28.  
 <offsetdefaults 18> Referenced in 23.  
 <openlogfile 69b> Referenced in 72.  
 <option 74> Referenced in 87a.  
 <outputdiagnostic 64b> Referenced in 72.  
 <outputdiagnosticw32 66> Referenced in 72.  
 <outputepf 68a> Referenced in 72.  
 <outputfiles 72> Referenced in 40a, 57b, 86b, 89, 91.  
 <outputfilesvars 60a> Referenced in 72.  
 <outputformats 71b> Referenced in 72.  
 <outputinp 69a> Referenced in 72.  
 <outputxye 68b> Referenced in 72.  
 <pdsout 64a> Referenced in 72.  
 <peakfunction 99a> Referenced in 101b.  
 <pkfit 101b> Referenced in 107.  
 <pointfilter 41> Referenced in 33, 86b.  
 <processline 26> Referenced in 28.  
 <processscan 33> Referenced in 40a, 57b, 89, 91.  
 <random 136a> Referenced in 28.  
 <rdnums 7b> Referenced in 11a.  
 <readheader 6> Referenced in 11a.  
 <rebin 28> Referenced in 31a, 32, 40a, 57b, 86b, 89, 91.  
 <renormset 83c> Referenced in 87a.  
 <renset 77> Referenced in 72.  
 <report 31b> Referenced in 31a.  
 <reporteffic 45> Referenced in 56.  
 <reportsums 39> Referenced in 33.  
 <rstset 76> Referenced in 72.  
 <scale 84c> Referenced in 87a.  
 <scalpk 55b> Referenced in 56.  
 <scaltot 55a> Referenced in 56.  
 <SCAN 19> Referenced in 28.  
 <scantitle 60b> Referenced in 61, 63, 64a.  
 <setmonitorcol 81c> Referenced in 87a.  
 <snblset 84b> Referenced in 87a.  
 <specfiles 11a> Referenced in 11b, 12, 13, 14b, 15, 16, 31a, 32, 40a, 57b, 86b, 89, 91.  
 <SPECVARS 3> Referenced in 11a.  
 <spfout 63> Referenced in 72.  
 <split 7a> Referenced in 11a.  
 <summation 56> Referenced in 40a, 57b, 86b, 89, 91.  
 <sumthem 47> Referenced in 56.  
 <superzapem 54> Referenced in 56.  
 <superzapset 83a> Referenced in 87a.  
 <tempreswrite 46> Referenced in 56.  
 <tidyup 98> Referenced in 40a, 57b, 89, 91.

<tthhb 20b> Referenced in 28.  
 <tthlb 21a> Referenced in 28.  
 <unitsoptions 29> Referenced in 87a.  
 <updatepars 102a> Referenced in 101b.  
 <useroptions 87a> Referenced in 40a, 57b, 86b, 89, 91.  
 <wdset 80b> Referenced in 87a.  
 <whichcolumn 10a> Referenced in 11a.  
 <window 37> Referenced in 28.  
 <windowset 81b> Referenced in 87a.  
 <wlogfile 70a> Referenced in 72.  
 <zapping 82b> Referenced in 87a.  
 <zerodata 92b> Referenced in 91.

User specified index entries. Fill these in please!

Things to do still:

glitch and zinger elimination flat plate corrections reprocess for error determinations and  $R_{merge}$  temperature stuff

## 15 Appendix 1

Larry Finger's peak function which corrects for low angle asymmetry.

```

"profval.f" 114≡
  C LEVEL 7          SUBROUTINE PROFVAL()
      real*4 function Profval( Eta , Gamma , S_L , D_L , TwoTH ,
1      TwoTH0 , dPRdT, dPRdG, dPRdE , dPRdS , dPRdD , Use_Asym )
  c
  c 115.19 bugfix jpw 15-oct-2001
  c
  c Returns value of Profile
  c Eta is the mixing coefficient between Gaussian and Lorentzian
  c Gamma is the FWHM
  c S_L is source width/detector distance
  c D_L is detector width/detector distance
  c TwoTH is point at which to evaluate the profile
  c TwoTH0 is two theta value for peak
  c dPRdT is derivative of profile wrt TwoTH0
  c dPRdG is derivative of profile wrt Gamma
  c dPRdE is derivative of profile wrt Eta
  c dPRdS is derivative of profile wrt S_L
  c dPRdD is derivative of profile wrt D_L
  c Use_Asym is true if asymmetry to be used
  c
  c
  c Asymmetry due to axial divergence using the method of Finger, Cox and
  c Jephcoat, J. Appl. Cryst. 27, 892, 1992.

      implicit none
      real*4 Eta , Gamma , S_L , D_L , TwoTH
      real*4 TwoTH0 , dPRdT, dPRdG, dPRdE , dPRdS , dPRdD
      logical Use_Asym
      integer*4 NTERMS(14)/6,10,20,40,60,80,100,150,200,300,400,
1      600,800,1000/
  C 115.19      integer*4 Fstterm(14)/0,3,8,18,38,68,108,158,233,333,483,
  C changed first term to 1 instead of zero.
      integer*4 Fstterm(14)/1,3,8,18,38,68,108,158,233,333,483,
1      683,983,1383/
      real*4 RAD/57.2957795/

      integer*4 ArrayNum , K , NGT, ngt2 , it, i
  
```

```

real*4 CsTH          ! cos(theta)
real*4 TTH           ! tan(theta)
real*4 SnTwoTH       ! sin(twoth)
real*4 CsTwoTH       ! cos(twoth)
real*4 ApB           ! (S + H)/L
real*4 AmB           ! (S - H)/L
real*4 ApB2          ! (ApB) **2
real*4 Einfl         ! 2phi value for inflection point
real*4 Emin          ! 2phi value for minimum
real*4 dEminDA       ! derivative of Emin wrt A
real*4 tmp , tmp1 , tmp2 ! intermediate values
real*4 WP(1883) , XP(1883)! Storage for Gauss-Legendre weights and intervals
real*4 Delta         ! Angle of integration for convolution
real*4 dDELTA dA     ! derivative of DELTA wrt A (S/L)
real*4 sinDELTA      ! sine of DELTA
real*4 cosDELTA      ! cosine of DELTA
real*4 tanDELTA      ! tangent of DELTA
real*4 RcosDELTA     ! 1/cos(DELTA)
real*4 F , dFdA
real*4 G , dGdA , dGdB , PsVoigt
real*4 sumWG , sumWRG , sumWdGdA , sumWRdGdA , sumWdGdB , sumWRdGdB
real*4 sumWGdRdG , sumWGdRdE , sumWGdRdA , sumWGdRdB , sumWGdRd2t
real*4 stepsize

!
! Values for the abscissas and weights of the Gauss-Legendre
! N-point quadrature formula have been precomputed using routine
! Gauleg from "Numerical Recipes" (Press, Flannery, Teukolsky
! and Vetterling, 1986, Cambridge University Press,
! ISBN 0 521 30811 9), and are stored in the DATA statements
! for XP and WP below.
!
data (xp(i),i= 1, 40)/
$.2386192E+00,.6612094E+00,.9324695E+00,.1488743E+00,.4333954E+00,
$.6794096E+00,.8650634E+00,.9739065E+00,.7652652E-01,.2277859E+00,
$.3737061E+00,.5108670E+00,.6360537E+00,.7463319E+00,.8391170E+00,
$.9122344E+00,.9639719E+00,.9931286E+00,.3877242E-01,.1160841E+00,
$.1926976E+00,.2681522E+00,.3419941E+00,.4137792E+00,.4830758E+00,
$.5494671E+00,.6125539E+00,.6719567E+00,.7273183E+00,.7783057E+00,
$.8246122E+00,.8659595E+00,.9020988E+00,.9328128E+00,.9579168E+00,
$.9772599E+00,.9907262E+00,.9982377E+00,.2595977E-01,.7780933E-01/
data (xp(i),i= 41, 80)/
$.1294491E+00,.1807400E+00,.2315436E+00,.2817229E+00,.3311428E+00,
$.3796701E+00,.4271737E+00,.4735258E+00,.5186014E+00,.5622789E+00,
$.6044406E+00,.6449728E+00,.6837663E+00,.7207165E+00,.7557238E+00,
$.7886937E+00,.8195375E+00,.8481720E+00,.8745199E+00,.8985103E+00,
$.9200785E+00,.9391663E+00,.9557223E+00,.9697018E+00,.9810672E+00,
$.9897879E+00,.9958405E+00,.9992101E+00,.1951138E-01,.5850444E-01,
$.9740840E-01,.1361640E+00,.1747123E+00,.2129945E+00,.2509524E+00,
$.2885281E+00,.3256644E+00,.3623048E+00,.3983934E+00,.4338754E+00/
data (xp(i),i= 81, 120)/
$.4686966E+00,.5028041E+00,.5361459E+00,.5686713E+00,.6003306E+00,
$.6310758E+00,.6608599E+00,.6896376E+00,.7173652E+00,.7440003E+00,
$.7695024E+00,.7938327E+00,.8169541E+00,.8388315E+00,.8594314E+00,
$.8787226E+00,.8966756E+00,.9132631E+00,.9284599E+00,.9422428E+00,
$.9545908E+00,.9654851E+00,.9749091E+00,.9828486E+00,.9892913E+00,
$.9942275E+00,.9976499E+00,.9995538E+00,.1562898E-01,.4687168E-01,
$.7806858E-01,.1091892E+00,.1402031E+00,.1710801E+00,.2017899E+00,
$.2323025E+00,.2625881E+00,.2926172E+00,.3223603E+00,.3517885E+00/
data (xp(i),i= 121, 160)/
$.3808730E+00,.4095853E+00,.4378974E+00,.4657816E+00,.4932108E+00,

```

```

$.5201580E+00,.5465970E+00,.5725019E+00,.5978475E+00,.6226089E+00,
$.6467619E+00,.6702830E+00,.6931492E+00,.7153381E+00,.7368281E+00,
$.7575981E+00,.7776279E+00,.7968979E+00,.8153892E+00,.8330839E+00,
$.8499645E+00,.8660147E+00,.8812187E+00,.8955616E+00,.9090296E+00,
$.9216093E+00,.9332885E+00,.9440559E+00,.9539008E+00,.9628137E+00,
$.9707858E+00,.9778094E+00,.9838775E+00,.9889844E+00,.9931249E+00,
$.9962951E+00,.9984920E+00,.9997137E+00,.1043694E-01,.3130627E-01/
data (xp(i),i= 161, 200)/
$.5216195E-01,.7299491E-01,.9379607E-01,.1145563E+00,.1352667E+00,
$.1559181E+00,.1765016E+00,.1970082E+00,.2174290E+00,.2377550E+00,
$.2579774E+00,.2780874E+00,.2980762E+00,.3179352E+00,.3376556E+00,
$.3572289E+00,.3766466E+00,.3959001E+00,.4149811E+00,.4338813E+00,
$.4525925E+00,.4711065E+00,.4894151E+00,.5075106E+00,.5253849E+00,
$.5430303E+00,.5604390E+00,.5776036E+00,.5945165E+00,.6111703E+00,
$.6275579E+00,.6436720E+00,.6595056E+00,.6750519E+00,.6903041E+00,
$.7052554E+00,.7198995E+00,.7342299E+00,.7482404E+00,.7619248E+00/
data (xp(i),i= 201, 240)/
$.7752773E+00,.7882919E+00,.8009631E+00,.8132853E+00,.8252531E+00,
$.8368613E+00,.8481049E+00,.8589789E+00,.8694787E+00,.8795996E+00,
$.8893372E+00,.8986874E+00,.9076460E+00,.9162090E+00,.9243729E+00,
$.9321340E+00,.9394890E+00,.9464346E+00,.9529678E+00,.9590857E+00,
$.9647858E+00,.9700655E+00,.9749225E+00,.9793548E+00,.9833603E+00,
$.9869373E+00,.9900843E+00,.9927999E+00,.9950829E+00,.9969323E+00,
$.9983473E+00,.9993274E+00,.9998723E+00,.7834291E-02,.2350095E-01,
$.3916184E-01,.5481311E-01,.7045093E-01,.8607145E-01,.1016708E+00/
data (xp(i),i= 241, 280)/
$.1172453E+00,.1327909E+00,.1483040E+00,.1637806E+00,.1792170E+00,
$.1946095E+00,.2099541E+00,.2252472E+00,.2404850E+00,.2556638E+00,
$.2707798E+00,.2858293E+00,.3008086E+00,.3157141E+00,.3305421E+00,
$.3452890E+00,.3599510E+00,.3745247E+00,.3890065E+00,.4033927E+00,
$.4176799E+00,.4318646E+00,.4459432E+00,.4599124E+00,.4737686E+00,
$.4875086E+00,.5011288E+00,.5146260E+00,.5279969E+00,.5412382E+00,
$.5543465E+00,.5673188E+00,.5801518E+00,.5928424E+00,.6053874E+00,
$.6177838E+00,.6300285E+00,.6421185E+00,.6540509E+00,.6658228E+00/
data (xp(i),i= 281, 320)/
$.6774311E+00,.6888732E+00,.7001461E+00,.7112472E+00,.7221736E+00,
$.7329227E+00,.7434919E+00,.7538786E+00,.7640801E+00,.7740941E+00,
$.7839181E+00,.7935496E+00,.8029862E+00,.8122257E+00,.8212659E+00,
$.8301044E+00,.8387391E+00,.8471679E+00,.8553887E+00,.8633995E+00,
$.8711983E+00,.8787832E+00,.8861524E+00,.8933041E+00,.9002364E+00,
$.9069477E+00,.9134364E+00,.9197008E+00,.9257394E+00,.9315507E+00,
$.9371333E+00,.9424859E+00,.9476071E+00,.9524956E+00,.9571503E+00,
$.9615700E+00,.9657536E+00,.9697002E+00,.9734086E+00,.9768781E+00/
data (xp(i),i= 321, 360)/
$.9801078E+00,.9830968E+00,.9858445E+00,.9883502E+00,.9906132E+00,
$.9926330E+00,.9944091E+00,.9959410E+00,.9972285E+00,.9982712E+00,
$.9990687E+00,.9996210E+00,.9999281E+00,.5227245E-02,.1568116E-01,
$.2613337E-01,.3658271E-01,.4702806E-01,.5746827E-01,.6790220E-01,
$.7832871E-01,.8874665E-01,.9915490E-01,.1095523E+00,.1199377E+00,
$.1303101E+00,.1406682E+00,.1510109E+00,.1613371E+00,.1716456E+00,
$.1819354E+00,.1922054E+00,.2024543E+00,.2126811E+00,.2228846E+00,
$.2330638E+00,.2432175E+00,.2533446E+00,.2634441E+00,.2735147E+00/
data (xp(i),i= 361, 400)/
$.2835555E+00,.2935652E+00,.3035429E+00,.3134874E+00,.3233976E+00,
$.3332725E+00,.3431110E+00,.3529120E+00,.3626744E+00,.3723971E+00,
$.3820792E+00,.3917194E+00,.4013169E+00,.4108705E+00,.4203792E+00,
$.4298420E+00,.4392578E+00,.4486255E+00,.4579443E+00,.4672130E+00,
$.4764306E+00,.4855961E+00,.4947086E+00,.5037670E+00,.5127704E+00,
$.5217177E+00,.5306079E+00,.5394402E+00,.5482135E+00,.5569269E+00,
$.5655795E+00,.5741702E+00,.5826982E+00,.5911624E+00,.5995621E+00,

```

```

$.6078963E+00,.6161639E+00,.6243643E+00,.6324964E+00,.6405594E+00/
data (xp(i),i= 401, 440)/
$.6485524E+00,.6564744E+00,.6643248E+00,.6721025E+00,.6798068E+00,
$.6874367E+00,.6949916E+00,.7024704E+00,.7098725E+00,.7171970E+00,
$.7244432E+00,.7316101E+00,.7386971E+00,.7457033E+00,.7526281E+00,
$.7594705E+00,.7662300E+00,.7729057E+00,.7794970E+00,.7860030E+00,
$.7924232E+00,.7987567E+00,.8050030E+00,.8111612E+00,.8172308E+00,
$.8232111E+00,.8291014E+00,.8349011E+00,.8406095E+00,.8462260E+00,
$.8517501E+00,.8571811E+00,.8625184E+00,.8677614E+00,.8729095E+00,
$.8779623E+00,.8829191E+00,.8877794E+00,.8925427E+00,.8972084E+00/
data (xp(i),i= 441, 480)/
$.9017761E+00,.9062452E+00,.9106152E+00,.9148857E+00,.9190563E+00,
$.9231263E+00,.9270955E+00,.9309634E+00,.9347295E+00,.9383934E+00,
$.9419548E+00,.9454132E+00,.9487683E+00,.9520197E+00,.9551671E+00,
$.9582100E+00,.9611482E+00,.9639814E+00,.9667092E+00,.9693313E+00,
$.9718476E+00,.9742575E+00,.9765610E+00,.9787578E+00,.9808476E+00,
$.9828302E+00,.9847054E+00,.9864729E+00,.9881326E+00,.9896844E+00,
$.9911279E+00,.9924632E+00,.9936899E+00,.9948081E+00,.9958175E+00,
$.9967181E+00,.9975097E+00,.9981923E+00,.9987659E+00,.9992302E+00/
data (xp(i),i= 481, 520)/
$.9995854E+00,.9998313E+00,.9999680E+00,.3922075E-02,.1176598E-01,
$.1960917E-01,.2745115E-01,.3529144E-01,.4312955E-01,.5096502E-01,
$.5879735E-01,.6662606E-01,.7445067E-01,.8227070E-01,.9008566E-01,
$.9789509E-01,.1056985E+00,.1134954E+00,.1212853E+00,.1290678E+00,
$.1368423E+00,.1446083E+00,.1523655E+00,.1601134E+00,.1678513E+00,
$.1755790E+00,.1832958E+00,.1910013E+00,.1986951E+00,.2063767E+00,
$.2140456E+00,.2217013E+00,.2293434E+00,.2369713E+00,.2445847E+00,
$.2521830E+00,.2597658E+00,.2673327E+00,.2748830E+00,.2824165E+00/
data (xp(i),i= 521, 560)/
$.2899326E+00,.2974308E+00,.3049108E+00,.3123719E+00,.3198139E+00,
$.3272362E+00,.3346383E+00,.3420199E+00,.3493804E+00,.3567194E+00,
$.3640365E+00,.3713311E+00,.3786029E+00,.3858515E+00,.3930762E+00,
$.4002768E+00,.4074528E+00,.4146037E+00,.4217291E+00,.4288285E+00,
$.4359016E+00,.4429478E+00,.4499667E+00,.4569580E+00,.4639212E+00,
$.4708558E+00,.4777615E+00,.4846377E+00,.4914841E+00,.4983003E+00,
$.5050859E+00,.5118403E+00,.5185633E+00,.5252543E+00,.5319131E+00,
$.5385391E+00,.5451319E+00,.5516912E+00,.5582166E+00,.5647076E+00/
data (xp(i),i= 561, 600)/
$.5711639E+00,.5775851E+00,.5839707E+00,.5903203E+00,.5966337E+00,
$.6029103E+00,.6091498E+00,.6153519E+00,.6215161E+00,.6276420E+00,
$.6337293E+00,.6397777E+00,.6457866E+00,.6517559E+00,.6576850E+00,
$.6635737E+00,.6694215E+00,.6752281E+00,.6809932E+00,.6867164E+00,
$.6923974E+00,.6980357E+00,.7036311E+00,.7091832E+00,.7146916E+00,
$.7201561E+00,.7255763E+00,.7309518E+00,.7362823E+00,.7415676E+00,
$.7468072E+00,.7520008E+00,.7571482E+00,.7622490E+00,.7673029E+00,
$.7723096E+00,.7772688E+00,.7821801E+00,.7870433E+00,.7918581E+00/
data (xp(i),i= 601, 640)/
$.7966241E+00,.8013412E+00,.8060089E+00,.8106271E+00,.8151953E+00,
$.8197134E+00,.8241811E+00,.8285980E+00,.8329640E+00,.8372787E+00,
$.8415419E+00,.8457533E+00,.8499127E+00,.8540198E+00,.8580743E+00,
$.8620760E+00,.8660247E+00,.8699201E+00,.8737620E+00,.8775501E+00,
$.8812842E+00,.8849641E+00,.8885896E+00,.8921603E+00,.8956762E+00,
$.8991369E+00,.9025424E+00,.9058923E+00,.9091864E+00,.9124246E+00,
$.9156067E+00,.9187324E+00,.9218016E+00,.9248141E+00,.9277697E+00,
$.9306682E+00,.9335094E+00,.9362932E+00,.9390194E+00,.9416878E+00/
data (xp(i),i= 641, 680)/
$.9442982E+00,.9468506E+00,.9493447E+00,.9517803E+00,.9541574E+00,
$.9564759E+00,.9587354E+00,.9609360E+00,.9630774E+00,.9651596E+00,
$.9671823E+00,.9691456E+00,.9710493E+00,.9728932E+00,.9746772E+00,
$.9764012E+00,.9780652E+00,.9796690E+00,.9812125E+00,.9826957E+00,

```

```

$.9841183E+00,.9854805E+00,.9867820E+00,.9880227E+00,.9892027E+00,
$.9903218E+00,.9913800E+00,.9923771E+00,.9933133E+00,.9941882E+00,
$.9950021E+00,.9957547E+00,.9964460E+00,.9970760E+00,.9976447E+00,
$.9981519E+00,.9985978E+00,.9989822E+00,.9993052E+00,.9995666E+00/
data (xp(i),i= 681, 720)/
$.9997666E+00,.9999050E+00,.9999820E+00,.2615810E-02,.7847359E-02,
$.1307869E-01,.1830967E-01,.2354014E-01,.2876997E-01,.3399902E-01,
$.3922713E-01,.4445417E-01,.4967999E-01,.5490445E-01,.6012741E-01,
$.6534873E-01,.7056825E-01,.7578585E-01,.8100137E-01,.8621467E-01,
$.9142561E-01,.9663405E-01,.1018398E+00,.1070429E+00,.1122429E+00,
$.1174399E+00,.1226337E+00,.1278242E+00,.1330111E+00,.1381944E+00,
$.1433739E+00,.1485495E+00,.1537210E+00,.1588884E+00,.1640513E+00,
$.1692098E+00,.1743636E+00,.1795127E+00,.1846569E+00,.1897960E+00/
data (xp(i),i= 721, 760)/
$.1949299E+00,.2000585E+00,.2051816E+00,.2102991E+00,.2154108E+00,
$.2205166E+00,.2256164E+00,.2307101E+00,.2357974E+00,.2408782E+00,
$.2459525E+00,.2510200E+00,.2560807E+00,.2611343E+00,.2661808E+00,
$.2712201E+00,.2762519E+00,.2812761E+00,.2862926E+00,.2913013E+00,
$.2963021E+00,.3012947E+00,.3062790E+00,.3112550E+00,.3162225E+00,
$.3211813E+00,.3261313E+00,.3310724E+00,.3360045E+00,.3409273E+00,
$.3458408E+00,.3507449E+00,.3556393E+00,.3605240E+00,.3653989E+00,
$.3702637E+00,.3751184E+00,.3799629E+00,.3847969E+00,.3896204E+00/
data (xp(i),i= 761, 800)/
$.3944333E+00,.3992353E+00,.4040264E+00,.4088065E+00,.4135754E+00,
$.4183329E+00,.4230790E+00,.4278136E+00,.4325364E+00,.4372474E+00,
$.4419464E+00,.4466333E+00,.4513080E+00,.4559703E+00,.4606202E+00,
$.4652574E+00,.4698819E+00,.4744936E+00,.4790923E+00,.4836778E+00,
$.4882502E+00,.4928091E+00,.4973546E+00,.5018864E+00,.5064046E+00,
$.5109088E+00,.5153991E+00,.5198753E+00,.5243372E+00,.5287848E+00,
$.5332179E+00,.5376364E+00,.5420402E+00,.5464292E+00,.5508032E+00,
$.5551622E+00,.5595059E+00,.5638343E+00,.5681473E+00,.5724448E+00/
data (xp(i),i= 801, 840)/
$.5767266E+00,.5809926E+00,.5852427E+00,.5894768E+00,.5936947E+00,
$.5978964E+00,.6020817E+00,.6062506E+00,.6104028E+00,.6145384E+00,
$.6186571E+00,.6227589E+00,.6268437E+00,.6309113E+00,.6349616E+00,
$.6389945E+00,.6430100E+00,.6470079E+00,.6509880E+00,.6549504E+00,
$.6588948E+00,.6628211E+00,.6667294E+00,.6706194E+00,.6744910E+00,
$.6783442E+00,.6821788E+00,.6859947E+00,.6897919E+00,.6935702E+00,
$.6973295E+00,.7010697E+00,.7047907E+00,.7084924E+00,.7121748E+00,
$.7158376E+00,.7194809E+00,.7231044E+00,.7267082E+00,.7302921E+00/
data (xp(i),i= 841, 880)/
$.7338560E+00,.7373998E+00,.7409234E+00,.7444268E+00,.7479097E+00,
$.7513722E+00,.7548142E+00,.7582355E+00,.7616360E+00,.7650157E+00,
$.7683744E+00,.7717121E+00,.7750287E+00,.7783241E+00,.7815982E+00,
$.7848508E+00,.7880821E+00,.7912917E+00,.7944797E+00,.7976459E+00,
$.8007903E+00,.8039128E+00,.8070132E+00,.8100916E+00,.8131479E+00,
$.8161818E+00,.8191934E+00,.8221826E+00,.8251493E+00,.8280935E+00,
$.8310149E+00,.8339136E+00,.8367895E+00,.8396425E+00,.8424725E+00,
$.8452794E+00,.8480632E+00,.8508238E+00,.8535611E+00,.8562750E+00/
data (xp(i),i= 881, 920)/
$.8589656E+00,.8616325E+00,.8642760E+00,.8668957E+00,.8694918E+00,
$.8720640E+00,.8746124E+00,.8771368E+00,.8796372E+00,.8821136E+00,
$.8845658E+00,.8869937E+00,.8893975E+00,.8917768E+00,.8941318E+00,
$.8964623E+00,.8987683E+00,.9010496E+00,.9033063E+00,.9055383E+00,
$.9077455E+00,.9099278E+00,.9120852E+00,.9142177E+00,.9163252E+00,
$.9184075E+00,.9204648E+00,.9224968E+00,.9245036E+00,.9264851E+00,
$.9284412E+00,.9303720E+00,.9322772E+00,.9341570E+00,.9360111E+00,
$.9378397E+00,.9396426E+00,.9414198E+00,.9431712E+00,.9448967E+00/
data (xp(i),i= 921, 960)/
$.9465965E+00,.9482703E+00,.9499181E+00,.9515400E+00,.9531358E+00,

```



```

$.9547056E+00,.9562492E+00,.9577666E+00,.9592578E+00,.9607228E+00,
$.9621615E+00,.9635738E+00,.9649597E+00,.9663193E+00,.9676524E+00,
$.9689590E+00,.9702391E+00,.9714927E+00,.9727196E+00,.9739199E+00,
$.9750936E+00,.9762406E+00,.9773609E+00,.9784544E+00,.9795211E+00,
$.9805610E+00,.9815741E+00,.9825603E+00,.9835197E+00,.9844521E+00,
$.9853575E+00,.9862360E+00,.9870876E+00,.9879120E+00,.9887095E+00,
$.9894799E+00,.9902232E+00,.9909394E+00,.9916285E+00,.9922904E+00/
data (xp(i),i= 961,1000)/
$.9929252E+00,.9935328E+00,.9941132E+00,.9946664E+00,.9951924E+00,
$.9956911E+00,.9961626E+00,.9966068E+00,.9970238E+00,.9974135E+00,
$.9977758E+00,.9981109E+00,.9984186E+00,.9986990E+00,.9989521E+00,
$.9991778E+00,.9993762E+00,.9995472E+00,.9996909E+00,.9998072E+00,
$.9998962E+00,.9999577E+00,.9999920E+00,.1962267E-02,.5886772E-02,
$.9811186E-02,.1373545E-01,.1765950E-01,.2158328E-01,.2550673E-01,
$.2942978E-01,.3335238E-01,.3727447E-01,.4119598E-01,.4511686E-01,
$.4903704E-01,.5295647E-01,.5687508E-01,.6079282E-01,.6470962E-01/
data (xp(i),i=1001,1040)/
$.6862542E-01,.7254017E-01,.7645380E-01,.8036625E-01,.8427746E-01,
$.8818738E-01,.9209594E-01,.9600308E-01,.9990874E-01,.1038129E+00,
$.1077154E+00,.1116162E+00,.1155154E+00,.1194128E+00,.1233083E+00,
$.1272019E+00,.1310936E+00,.1349832E+00,.1388708E+00,.1427562E+00,
$.1466395E+00,.1505204E+00,.1543991E+00,.1582754E+00,.1621492E+00,
$.1660205E+00,.1698893E+00,.1737555E+00,.1776190E+00,.1814798E+00,
$.1853377E+00,.1891928E+00,.1930450E+00,.1968942E+00,.2007404E+00,
$.2045835E+00,.2084235E+00,.2122602E+00,.2160937E+00,.2199238E+00/
data (xp(i),i=1041,1080)/
$.2237505E+00,.2275738E+00,.2313936E+00,.2352099E+00,.2390225E+00,
$.2428314E+00,.2466366E+00,.2504380E+00,.2542355E+00,.2580292E+00,
$.2618188E+00,.2656044E+00,.2693859E+00,.2731633E+00,.2769365E+00,
$.2807054E+00,.2844699E+00,.2882301E+00,.2919859E+00,.2957372E+00,
$.2994839E+00,.3032259E+00,.3069634E+00,.3106961E+00,.3144240E+00,
$.3181470E+00,.3218652E+00,.3255784E+00,.3292866E+00,.3329897E+00,
$.3366877E+00,.3403805E+00,.3440681E+00,.3477503E+00,.3514272E+00,
$.3550987E+00,.3587648E+00,.3624253E+00,.3660802E+00,.3697295E+00/
data (xp(i),i=1081,1120)/
$.3733731E+00,.3770109E+00,.3806429E+00,.3842691E+00,.3878893E+00,
$.3915036E+00,.3951118E+00,.3987140E+00,.4023100E+00,.4058998E+00,
$.4094834E+00,.4130607E+00,.4166316E+00,.4201960E+00,.4237541E+00,
$.4273055E+00,.4308504E+00,.4343887E+00,.4379203E+00,.4414451E+00,
$.4449632E+00,.4484743E+00,.4519786E+00,.4554759E+00,.4589662E+00,
$.4624494E+00,.4659255E+00,.4693945E+00,.4728562E+00,.4763106E+00,
$.4797577E+00,.4831973E+00,.4866296E+00,.4900543E+00,.4934715E+00,
$.4968812E+00,.5002831E+00,.5036774E+00,.5070638E+00,.5104425E+00/
data (xp(i),i=1121,1160)/
$.5138133E+00,.5171762E+00,.5205312E+00,.5238781E+00,.5272169E+00,
$.5305477E+00,.5338702E+00,.5371846E+00,.5404906E+00,.5437884E+00,
$.5470777E+00,.5503587E+00,.5536311E+00,.5568951E+00,.5601504E+00,
$.5633972E+00,.5666352E+00,.5698645E+00,.5730851E+00,.5762968E+00,
$.5794996E+00,.5826936E+00,.5858785E+00,.5890544E+00,.5922213E+00,
$.5953790E+00,.5985276E+00,.6016669E+00,.6047970E+00,.6079177E+00,
$.6110291E+00,.6141311E+00,.6172236E+00,.6203066E+00,.6233801E+00,
$.6264440E+00,.6294982E+00,.6325427E+00,.6355775E+00,.6386024E+00/
data (xp(i),i=1161,1200)/
$.6416176E+00,.6446229E+00,.6476182E+00,.6506036E+00,.6535789E+00,
$.6565442E+00,.6594994E+00,.6624444E+00,.6653792E+00,.6683037E+00,
$.6712180E+00,.6741219E+00,.6770155E+00,.6798986E+00,.6827712E+00,
$.6856333E+00,.6884849E+00,.6913259E+00,.6941562E+00,.6969758E+00,
$.6997847E+00,.7025828E+00,.7053701E+00,.7081465E+00,.7109120E+00,
$.7136666E+00,.7164102E+00,.7191427E+00,.7218642E+00,.7245746E+00,
$.7272737E+00,.7299617E+00,.7326385E+00,.7353039E+00,.7379581E+00,

```

```

$.7406008E+00,.7432322E+00,.7458521E+00,.7484606E+00,.7510575E+00/
data (xp(i),i=1201,1240)/
$.7536428E+00,.7562165E+00,.7587786E+00,.7613290E+00,.7638676E+00,
$.7663945E+00,.7689096E+00,.7714129E+00,.7739043E+00,.7763837E+00,
$.7788512E+00,.7813067E+00,.7837502E+00,.7861816E+00,.7886009E+00,
$.7910080E+00,.7934030E+00,.7957857E+00,.7981562E+00,.8005144E+00,
$.8028602E+00,.8051937E+00,.8075148E+00,.8098234E+00,.8121196E+00,
$.8144033E+00,.8166744E+00,.8189330E+00,.8211789E+00,.8234122E+00,
$.8256328E+00,.8278407E+00,.8300358E+00,.8322182E+00,.8343877E+00,
$.8365444E+00,.8386882E+00,.8408191E+00,.8429370E+00,.8450420E+00/
data (xp(i),i=1241,1280)/
$.8471339E+00,.8492128E+00,.8512786E+00,.8533313E+00,.8553709E+00,
$.8573972E+00,.8594104E+00,.8614104E+00,.8633970E+00,.8653704E+00,
$.8673304E+00,.8692771E+00,.8712104E+00,.8731303E+00,.8750367E+00,
$.8769297E+00,.8788091E+00,.8806750E+00,.8825274E+00,.8843662E+00,
$.8861913E+00,.8880028E+00,.8898006E+00,.8915847E+00,.8933550E+00,
$.8951117E+00,.8968545E+00,.8985835E+00,.9002987E+00,.9020000E+00,
$.9036874E+00,.9053609E+00,.9070204E+00,.9086660E+00,.9102976E+00,
$.9119152E+00,.9135187E+00,.9151081E+00,.9166835E+00,.9182447E+00/
data (xp(i),i=1281,1320)/
$.9197918E+00,.9213247E+00,.9228435E+00,.9243480E+00,.9258383E+00,
$.9273143E+00,.9287760E+00,.9302235E+00,.9316566E+00,.9330754E+00,
$.9344797E+00,.9358697E+00,.9372453E+00,.9386065E+00,.9399532E+00,
$.9412854E+00,.9426031E+00,.9439063E+00,.9451950E+00,.9464691E+00,
$.9477286E+00,.9489735E+00,.9502038E+00,.9514195E+00,.9526205E+00,
$.9538069E+00,.9549785E+00,.9561355E+00,.9572777E+00,.9584052E+00,
$.9595179E+00,.9606159E+00,.9616990E+00,.9627673E+00,.9638208E+00,
$.9648595E+00,.9658833E+00,.9668922E+00,.9678863E+00,.9688654E+00/
data (xp(i),i=1321,1360)/
$.9698296E+00,.9707788E+00,.9717132E+00,.9726325E+00,.9735369E+00,
$.9744262E+00,.9753006E+00,.9761600E+00,.9770043E+00,.9778335E+00,
$.9786477E+00,.9794468E+00,.9802309E+00,.9809998E+00,.9817537E+00,
$.9824924E+00,.9832160E+00,.9839244E+00,.9846177E+00,.9852958E+00,
$.9859587E+00,.9866065E+00,.9872390E+00,.9878564E+00,.9884585E+00,
$.9890455E+00,.9896171E+00,.9901736E+00,.9907148E+00,.9912407E+00,
$.9917514E+00,.9922468E+00,.9927269E+00,.9931917E+00,.9936412E+00,
$.9940754E+00,.9944943E+00,.9948979E+00,.9952862E+00,.9956591E+00/
data (xp(i),i=1361,1400)/
$.9960167E+00,.9963590E+00,.9966859E+00,.9969974E+00,.9972936E+00,
$.9975745E+00,.9978400E+00,.9980901E+00,.9983248E+00,.9985442E+00,
$.9987482E+00,.9989368E+00,.9991100E+00,.9992678E+00,.9994103E+00,
$.9995373E+00,.9996489E+00,.9997452E+00,.9998261E+00,.9998915E+00,
$.9999416E+00,.9999762E+00,.9999955E+00,.1570010E-02,.4710016E-02,
$.7849975E-02,.1098986E-01,.1412963E-01,.1726926E-01,.2040873E-01,
$.2354799E-01,.2668702E-01,.2982579E-01,.3296426E-01,.3610241E-01,
$.3924020E-01,.4237761E-01,.4551459E-01,.4865113E-01,.5178719E-01/
data (xp(i),i=1401,1440)/
$.5492274E-01,.5805775E-01,.6119218E-01,.6432601E-01,.6745921E-01,
$.7059174E-01,.7372358E-01,.7685468E-01,.7998504E-01,.8311460E-01,
$.8624334E-01,.8937123E-01,.9249824E-01,.9562434E-01,.9874950E-01,
$.1018737E+00,.1049969E+00,.1081190E+00,.1112401E+00,.1143601E+00,
$.1174789E+00,.1205966E+00,.1237131E+00,.1268284E+00,.1299424E+00,
$.1330552E+00,.1361666E+00,.1392767E+00,.1423855E+00,.1454928E+00,
$.1485987E+00,.1517031E+00,.1548060E+00,.1579074E+00,.1610073E+00,
$.1641055E+00,.1672022E+00,.1702971E+00,.1733905E+00,.1764821E+00/
data (xp(i),i=1441,1480)/
$.1795719E+00,.1826600E+00,.1857463E+00,.1888308E+00,.1919134E+00,
$.1949941E+00,.1980728E+00,.2011497E+00,.2042245E+00,.2072973E+00,
$.2103681E+00,.2134368E+00,.2165035E+00,.2195679E+00,.2226302E+00,
$.2256904E+00,.2287482E+00,.2318039E+00,.2348572E+00,.2379083E+00,

```

```

$.2409569E+00,.2440033E+00,.2470472E+00,.2500886E+00,.2531276E+00,
$.2561641E+00,.2591981E+00,.2622295E+00,.2652584E+00,.2682846E+00,
$.2713082E+00,.2743291E+00,.2773473E+00,.2803628E+00,.2833755E+00,
$.2863854E+00,.2893925E+00,.2923967E+00,.2953980E+00,.2983965E+00/
data (xp(i),i=1481,1520)/
$.3013920E+00,.3043845E+00,.3073740E+00,.3103605E+00,.3133439E+00,
$.3163243E+00,.3193015E+00,.3222756E+00,.3252465E+00,.3282141E+00,
$.3311786E+00,.3341398E+00,.3370977E+00,.3400522E+00,.3430035E+00,
$.3459513E+00,.3488957E+00,.3518367E+00,.3547742E+00,.3577082E+00,
$.3606387E+00,.3635657E+00,.3664890E+00,.3694087E+00,.3723248E+00,
$.3752373E+00,.3781460E+00,.3810510E+00,.3839522E+00,.3868497E+00,
$.3897433E+00,.3926331E+00,.3955190E+00,.3984010E+00,.4012791E+00,
$.4041533E+00,.4070234E+00,.4098896E+00,.4127517E+00,.4156097E+00/
data (xp(i),i=1521,1560)/
$.4184636E+00,.4213134E+00,.4241591E+00,.4270006E+00,.4298378E+00,
$.4326708E+00,.4354996E+00,.4383241E+00,.4411442E+00,.4439600E+00,
$.4467714E+00,.4495784E+00,.4523810E+00,.4551791E+00,.4579727E+00,
$.4607618E+00,.4635464E+00,.4663264E+00,.4691018E+00,.4718726E+00,
$.4746387E+00,.4774001E+00,.4801569E+00,.4829089E+00,.4856561E+00,
$.4883986E+00,.4911362E+00,.4938690E+00,.4965969E+00,.4993199E+00,
$.5020381E+00,.5047512E+00,.5074594E+00,.5101626E+00,.5128607E+00,
$.5155538E+00,.5182418E+00,.5209247E+00,.5236025E+00,.5262750E+00/
data (xp(i),i=1561,1600)/
$.5289425E+00,.5316046E+00,.5342616E+00,.5369133E+00,.5395597E+00,
$.5422007E+00,.5448365E+00,.5474668E+00,.5500918E+00,.5527113E+00,
$.5553254E+00,.5579340E+00,.5605371E+00,.5631347E+00,.5657267E+00,
$.5683131E+00,.5708940E+00,.5734692E+00,.5760387E+00,.5786026E+00,
$.5811608E+00,.5837132E+00,.5862599E+00,.5888008E+00,.5913360E+00,
$.5938652E+00,.5963886E+00,.5989062E+00,.6014178E+00,.6039235E+00,
$.6064233E+00,.6089170E+00,.6114048E+00,.6138865E+00,.6163622E+00,
$.6188318E+00,.6212953E+00,.6237527E+00,.6262040E+00,.6286490E+00/
data (xp(i),i=1601,1640)/
$.6310879E+00,.6335205E+00,.6359469E+00,.6383670E+00,.6407808E+00,
$.6431883E+00,.6455895E+00,.6479843E+00,.6503727E+00,.6527547E+00,
$.6551303E+00,.6574994E+00,.6598620E+00,.6622181E+00,.6645677E+00,
$.6669107E+00,.6692472E+00,.6715770E+00,.6739003E+00,.6762169E+00,
$.6785268E+00,.6808300E+00,.6831266E+00,.6854163E+00,.6876994E+00,
$.6899756E+00,.6922451E+00,.6945077E+00,.6967635E+00,.6990124E+00,
$.7012544E+00,.7034895E+00,.7057176E+00,.7079388E+00,.7101530E+00,
$.7123603E+00,.7145605E+00,.7167536E+00,.7189397E+00,.7211187E+00/
data (xp(i),i=1641,1680)/
$.7232906E+00,.7254553E+00,.7276129E+00,.7297634E+00,.7319066E+00,
$.7340426E+00,.7361714E+00,.7382929E+00,.7404071E+00,.7425141E+00,
$.7446137E+00,.7467060E+00,.7487909E+00,.7508684E+00,.7529385E+00,
$.7550013E+00,.7570565E+00,.7591043E+00,.7611446E+00,.7631774E+00,
$.7652027E+00,.7672204E+00,.7692306E+00,.7712332E+00,.7732282E+00,
$.7752155E+00,.7771953E+00,.7791673E+00,.7811317E+00,.7830884E+00,
$.7850373E+00,.7869785E+00,.7889120E+00,.7908376E+00,.7927555E+00,
$.7946656E+00,.7965678E+00,.7984622E+00,.8003486E+00,.8022273E+00/
data (xp(i),i=1681,1720)/
$.8040979E+00,.8059607E+00,.8078155E+00,.8096624E+00,.8115013E+00,
$.8133321E+00,.8151550E+00,.8169698E+00,.8187765E+00,.8205752E+00,
$.8223658E+00,.8241483E+00,.8259227E+00,.8276889E+00,.8294470E+00,
$.8311968E+00,.8329385E+00,.8346720E+00,.8363972E+00,.8381142E+00,
$.8398229E+00,.8415234E+00,.8432156E+00,.8448994E+00,.8465749E+00,
$.8482421E+00,.8499009E+00,.8515513E+00,.8531933E+00,.8548269E+00,
$.8564521E+00,.8580688E+00,.8596771E+00,.8612769E+00,.8628682E+00,
$.8644510E+00,.8660253E+00,.8675910E+00,.8691482E+00,.8706968E+00/
data (xp(i),i=1721,1760)/
$.8722369E+00,.8737683E+00,.8752911E+00,.8768053E+00,.8783108E+00,

```

```

$.8798077E+00,.8812959E+00,.8827754E+00,.8842463E+00,.8857083E+00,
$.8871617E+00,.8886063E+00,.8900422E+00,.8914692E+00,.8928875E+00,
$.8942970E+00,.8956977E+00,.8970895E+00,.8984725E+00,.8998466E+00,
$.9012119E+00,.9025683E+00,.9039157E+00,.9052543E+00,.9065839E+00,
$.9079046E+00,.9092164E+00,.9105192E+00,.9118130E+00,.9130978E+00,
$.9143736E+00,.9156404E+00,.9168982E+00,.9181469E+00,.9193866E+00,
$.9206172E+00,.9218387E+00,.9230511E+00,.9242545E+00,.9254487E+00/
data (xp(i),i=1761,1800)/
$.9266338E+00,.9278098E+00,.9289766E+00,.9301343E+00,.9312828E+00,
$.9324221E+00,.9335522E+00,.9346731E+00,.9357848E+00,.9368872E+00,
$.9379805E+00,.9390645E+00,.9401392E+00,.9412046E+00,.9422608E+00,
$.9433077E+00,.9443453E+00,.9453735E+00,.9463925E+00,.9474021E+00,
$.9484024E+00,.9493933E+00,.9503749E+00,.9513471E+00,.9523099E+00,
$.9532633E+00,.9542073E+00,.9551420E+00,.9560672E+00,.9569829E+00,
$.9578893E+00,.9587862E+00,.9596736E+00,.9605516E+00,.9614201E+00,
$.9622791E+00,.9631287E+00,.9639687E+00,.9647992E+00,.9656203E+00/
data (xp(i),i=1801,1840)/
$.9664318E+00,.9672338E+00,.9680262E+00,.9688091E+00,.9695824E+00,
$.9703462E+00,.9711004E+00,.9718451E+00,.9725801E+00,.9733056E+00,
$.9740215E+00,.9747278E+00,.9754244E+00,.9761115E+00,.9767889E+00,
$.9774567E+00,.9781148E+00,.9787633E+00,.9794022E+00,.9800314E+00,
$.9806509E+00,.9812608E+00,.9818610E+00,.9824515E+00,.9830323E+00,
$.9836035E+00,.9841649E+00,.9847166E+00,.9852586E+00,.9857909E+00,
$.9863135E+00,.9868264E+00,.9873295E+00,.9878229E+00,.9883066E+00,
$.9887805E+00,.9892447E+00,.9896991E+00,.9901437E+00,.9905786E+00/
data (xp(i),i=1841,1880)/
$.9910037E+00,.9914191E+00,.9918247E+00,.9922205E+00,.9926065E+00,
$.9929827E+00,.9933492E+00,.9937058E+00,.9940527E+00,.9943897E+00,
$.9947169E+00,.9950344E+00,.9953420E+00,.9956398E+00,.9959278E+00,
$.9962060E+00,.9964743E+00,.9967328E+00,.9969815E+00,.9972204E+00,
$.9974494E+00,.9976686E+00,.9978780E+00,.9980775E+00,.9982672E+00,
$.9984471E+00,.9986171E+00,.9987772E+00,.9989275E+00,.9990680E+00,
$.9991986E+00,.9993193E+00,.9994302E+00,.9995313E+00,.9996225E+00,
$.9997038E+00,.9997753E+00,.9998369E+00,.9998886E+00,.9999306E+00/
data (xp(i),i=1881,1883)/
$.9999626E+00,.9999848E+00,.9999971E+00/
data (wp(i),i= 1, 40)/
$.4679139E+00,.3607616E+00,.1713245E+00,.2955242E+00,.2692667E+00,
$.2190864E+00,.1494513E+00,.6667134E-01,.1527534E+00,.1491730E+00,
$.1420961E+00,.1316886E+00,.1181945E+00,.1019301E+00,.8327674E-01,
$.6267205E-01,.4060143E-01,.1761401E-01,.7750595E-01,.7703982E-01,
$.7611036E-01,.7472317E-01,.7288658E-01,.7061165E-01,.6791205E-01,
$.6480401E-01,.6130624E-01,.5743977E-01,.5322785E-01,.4869581E-01,
$.4387091E-01,.3878217E-01,.3346020E-01,.2793701E-01,.2224585E-01,
$.1642106E-01,.1049828E-01,.4521277E-02,.5190788E-01,.5176794E-01/
data (wp(i),i= 41, 80)/
$.5148845E-01,.5107016E-01,.5051418E-01,.4982204E-01,.4899558E-01,
$.4803703E-01,.4694899E-01,.4573438E-01,.4439648E-01,.4293889E-01,
$.4136555E-01,.3968070E-01,.3788887E-01,.3599490E-01,.3400389E-01,
$.3192122E-01,.2975249E-01,.2750356E-01,.2518048E-01,.2278952E-01,
$.2033712E-01,.1782990E-01,.1527462E-01,.1267817E-01,.1004756E-01,
$.7389931E-02,.4712730E-02,.2026812E-02,.3901781E-01,.3895840E-01,
$.3883965E-01,.3866176E-01,.3842499E-01,.3812971E-01,.3777636E-01,
$.3736549E-01,.3689771E-01,.3637375E-01,.3579439E-01,.3516053E-01/
data (wp(i),i= 81, 120)/
$.3447312E-01,.3373321E-01,.3294194E-01,.3210050E-01,.3121017E-01,
$.3027232E-01,.2928837E-01,.2825982E-01,.2718823E-01,.2607524E-01,
$.2492254E-01,.2373188E-01,.2250509E-01,.2124403E-01,.1995061E-01,
$.1862681E-01,.1727465E-01,.1589618E-01,.1449351E-01,.1306876E-01,
$.1162411E-01,.1016177E-01,.8683945E-02,.7192905E-02,.5690922E-02,

```

```

$.4180313E-02,.2663534E-02,.1144950E-02,.3125542E-01,.3122488E-01,
$.3116384E-01,.3107234E-01,.3095048E-01,.3079838E-01,.3061619E-01,
$.3040408E-01,.3016227E-01,.2989098E-01,.2959049E-01,.2926108E-01/
data (wp(i),i= 121, 160)/
$.2890309E-01,.2851685E-01,.2810276E-01,.2766120E-01,.2719261E-01,
$.2669746E-01,.2617622E-01,.2562940E-01,.2505754E-01,.2446120E-01,
$.2384096E-01,.2319742E-01,.2253122E-01,.2184300E-01,.2113344E-01,
$.2040323E-01,.1965309E-01,.1888374E-01,.1809594E-01,.1729046E-01,
$.1646809E-01,.1562962E-01,.1477588E-01,.1390771E-01,.1302595E-01,
$.1213146E-01,.1122511E-01,.1030780E-01,.9380420E-02,.8443871E-02,
$.7499073E-02,.6546948E-02,.5588428E-02,.4624450E-02,.3655961E-02,
$.2683925E-02,.1709393E-02,.7346345E-03,.2087312E-01,.2086402E-01/
data (wp(i),i= 161, 200)/
$.2084584E-01,.2081857E-01,.2078223E-01,.2073683E-01,.2068240E-01,
$.2061896E-01,.2054653E-01,.2046515E-01,.2037486E-01,.2027568E-01,
$.2016767E-01,.2005088E-01,.1992534E-01,.1979113E-01,.1964829E-01,
$.1949689E-01,.1933700E-01,.1916867E-01,.1899200E-01,.1880705E-01,
$.1861391E-01,.1841266E-01,.1820338E-01,.1798617E-01,.1776113E-01,
$.1752835E-01,.1728792E-01,.1703997E-01,.1678459E-01,.1652190E-01,
$.1625201E-01,.1597504E-01,.1569110E-01,.1540033E-01,.1510285E-01,
$.1479879E-01,.1448828E-01,.1417146E-01,.1384846E-01,.1351943E-01/
data (wp(i),i= 201, 240)/
$.1318451E-01,.1284384E-01,.1249758E-01,.1214587E-01,.1178887E-01,
$.1142673E-01,.1105962E-01,.1068768E-01,.1031109E-01,.9930004E-02,
$.9544593E-02,.9155022E-02,.8761463E-02,.8364086E-02,.7963064E-02,
$.7558573E-02,.7150788E-02,.6739888E-02,.6326051E-02,.5909457E-02,
$.5490289E-02,.5068728E-02,.4644959E-02,.4219166E-02,.3791535E-02,
$.3362252E-02,.2931504E-02,.2499479E-02,.2066366E-02,.1632357E-02,
$.1197647E-02,.7624721E-03,.3276087E-03,.1566826E-01,.1566442E-01,
$.1565672E-01,.1564519E-01,.1562981E-01,.1561059E-01,.1558755E-01/
data (wp(i),i= 241, 280)/
$.1556067E-01,.1552998E-01,.1549547E-01,.1545716E-01,.1541506E-01,
$.1536917E-01,.1531950E-01,.1526608E-01,.1520891E-01,.1514800E-01,
$.1508338E-01,.1501505E-01,.1494303E-01,.1486735E-01,.1478802E-01,
$.1470505E-01,.1461848E-01,.1452832E-01,.1443459E-01,.1433731E-01,
$.1423652E-01,.1413223E-01,.1402447E-01,.1391327E-01,.1379866E-01,
$.1368065E-01,.1355929E-01,.1343460E-01,.1330661E-01,.1317535E-01,
$.1304086E-01,.1290316E-01,.1276230E-01,.1261831E-01,.1247122E-01,
$.1232106E-01,.1216788E-01,.1201172E-01,.1185260E-01,.1169058E-01/
data (wp(i),i= 281, 320)/
$.1152568E-01,.1135796E-01,.1118744E-01,.1101418E-01,.1083822E-01,
$.1065959E-01,.1047835E-01,.1029454E-01,.1010820E-01,.9919373E-02,
$.9728115E-02,.9534468E-02,.9338480E-02,.9140200E-02,.8939676E-02,
$.8736957E-02,.8532093E-02,.8325134E-02,.8116132E-02,.7905137E-02,
$.7692201E-02,.7477377E-02,.7260717E-02,.7042274E-02,.6822103E-02,
$.6600256E-02,.6376790E-02,.6151757E-02,.5925215E-02,.5697218E-02,
$.5467822E-02,.5237083E-02,.5005059E-02,.4771806E-02,.4537382E-02,
$.4301844E-02,.4065249E-02,.3827657E-02,.3589125E-02,.3349711E-02/
data (wp(i),i= 321, 360)/
$.3109476E-02,.2868477E-02,.2626773E-02,.2384425E-02,.2141492E-02,
$.1898033E-02,.1654108E-02,.1409777E-02,.1165101E-02,.9201405E-03,
$.6749606E-03,.4296466E-03,.1845901E-03,.1045439E-01,.1045325E-01,
$.1045097E-01,.1044754E-01,.1044297E-01,.1043726E-01,.1043041E-01,
$.1042242E-01,.1041329E-01,.1040302E-01,.1039161E-01,.1037907E-01,
$.1036539E-01,.1035058E-01,.1033464E-01,.1031758E-01,.1029938E-01,
$.1028006E-01,.1025961E-01,.1023804E-01,.1021535E-01,.1019155E-01,
$.1016663E-01,.1014060E-01,.1011347E-01,.1008523E-01,.1005588E-01/
data (wp(i),i= 361, 400)/
$.1002544E-01,.9993899E-02,.9961267E-02,.9927547E-02,.9892741E-02,
$.9856855E-02,.9819891E-02,.9781854E-02,.9742747E-02,.9702576E-02,

```

```

$.9661345E-02,.9619057E-02,.9575718E-02,.9531333E-02,.9485905E-02,
$.9439441E-02,.9391946E-02,.9343424E-02,.9293880E-02,.9243321E-02,
$.9191751E-02,.9139177E-02,.9085604E-02,.9031038E-02,.8975485E-02,
$.8918951E-02,.8861442E-02,.8802965E-02,.8743525E-02,.8683130E-02,
$.8621786E-02,.8559499E-02,.8496277E-02,.8432127E-02,.8367054E-02,
$.8301068E-02,.8234174E-02,.8166380E-02,.8097693E-02,.8028121E-02/
data (wp(i),i= 401, 440)/
$.7957672E-02,.7886353E-02,.7814173E-02,.7741138E-02,.7667257E-02,
$.7592538E-02,.7516989E-02,.7440619E-02,.7363435E-02,.7285447E-02,
$.7206662E-02,.7127090E-02,.7046739E-02,.6965617E-02,.6883734E-02,
$.6801099E-02,.6717721E-02,.6633608E-02,.6548770E-02,.6463217E-02,
$.6376957E-02,.6290000E-02,.6202356E-02,.6114033E-02,.6025043E-02,
$.5935394E-02,.5845096E-02,.5754159E-02,.5662594E-02,.5570409E-02,
$.5477616E-02,.5384224E-02,.5290244E-02,.5195685E-02,.5100559E-02,
$.5004875E-02,.4908644E-02,.4811876E-02,.4714583E-02,.4616774E-02/
data (wp(i),i= 441, 480)/
$.4518461E-02,.4419654E-02,.4320364E-02,.4220601E-02,.4120378E-02,
$.4019704E-02,.3918590E-02,.3817049E-02,.3715090E-02,.3612725E-02,
$.3509965E-02,.3406822E-02,.3303306E-02,.3199429E-02,.3095203E-02,
$.2990638E-02,.2885746E-02,.2780539E-02,.2675029E-02,.2569225E-02,
$.2463141E-02,.2356788E-02,.2250177E-02,.2143320E-02,.2036229E-02,
$.1928915E-02,.1821391E-02,.1713667E-02,.1605756E-02,.1497670E-02,
$.1389420E-02,.1281018E-02,.1172476E-02,.1063806E-02,.9550200E-03,
$.8461294E-03,.7371464E-03,.6280830E-03,.5189512E-03,.4097636E-03/
data (wp(i),i= 481, 520)/
$.3005340E-03,.1912855E-03,.8217779E-04,.7844110E-02,.7843627E-02,
$.7842662E-02,.7841214E-02,.7839284E-02,.7836871E-02,.7833976E-02,
$.7830599E-02,.7826741E-02,.7822400E-02,.7817579E-02,.7812276E-02,
$.7806493E-02,.7800229E-02,.7793485E-02,.7786262E-02,.7778560E-02,
$.7770379E-02,.7761720E-02,.7752583E-02,.7742970E-02,.7732880E-02,
$.7722314E-02,.7711273E-02,.7699757E-02,.7687768E-02,.7675306E-02,
$.7662371E-02,.7648965E-02,.7635088E-02,.7620742E-02,.7605926E-02,
$.7590643E-02,.7574892E-02,.7558676E-02,.7541994E-02,.7524848E-02/
data (wp(i),i= 521, 560)/
$.7507240E-02,.7489169E-02,.7470638E-02,.7451646E-02,.7432197E-02,
$.7412290E-02,.7391927E-02,.7371109E-02,.7349838E-02,.7328114E-02,
$.7305939E-02,.7283315E-02,.7260243E-02,.7236724E-02,.7212760E-02,
$.7188352E-02,.7163501E-02,.7138210E-02,.7112480E-02,.7086312E-02,
$.7059708E-02,.7032669E-02,.7005198E-02,.6977296E-02,.6948964E-02,
$.6920205E-02,.6891020E-02,.6861412E-02,.6831380E-02,.6800929E-02,
$.6770059E-02,.6738773E-02,.6707072E-02,.6674958E-02,.6642433E-02,
$.6609500E-02,.6576160E-02,.6542416E-02,.6508269E-02,.6473721E-02/
data (wp(i),i= 561, 600)/
$.6438775E-02,.6403433E-02,.6367697E-02,.6331569E-02,.6295052E-02,
$.6258147E-02,.6220857E-02,.6183184E-02,.6145131E-02,.6106700E-02,
$.6067893E-02,.6028713E-02,.5989161E-02,.5949241E-02,.5908956E-02,
$.5868306E-02,.5827296E-02,.5785926E-02,.5744201E-02,.5702123E-02,
$.5659693E-02,.5616916E-02,.5573792E-02,.5530326E-02,.5486520E-02,
$.5442376E-02,.5397897E-02,.5353085E-02,.5307945E-02,.5262478E-02,
$.5216687E-02,.5170575E-02,.5124145E-02,.5077400E-02,.5030342E-02,
$.4982975E-02,.4935301E-02,.4887323E-02,.4839045E-02,.4790469E-02/
data (wp(i),i= 601, 640)/
$.4741598E-02,.4692436E-02,.4642984E-02,.4593247E-02,.4543228E-02,
$.4492929E-02,.4442353E-02,.4391504E-02,.4340385E-02,.4288999E-02,
$.4237349E-02,.4185438E-02,.4133270E-02,.4080847E-02,.4028173E-02,
$.3975251E-02,.3922085E-02,.3868678E-02,.3815032E-02,.3761152E-02,
$.3707040E-02,.3652700E-02,.3598135E-02,.3543349E-02,.3488345E-02,
$.3433126E-02,.3377697E-02,.3322059E-02,.3266217E-02,.3210173E-02,
$.3153933E-02,.3097498E-02,.3040873E-02,.2984060E-02,.2927064E-02,
$.2869888E-02,.2812535E-02,.2755009E-02,.2697314E-02,.2639453E-02/

```

```

data (wp(i),i= 641, 680)/
$.2581429E-02,.2523247E-02,.2464909E-02,.2406419E-02,.2347782E-02,
$.2289000E-02,.2230077E-02,.2171017E-02,.2111823E-02,.2052500E-02,
$.1993050E-02,.1933477E-02,.1873786E-02,.1813979E-02,.1754061E-02,
$.1694034E-02,.1633904E-02,.1573673E-02,.1513345E-02,.1452924E-02,
$.1392413E-02,.1331817E-02,.1271139E-02,.1210383E-02,.1149552E-02,
$.1088651E-02,.1027682E-02,.9666507E-03,.9055595E-03,.8444126E-03,
$.7832138E-03,.7219667E-03,.6606753E-03,.5993432E-03,.5379742E-03,
$.4765722E-03,.4151409E-03,.3536841E-03,.2922057E-03,.2307099E-03/
data (wp(i),i= 681, 720)/
$.1692014E-03,.1076904E-03,.4626372E-04,.5231608E-02,.5231465E-02,
$.5231179E-02,.5230749E-02,.5230177E-02,.5229461E-02,.5228602E-02,
$.5227600E-02,.5226454E-02,.5225166E-02,.5223735E-02,.5222161E-02,
$.5220444E-02,.5218584E-02,.5216581E-02,.5214435E-02,.5212147E-02,
$.5209716E-02,.5207142E-02,.5204426E-02,.5201567E-02,.5198567E-02,
$.5195423E-02,.5192138E-02,.5188710E-02,.5185141E-02,.5181429E-02,
$.5177576E-02,.5173581E-02,.5169445E-02,.5165167E-02,.5160747E-02,
$.5156186E-02,.5151485E-02,.5146642E-02,.5141658E-02,.5136534E-02/
data (wp(i),i= 721, 760)/
$.5131269E-02,.5125863E-02,.5120318E-02,.5114632E-02,.5108806E-02,
$.5102840E-02,.5096735E-02,.5090490E-02,.5084106E-02,.5077583E-02,
$.5070920E-02,.5064119E-02,.5057180E-02,.5050102E-02,.5042885E-02,
$.5035531E-02,.5028039E-02,.5020409E-02,.5012642E-02,.5004738E-02,
$.4996696E-02,.4988518E-02,.4980203E-02,.4971752E-02,.4963165E-02,
$.4954443E-02,.4945584E-02,.4936590E-02,.4927461E-02,.4918197E-02,
$.4908799E-02,.4899266E-02,.4889599E-02,.4879799E-02,.4869864E-02,
$.4859797E-02,.4849596E-02,.4839263E-02,.4828797E-02,.4818199E-02/
data (wp(i),i= 761, 800)/
$.4807470E-02,.4796608E-02,.4785616E-02,.4774492E-02,.4763238E-02,
$.4751853E-02,.4740338E-02,.4728694E-02,.4716920E-02,.4705017E-02,
$.4692985E-02,.4680825E-02,.4668537E-02,.4656121E-02,.4643577E-02,
$.4630907E-02,.4618109E-02,.4605185E-02,.4592136E-02,.4578960E-02,
$.4565659E-02,.4552233E-02,.4538683E-02,.4525008E-02,.4511210E-02,
$.4497288E-02,.4483243E-02,.4469075E-02,.4454785E-02,.4440373E-02,
$.4425840E-02,.4411185E-02,.4396410E-02,.4381514E-02,.4366498E-02,
$.4351363E-02,.4336109E-02,.4320736E-02,.4305245E-02,.4289636E-02/
data (wp(i),i= 801, 840)/
$.4273910E-02,.4258066E-02,.4242106E-02,.4226030E-02,.4209839E-02,
$.4193532E-02,.4177110E-02,.4160574E-02,.4143924E-02,.4127161E-02,
$.4110284E-02,.4093296E-02,.4076195E-02,.4058982E-02,.4041659E-02,
$.4024225E-02,.4006681E-02,.3989027E-02,.3971264E-02,.3953392E-02,
$.3935412E-02,.3917324E-02,.3899129E-02,.3880828E-02,.3862420E-02,
$.3843906E-02,.3825288E-02,.3806564E-02,.3787737E-02,.3768805E-02,
$.3749771E-02,.3730634E-02,.3711395E-02,.3692054E-02,.3672612E-02,
$.3653070E-02,.3633427E-02,.3613685E-02,.3593845E-02,.3573906E-02/
data (wp(i),i= 841, 880)/
$.3553869E-02,.3533735E-02,.3513504E-02,.3493177E-02,.3472754E-02,
$.3452237E-02,.3431624E-02,.3410918E-02,.3390119E-02,.3369227E-02,
$.3348242E-02,.3327166E-02,.3305999E-02,.3284741E-02,.3263394E-02,
$.3241957E-02,.3220431E-02,.3198818E-02,.3177116E-02,.3155328E-02,
$.3133454E-02,.3111493E-02,.3089448E-02,.3067318E-02,.3045104E-02,
$.3022806E-02,.3000426E-02,.2977964E-02,.2955420E-02,.2932796E-02,
$.2910091E-02,.2887306E-02,.2864443E-02,.2841501E-02,.2818481E-02,
$.2795384E-02,.2772211E-02,.2748961E-02,.2725637E-02,.2702238E-02/
data (wp(i),i= 881, 920)/
$.2678765E-02,.2655218E-02,.2631599E-02,.2607908E-02,.2584146E-02,
$.2560312E-02,.2536409E-02,.2512437E-02,.2488395E-02,.2464286E-02,
$.2440109E-02,.2415865E-02,.2391555E-02,.2367179E-02,.2342739E-02,
$.2318235E-02,.2293667E-02,.2269037E-02,.2244344E-02,.2219590E-02,
$.2194775E-02,.2169901E-02,.2144966E-02,.2119973E-02,.2094922E-02,

```

```

$.2069814E-02,.2044649E-02,.2019428E-02,.1994152E-02,.1968821E-02,
$.1943437E-02,.1917999E-02,.1892508E-02,.1866966E-02,.1841373E-02,
$.1815729E-02,.1790036E-02,.1764294E-02,.1738503E-02,.1712665E-02/
data (wp(i),i= 921, 960)/
$.1686780E-02,.1660848E-02,.1634872E-02,.1608850E-02,.1582785E-02,
$.1556676E-02,.1530525E-02,.1504331E-02,.1478097E-02,.1451822E-02,
$.1425507E-02,.1399154E-02,.1372762E-02,.1346332E-02,.1319866E-02,
$.1293363E-02,.1266825E-02,.1240253E-02,.1213646E-02,.1187006E-02,
$.1160334E-02,.1133630E-02,.1106895E-02,.1080130E-02,.1053335E-02,
$.1026511E-02,.9996593E-03,.9727801E-03,.9458743E-03,.9189426E-03,
$.8919858E-03,.8650045E-03,.8379996E-03,.8109717E-03,.7839217E-03,
$.7568502E-03,.7297579E-03,.7026457E-03,.6755143E-03,.6483644E-03/
data (wp(i),i= 961,1000)/
$.6211967E-03,.5940120E-03,.5668111E-03,.5395947E-03,.5123635E-03,
$.4851182E-03,.4578597E-03,.4305887E-03,.4033058E-03,.3760120E-03,
$.3487078E-03,.3213941E-03,.2940716E-03,.2667411E-03,.2394033E-03,
$.2120589E-03,.1847087E-03,.1573535E-03,.1299941E-03,.1026314E-03,
$.7526651E-04,.4790311E-04,.2057885E-04,.3924530E-02,.3924469E-02,
$.3924348E-02,.3924167E-02,.3923925E-02,.3923623E-02,.3923260E-02,
$.3922837E-02,.3922354E-02,.3921810E-02,.3921206E-02,.3920541E-02,
$.3919816E-02,.3919030E-02,.3918185E-02,.3917278E-02,.3916312E-02/
data (wp(i),i=1001,1040)/
$.3915285E-02,.3914198E-02,.3913051E-02,.3911843E-02,.3910575E-02,
$.3909247E-02,.3907858E-02,.3906410E-02,.3904901E-02,.3903332E-02,
$.3901703E-02,.3900014E-02,.3898265E-02,.3896456E-02,.3894587E-02,
$.3892658E-02,.3890668E-02,.3888619E-02,.3886510E-02,.3884342E-02,
$.3882113E-02,.3879825E-02,.3877476E-02,.3875068E-02,.3872601E-02,
$.3870074E-02,.3867487E-02,.3864840E-02,.3862134E-02,.3859369E-02,
$.3856544E-02,.3853660E-02,.3850716E-02,.3847713E-02,.3844651E-02,
$.3841530E-02,.3838349E-02,.3835109E-02,.3831811E-02,.3828453E-02/
data (wp(i),i=1041,1080)/
$.3825036E-02,.3821561E-02,.3818026E-02,.3814433E-02,.3810781E-02,
$.3807070E-02,.3803300E-02,.3799473E-02,.3795586E-02,.3791641E-02,
$.3787638E-02,.3783576E-02,.3779456E-02,.3775278E-02,.3771042E-02,
$.3766747E-02,.3762395E-02,.3757984E-02,.3753516E-02,.3748990E-02,
$.3744406E-02,.3739765E-02,.3735066E-02,.3730309E-02,.3725495E-02,
$.3720624E-02,.3715695E-02,.3710709E-02,.3705666E-02,.3700566E-02,
$.3695408E-02,.3690194E-02,.3684923E-02,.3679596E-02,.3674211E-02,
$.3668770E-02,.3663273E-02,.3657719E-02,.3652109E-02,.3646442E-02/
data (wp(i),i=1081,1120)/
$.3640720E-02,.3634941E-02,.3629106E-02,.3623216E-02,.3617269E-02,
$.3611267E-02,.3605209E-02,.3599096E-02,.3592927E-02,.3586703E-02,
$.3580424E-02,.3574090E-02,.3567700E-02,.3561256E-02,.3554757E-02,
$.3548203E-02,.3541594E-02,.3534931E-02,.3528213E-02,.3521441E-02,
$.3514615E-02,.3507734E-02,.3500800E-02,.3493812E-02,.3486770E-02,
$.3479674E-02,.3472524E-02,.3465321E-02,.3458065E-02,.3450756E-02,
$.3443393E-02,.3435977E-02,.3428508E-02,.3420987E-02,.3413413E-02,
$.3405786E-02,.3398107E-02,.3390375E-02,.3382592E-02,.3374756E-02/
data (wp(i),i=1121,1160)/
$.3366868E-02,.3358928E-02,.3350937E-02,.3342894E-02,.3334800E-02,
$.3326654E-02,.3318457E-02,.3310208E-02,.3301909E-02,.3293559E-02,
$.3285158E-02,.3276707E-02,.3268205E-02,.3259653E-02,.3251051E-02,
$.3242398E-02,.3233696E-02,.3224944E-02,.3216142E-02,.3207290E-02,
$.3198390E-02,.3189440E-02,.3180440E-02,.3171392E-02,.3162295E-02,
$.3153149E-02,.3143955E-02,.3134712E-02,.3125421E-02,.3116082E-02,
$.3106695E-02,.3097260E-02,.3087778E-02,.3078247E-02,.3068670E-02,
$.3059045E-02,.3049373E-02,.3039654E-02,.3029888E-02,.3020075E-02/
data (wp(i),i=1161,1200)/
$.3010217E-02,.3000311E-02,.2990360E-02,.2980362E-02,.2970318E-02,
$.2960229E-02,.2950094E-02,.2939914E-02,.2929688E-02,.2919418E-02,

```



```

$.2909102E-02,.2898742E-02,.2888336E-02,.2877887E-02,.2867393E-02,
$.2856855E-02,.2846273E-02,.2835647E-02,.2824977E-02,.2814264E-02,
$.2803508E-02,.2792708E-02,.2781865E-02,.2770980E-02,.2760052E-02,
$.2749081E-02,.2738068E-02,.2727013E-02,.2715915E-02,.2704776E-02,
$.2693596E-02,.2682374E-02,.2671110E-02,.2659805E-02,.2648460E-02,
$.2637073E-02,.2625646E-02,.2614179E-02,.2602671E-02,.2591123E-02/
data (wp(i),i=1201,1240)/
$.2579536E-02,.2567908E-02,.2556241E-02,.2544535E-02,.2532789E-02,
$.2521005E-02,.2509181E-02,.2497319E-02,.2485419E-02,.2473480E-02,
$.2461503E-02,.2449488E-02,.2437436E-02,.2425346E-02,.2413218E-02,
$.2401054E-02,.2388852E-02,.2376614E-02,.2364339E-02,.2352028E-02,
$.2339680E-02,.2327296E-02,.2314877E-02,.2302422E-02,.2289931E-02,
$.2277405E-02,.2264844E-02,.2252249E-02,.2239618E-02,.2226953E-02,
$.2214254E-02,.2201520E-02,.2188753E-02,.2175952E-02,.2163117E-02,
$.2150249E-02,.2137349E-02,.2124415E-02,.2111448E-02,.2098449E-02/
data (wp(i),i=1241,1280)/
$.2085417E-02,.2072354E-02,.2059258E-02,.2046131E-02,.2032972E-02,
$.2019782E-02,.2006561E-02,.1993309E-02,.1980026E-02,.1966713E-02,
$.1953370E-02,.1939996E-02,.1926592E-02,.1913159E-02,.1899697E-02,
$.1886205E-02,.1872684E-02,.1859134E-02,.1845555E-02,.1831949E-02,
$.1818314E-02,.1804650E-02,.1790960E-02,.1777241E-02,.1763495E-02,
$.1749722E-02,.1735922E-02,.1722096E-02,.1708242E-02,.1694363E-02,
$.1680457E-02,.1666526E-02,.1652569E-02,.1638586E-02,.1624578E-02,
$.1610545E-02,.1596488E-02,.1582405E-02,.1568299E-02,.1554168E-02/
data (wp(i),i=1281,1320)/
$.1540013E-02,.1525835E-02,.1511633E-02,.1497407E-02,.1483159E-02,
$.1468888E-02,.1454594E-02,.1440278E-02,.1425940E-02,.1411579E-02,
$.1397197E-02,.1382794E-02,.1368369E-02,.1353923E-02,.1339456E-02,
$.1324969E-02,.1310461E-02,.1295933E-02,.1281385E-02,.1266817E-02,
$.1252230E-02,.1237623E-02,.1222998E-02,.1208353E-02,.1193690E-02,
$.1179009E-02,.1164309E-02,.1149592E-02,.1134857E-02,.1120104E-02,
$.1105334E-02,.1090547E-02,.1075743E-02,.1060923E-02,.1046086E-02,
$.1031234E-02,.1016365E-02,.1001481E-02,.9865808E-03,.9716658E-03/
data (wp(i),i=1321,1360)/
$.9567359E-03,.9417913E-03,.9268321E-03,.9118587E-03,.8968712E-03,
$.8818700E-03,.8668551E-03,.8518269E-03,.8367855E-03,.8217313E-03,
$.8066644E-03,.7915851E-03,.7764936E-03,.7613902E-03,.7462750E-03,
$.7311483E-03,.7160104E-03,.7008614E-03,.6857017E-03,.6705313E-03,
$.6553507E-03,.6401600E-03,.6249593E-03,.6097491E-03,.5945295E-03,
$.5793007E-03,.5640630E-03,.5488166E-03,.5335618E-03,.5182987E-03,
$.5030277E-03,.4877489E-03,.4724626E-03,.4571690E-03,.4418684E-03,
$.4265610E-03,.4112470E-03,.3959267E-03,.3806003E-03,.3652680E-03/
data (wp(i),i=1361,1400)/
$.3499301E-03,.3345867E-03,.3192383E-03,.3038849E-03,.2885269E-03,
$.2731644E-03,.2577977E-03,.2424270E-03,.2270526E-03,.2116747E-03,
$.1962935E-03,.1809093E-03,.1655224E-03,.1501328E-03,.1347410E-03,
$.1193471E-03,.1039514E-03,.8855408E-04,.7315545E-04,.5775582E-04,
$.4235569E-04,.2695689E-04,.1158044E-04,.3140018E-02,.3139987E-02,
$.3139926E-02,.3139833E-02,.3139709E-02,.3139554E-02,.3139368E-02,
$.3139152E-02,.3138904E-02,.3138625E-02,.3138316E-02,.3137975E-02,
$.3137604E-02,.3137201E-02,.3136768E-02,.3136304E-02,.3135809E-02/
data (wp(i),i=1401,1440)/
$.3135283E-02,.3134726E-02,.3134138E-02,.3133519E-02,.3132869E-02,
$.3132189E-02,.3131477E-02,.3130735E-02,.3129962E-02,.3129158E-02,
$.3128323E-02,.3127457E-02,.3126560E-02,.3125633E-02,.3124675E-02,
$.3123686E-02,.3122666E-02,.3121615E-02,.3120534E-02,.3119422E-02,
$.3118279E-02,.3117105E-02,.3115901E-02,.3114666E-02,.3113400E-02,
$.3112103E-02,.3110776E-02,.3109418E-02,.3108029E-02,.3106610E-02,
$.3105160E-02,.3103680E-02,.3102169E-02,.3100627E-02,.3099055E-02,
$.3097452E-02,.3095819E-02,.3094155E-02,.3092461E-02,.3090736E-02/

```

```

data (wp(i),i=1441,1480)/
$.3088981E-02,.3087195E-02,.3085379E-02,.3083532E-02,.3081655E-02,
$.3079748E-02,.3077810E-02,.3075842E-02,.3073843E-02,.3071815E-02,
$.3069756E-02,.3067666E-02,.3065547E-02,.3063397E-02,.3061217E-02,
$.3059007E-02,.3056766E-02,.3054496E-02,.3052195E-02,.3049865E-02,
$.3047504E-02,.3045113E-02,.3042692E-02,.3040242E-02,.3037761E-02,
$.3035250E-02,.3032709E-02,.3030139E-02,.3027538E-02,.3024908E-02,
$.3022248E-02,.3019558E-02,.3016838E-02,.3014089E-02,.3011310E-02,
$.3008501E-02,.3005662E-02,.3002794E-02,.2999896E-02,.2996969E-02/
data (wp(i),i=1481,1520)/
$.2994012E-02,.2991026E-02,.2988010E-02,.2984965E-02,.2981890E-02,
$.2978786E-02,.2975652E-02,.2972489E-02,.2969297E-02,.2966075E-02,
$.2962825E-02,.2959545E-02,.2956236E-02,.2952897E-02,.2949530E-02,
$.2946134E-02,.2942708E-02,.2939254E-02,.2935770E-02,.2932258E-02,
$.2928716E-02,.2925146E-02,.2921547E-02,.2917919E-02,.2914262E-02,
$.2910577E-02,.2906863E-02,.2903120E-02,.2899349E-02,.2895549E-02,
$.2891720E-02,.2887863E-02,.2883978E-02,.2880064E-02,.2876122E-02,
$.2872151E-02,.2868152E-02,.2864125E-02,.2860069E-02,.2855985E-02/
data (wp(i),i=1521,1560)/
$.2851873E-02,.2847734E-02,.2843565E-02,.2839369E-02,.2835145E-02,
$.2830893E-02,.2826613E-02,.2822305E-02,.2817970E-02,.2813606E-02,
$.2809215E-02,.2804796E-02,.2800350E-02,.2795875E-02,.2791374E-02,
$.2786844E-02,.2782288E-02,.2777704E-02,.2773092E-02,.2768453E-02,
$.2763787E-02,.2759094E-02,.2754373E-02,.2749625E-02,.2744850E-02,
$.2740048E-02,.2735219E-02,.2730363E-02,.2725480E-02,.2720571E-02,
$.2715634E-02,.2710671E-02,.2705681E-02,.2700664E-02,.2695621E-02,
$.2690551E-02,.2685454E-02,.2680331E-02,.2675182E-02,.2670006E-02/
data (wp(i),i=1561,1600)/
$.2664804E-02,.2659576E-02,.2654321E-02,.2649040E-02,.2643733E-02,
$.2638400E-02,.2633041E-02,.2627657E-02,.2622246E-02,.2616809E-02,
$.2611347E-02,.2605858E-02,.2600344E-02,.2594805E-02,.2589240E-02,
$.2583649E-02,.2578033E-02,.2572391E-02,.2566724E-02,.2561032E-02,
$.2555315E-02,.2549572E-02,.2543804E-02,.2538011E-02,.2532193E-02,
$.2526350E-02,.2520483E-02,.2514590E-02,.2508672E-02,.2502730E-02,
$.2496763E-02,.2490772E-02,.2484756E-02,.2478715E-02,.2472650E-02,
$.2466561E-02,.2460447E-02,.2454309E-02,.2448147E-02,.2441961E-02/
data (wp(i),i=1601,1640)/
$.2435751E-02,.2429516E-02,.2423258E-02,.2416976E-02,.2410670E-02,
$.2404340E-02,.2397986E-02,.2391609E-02,.2385209E-02,.2378784E-02,
$.2372337E-02,.2365866E-02,.2359371E-02,.2352853E-02,.2346312E-02,
$.2339748E-02,.2333161E-02,.2326551E-02,.2319918E-02,.2313262E-02,
$.2306584E-02,.2299882E-02,.2293158E-02,.2286411E-02,.2279642E-02,
$.2272850E-02,.2266036E-02,.2259200E-02,.2252341E-02,.2245460E-02,
$.2238557E-02,.2231631E-02,.2224684E-02,.2217715E-02,.2210724E-02,
$.2203711E-02,.2196677E-02,.2189620E-02,.2182543E-02,.2175443E-02/
data (wp(i),i=1641,1680)/
$.2168323E-02,.2161180E-02,.2154017E-02,.2146832E-02,.2139626E-02,
$.2132400E-02,.2125152E-02,.2117883E-02,.2110593E-02,.2103282E-02,
$.2095951E-02,.2088599E-02,.2081226E-02,.2073833E-02,.2066420E-02,
$.2058986E-02,.2051531E-02,.2044057E-02,.2036562E-02,.2029047E-02,
$.2021513E-02,.2013958E-02,.2006384E-02,.1998789E-02,.1991175E-02,
$.1983542E-02,.1975888E-02,.1968216E-02,.1960524E-02,.1952812E-02,
$.1945082E-02,.1937332E-02,.1929563E-02,.1921775E-02,.1913968E-02,
$.1906142E-02,.1898298E-02,.1890434E-02,.1882552E-02,.1874652E-02/
data (wp(i),i=1681,1720)/
$.1866733E-02,.1858795E-02,.1850840E-02,.1842866E-02,.1834874E-02,
$.1826863E-02,.1818835E-02,.1810789E-02,.1802725E-02,.1794643E-02,
$.1786544E-02,.1778427E-02,.1770292E-02,.1762140E-02,.1753970E-02,
$.1745784E-02,.1737580E-02,.1729359E-02,.1721120E-02,.1712865E-02,
$.1704593E-02,.1696304E-02,.1687999E-02,.1679677E-02,.1671338E-02,

```

```

$.1662983E-02,.1654611E-02,.1646223E-02,.1637819E-02,.1629399E-02,
$.1620962E-02,.1612510E-02,.1604042E-02,.1595557E-02,.1587058E-02,
$.1578542E-02,.1570011E-02,.1561465E-02,.1552903E-02,.1544325E-02/
data (wp(i),i=1721,1760)/
$.1535733E-02,.1527125E-02,.1518503E-02,.1509865E-02,.1501213E-02,
$.1492545E-02,.1483863E-02,.1475167E-02,.1466456E-02,.1457730E-02,
$.1448990E-02,.1440236E-02,.1431467E-02,.1422684E-02,.1413888E-02,
$.1405077E-02,.1396253E-02,.1387414E-02,.1378563E-02,.1369697E-02,
$.1360818E-02,.1351926E-02,.1343020E-02,.1334101E-02,.1325169E-02,
$.1316224E-02,.1307265E-02,.1298294E-02,.1289310E-02,.1280314E-02,
$.1271305E-02,.1262283E-02,.1253249E-02,.1244202E-02,.1235143E-02,
$.1226072E-02,.1216989E-02,.1207894E-02,.1198787E-02,.1189668E-02/
data (wp(i),i=1761,1800)/
$.1180538E-02,.1171396E-02,.1162242E-02,.1153077E-02,.1143900E-02,
$.1134712E-02,.1125513E-02,.1116303E-02,.1107082E-02,.1097850E-02,
$.1088607E-02,.1079354E-02,.1070089E-02,.1060815E-02,.1051529E-02,
$.1042234E-02,.1032928E-02,.1023612E-02,.1014286E-02,.1004950E-02,
$.9956034E-03,.9862475E-03,.9768819E-03,.9675067E-03,.9581219E-03,
$.9487276E-03,.9393240E-03,.9299112E-03,.9204892E-03,.9110581E-03,
$.9016180E-03,.8921690E-03,.8827112E-03,.8732448E-03,.8637697E-03,
$.8542861E-03,.8447941E-03,.8352937E-03,.8257851E-03,.8162684E-03/
data (wp(i),i=1801,1840)/
$.8067436E-03,.7972109E-03,.7876703E-03,.7781220E-03,.7685659E-03,
$.7590023E-03,.7494312E-03,.7398528E-03,.7302670E-03,.7206740E-03,
$.7110739E-03,.7014668E-03,.6918528E-03,.6822320E-03,.6726045E-03,
$.6629703E-03,.6533295E-03,.6436824E-03,.6340289E-03,.6243691E-03,
$.6147032E-03,.6050312E-03,.5953533E-03,.5856694E-03,.5759799E-03,
$.5662846E-03,.5565837E-03,.5468774E-03,.5371657E-03,.5274486E-03,
$.5177264E-03,.5079991E-03,.4982667E-03,.4885295E-03,.4787874E-03,
$.4690406E-03,.4592892E-03,.4495332E-03,.4397729E-03,.4300082E-03/
data (wp(i),i=1841,1880)/
$.4202392E-03,.4104661E-03,.4006890E-03,.3909079E-03,.3811229E-03,
$.3713342E-03,.3615418E-03,.3517459E-03,.3419465E-03,.3321437E-03,
$.3223377E-03,.3125285E-03,.3027162E-03,.2929009E-03,.2830827E-03,
$.2732617E-03,.2634380E-03,.2536118E-03,.2437830E-03,.2339519E-03,
$.2241184E-03,.2142827E-03,.2044449E-03,.1946051E-03,.1847634E-03,
$.1749198E-03,.1650745E-03,.1552276E-03,.1453792E-03,.1355293E-03,
$.1256781E-03,.1158257E-03,.1059721E-03,.9611747E-04,.8626190E-04,
$.7640548E-04,.6654832E-04,.5669051E-04,.4683217E-04,.3697344E-04/
data (wp(i),i=1881,1883)/
$.2711461E-04,.1725677E-04,.7413338E-05/

```

```

CsTH = cos(TwoTH0 * 0.5/RAD)
if (abs(CsTH) .lt. 1.0e-15) CsTH = 1.0e-15
TTH = sin(TwoTH0 * 0.5/RAD)/CsTH
CsTwoTH = cos(TwoTH0/RAD)
SnTwoTH = sin(TwoTH0/RAD)
ApB = S_L + D_L
AmB = S_L - D_L
ApB2 = ApB**2
if (((S_L .ne. 0.0) .or. (D_L .ne. 0.0)) .and. Use_Asym) then
  tmp = sqrt(1.0 + AmB**2)*CsTwoTH
  if (abs(tmp) .gt. 1.0) then
    Einfl = acos(CsTwoTH)*RAD
  else
    Einfl = acos(tmp)*RAD
  endif
tmp2 = 1.0 + ApB2
tmp = sqrt(tmp2) * CsTwoTH

```

```

c If S_L or D_L are zero, set Einfl = 2theta

      if ((S_L .eq. 0.0) .or. (D_L .eq. 0.0)) Einfl = TwoTH0
      if (abs(tmp) .le. 1.0) then
        Emin = acos(tmp) * RAD
        tmp1 = tmp2 * (1.0 - tmp2 * CsTwoTH**2)
      else
        tmp1 = 0.0
        if (tmp .gt. 0.0) then
          Emin = 0.0
        else
          Emin = 180.0
        endif
      endif
      if ((tmp1 .gt. 0.0) .and. (abs(tmp) .le. 1.0)) then
        dEmindA = -ApB * CsTwoTH/sqrt(tmp1)
      else
        dEmindA = 0.0
      endif
      ArrayNum = 1
      K = 400.0 * (TwoTH0 - Emin) ! Calculate number of terms needed
C From LWF - number of terms must be such that the interval between 2phi(min) and 2theta
C is in steps no larger than 0.005
C      stepsize = (twoth0-emin)/float(K)
      stepsize = 1.0/400.0 ! This seems to be always the case apart from rounding K
      if(stepsize .gt. gamma/10.0) then
        stepsize = gamma/10.0
        K = (twoth0-emin)/stepsize
      endif
      do while ((ArrayNum .lt. 14) .and. (K .gt. NTERMS(ArrayNum)))
        ArrayNum = ArrayNum + 1
      enddo
      NGT = nterms(ArrayNum) ! Save number of terms
      ngt2 = ngt / 2
c Clear terms needed for summations
      sumWG = 0.0
      sumWRG = 0.0
      sumWdGdA = 0.0
      sumWRdGdA = 0.0
      sumWdGdB = 0.0
      sumWRdGdB = 0.0
      sumWGdRd2t = 0.0
      sumWGdRdG = 0.0
      sumWGdRdE = 0.0
      sumWGdRdA = 0.0
      sumWGdRdB = 0.0
c Compute the convolution integral
      it = fstterm(arraynum)-ngt2
      do K = ngt2 , NGT
        delta = Emin + (TwoTH0 - Emin) * xp(K + it)
        dDeltadA = (1.0 - xp(k+it)) * dEmindA
        sinDELTA = sin(Delta/RAD)
        cosDELTA = cos(Delta/RAD)
        if (abs(cosDELTA) .lt. 1.0e-15) cosDELTA = 1.0e-15
        RcosDELTA = 1.0 / cosDELTA
        tanDELTA = tan(Delta/RAD)
        tmp = cosDELTA**2 - CsTwoTH**2
        if (tmp .gt. 0.0) then
          tmp1 = sqrt(tmp)
          F = abs(CsTwoTH) / tmp1
        endif
      enddo

```

```

        dFdA = cosDELTA * CsTwoTH * sinDELTA * dDELTAAdA
1      / (tmp1 * tmp1 * tmp1)
    else
        F = 0.0
        dFdA = 0.0
    endif
c calculate G(Delta,2theta) , FCJ eq. 7a and 7b
    if ( abs(Delta - Emin) .gt. abs(Einfl - Emin)) then
        if (S_L .gt. D_L) then
!
! N.B. this is the only place where d()/dA <> d()/dB
!
            G = 2.0 * D_L * F * RcosDELTA
            dGdA = 2.0 * D_L * RcosDELTA * (dFdA +
1          F*tanDELTA*dDELTAAdA)
            dGdB = dGdA + 2.0 * F * RcosDELTA
        else
            G = 2.0 * S_L * F * RcosDELTA
            dGdB = 2.0 * S_L * RcosDELTA
1          *(dFdA + F * tanDELTA * dDELTAAdA)
            dGdA = dGdB + 2.0 * F * RcosDELTA
        endif
    else
        G = (-1.0 + ApB * F) * RcosDELTA
        dGdA = RcosDELTA * (F - tanDELTA * dDELTAAdA + ApB * dFdA
1      + ApB * F * tanDELTA * dDELTAAdA)
        dGdB = dGdA
    endif
    tmp = PsVoigt(TwoTh-DELTA+TwoTHO,TwoTHO,eta,Gamma,dPRdT
1      ,dPRdG,dPRdE)
    sumWG = sumWG + wp(k+it) * G
    sumWRG = sumWRG + wp(k+it) * G * tmp
    sumWdGdA = sumWdGdA + wp(k+it) * dGdA
    sumWdGdB = sumWdGdB + wp(k+it) * dGdB
    sumWRdGdA = sumWRdGdA + wp(k+it) * dGdA * tmp
    sumWRdGdB = sumWRdGdB + wp(k+it) * dGdB * tmp
    sumWGdRd2t = sumWGdRd2t + wp(k+it) * G * dPRdT
    sumWGdRdG = sumWGdRdG + wp(k+it) * G * dPRdG
    sumWGdRdE = sumWGdRdE + wp(k+it) * G * dPRdE
    sumWGdRdA = sumWGdRdA + wp(k+it) * G * dPRdT * dDELTAAdA * RAD
    enddo
    if (sumWG .eq. 0.0) sumWG = 1.0
    Profval = sumWRG / sumWG
    dPRdT = sumWGdRd2t / sumWG
    dPRdG = sumWGdRdG / sumWG
    dPRdE = sumWGdRdE / sumWG
    dPRdS = (sumWRdGdA + sumWGdRdA) / sumWG - sumWRG *
1      sumWdGdA / sumWG**2

    dPRdD = (sumWRdGdB + sumWGdRdA) / sumWG - sumWRG *
1      sumWdGdB / sumWG**2
    else ! here for no asymmetry }
        tmp = PsVoigt(TwoTH,TwoTHO,eta,Gamma,dPRdT,dPRdG,dPRdE)
        Profval = tmp
        dPRdS = 0.0
        dPRdD = 0.0
    endif
    return
end

```

```

      real*4 function Gauss(Pos , Pos0 , Gamma , dGdT , dGdG )

c  Return value of Gaussian at 'Pos' for peak at 'Pos0' and 'Gamma'.
c  dGdT is derivative of G wrt Pos0.
c  dGdG is derivative of G wrt Gamma.

      implicit none
      real*4 Pos , Pos0 , Gamma , dGdT , dGdG
      real*4  c / 1.6651092/
      real*4  cg / 0.939437279/
      real*4  delp , temp

      delp = Pos - Pos0
      if (abs(delp)/Gamma .gt. 6) then
        Gauss = 0.0
        dGdT = 0.0
        dGdG = 0.0
      else
        temp = cg * exp(-(delp * c /Gamma)**2)/Gamma
        Gauss = temp
        dGdG = temp * ( -1.0 + 2.0 * (delp * c/Gamma)**2) / Gamma
        dGdT = 2.0 * c**2 * delp * temp/Gamma**2
      endif
      return
      end

      real*4 function Lorentz(Pos , Pos0 , Gamma , dLdT , dLdG )

c  Return value of Lorentzian at 'Pos' for peak at 'Pos0' and 'Gamma'.
c  dLdT is derivative of L wrt Pos0.
c  dLdG is derivative of L wrt Gamma.
      implicit none
      real*4 Pos , Pos0 , Gamma , dLdT , dLdG
      real*4  cl / 0.636619772/

      real*4  delp , denom

      delp = Pos - Pos0
      denom = 4.0 * delp**2 + Gamma**2
      Lorentz = cl * Gamma / denom
      dLdT = 8.0 * cl * Gamma * delp / denom**2
      dLdG = cl * (4.0 * delp**2 - Gamma**2) / denom**2
      return
      end

C
      real*4 function PsVoigt(TwoTH , TwoTH0 , eta , Gamma,
1          dPRdT , dPRdG , dPRdE )

c
c  Returns value of Pseudo Voigt
c  Eta is the mixing coefficient between Gaussian and Lorentzian
c  Gamma is the FWHM
c  TwoTH is point at which to evaluate the profile
c  TwoTH0 is two theta value for peak
c  dPRdT is derivative of profile wrt TwoTH0
c  dPRdG is derivative of profile wrt Gamma
c  dPRdE is derivative of profile wrt Eta
      implicit none
      real*4 TwoTH , TwoTH0 , eta , Gamma
      real*4 dPRdT , dPRdG , dPRdE

```

```

real*4  G,Gauss          ! Gaussian part
real*4  L,Lorentz        ! Lorentzian part
real*4  dGdT , dGdG , dLdT , dLdG
G = Gauss(TwoTH , TwoTH0 , Gamma , dGdT , dGdG )
L = Lorentz(TwoTH , TwoTH0 , Gamma , dLdT , dLdG )
PsVoigt = Eta * L + (1.0 - Eta) * G
dPRdT = Eta * dLdT + (1.0 - Eta) * dGdT
dPRdG = Eta * dLdG + (1.0 - Eta) * dGdG
dPRdE = L - G
return
end

```

◇

## 16 One Page Guide

The binning programs are **id31sum** and **id31sumall**. The former combines a set of scans from one specfile into an single rebinned scan, the latter sums up each scan separately.

Fill in the rest of the details here!

```

aplotmtv      : Makes plots
axmgr         : Ditto
bindump       : Part of test suite for binning programs
binit.pdf     : Full source and docs for binning programs
c2xye        : Pulls columns out of specfiles, used in fitting scripts
columns       : Tells you which columns are available in a specfile
cplotdemo     : ?
cyberfit      : Olivier's fitting of cyberscans
fitcol        : Script to fit columns in specfiles
fitem         : Script to fit columns for many scans in specfiles
fitit         : Fitting executable program
foot.mtv      : Needed for plotmtv
guesswork     : Makes eigenvalues from mprodds Pawley style fits
head.mtv      : Needed for plotmtc
id18plot      : ?
id18plot~     : ?
id31check     : Needed for id31sumcheck, this computes the chi^2 stuff
id31offsets   : Produces offsets for the multianalyser
id31sum       : Sums up scans and puts them in an xye file
id31sumall    : Sums up scans individually and puts them in .inp files
id31sumcheck  : Script which makes xye file, the inp files and checks
               : they are all the same
mprodd       : Rietveld program
pld           : script to plot diag.mtv for you (just type pld)
plepf         : script to plot epf file, you need to give it a filename
plmany        : script to plot many x,y,e files
plmesh        : script to plot mesh scans (typically optics MC2G vs pitch)
plmeshcontour : ditto, but with contours instead of lines
plmtv         : plots .mtv files
plotit        : Executable used in plotting scripts
plotmesh      : Executable used in plmesh
plotpr        : Plots Rietveld fits from mprodd
pls           : Script to plot data from specfiles
revi          : script to revise input file for mprodd
save\_exes\_18\_9\_02.tar.gz : exe's before update on 18/9/02
sifit         : Fits wavelength and zero to a set of silicon peaks
sipks         : Fits raw data for silicon to get si peaks and then almbda
sitexp        : Predicts temperature variation of silicon lattice par
smesh         : Plots a series of .inp files as a 3d plot
smeshcont     : ditto but with contours
smeshem       : Produces NeXuS files for lamp, as with smesh
solver        : Olivier's peak fitting routine
temp.mtv      : rubbish
testbins1     : part of binning testsuite
testbins2     : ditto
testnextscan  : exe, needed for id31sumcheck
testspect     : dumps out interesting information about specfiles

```

Some additional help for id31sum is gotten by just typing "id31sum", eg



example: `id31sum file.dat 0.01 1 10`  
 will process the scans 1 to 10 from file.dat with binsize 0.01

Optional arguments:

ed=n1,n2 to exclude detectors n1 and n2 (default=none)  
 es=m1,m2 to exclude scans m1 and m2 (default=none)  
 lowtth=xx.xx to set min two theta to use (default=-30.0)  
 hightth=xx.xx to set max two theta to use (default=160.0)  
 step=x.xxxx to force a stepsize (if very small)  
 wd=xx to set a limit for outliers on diagnostic file, diag.mtv  
 alp=xx for esd=sqrt(cts+alp), default is 0.5  
 mm=xx for to set minimum monitor counts threshold (default=5)  
 scalpk|scalmon|scaltot for pattern scaling (default=scalpk)  
 scalinp=xxx for counts per monitor count (scalmon, default=10000)  
 renorm to use current effics instead of values from temp.res  
 nodiag to prevent creation of diagnostic file, diag.mtv  
 gsas to output a .gsa file for gsas  
 spf to output a .spf file for profil  
 pds to output a .pds file for a pds file

id31sumall:  
 id31sumall version 0.1

example: `id31sumall file.dat 0.01 1 10`  
 will process the scans 1 to 10 from file.dat with binsize 0.01

Optional arguments:

ed=n1,n2 to exclude detectors n1 and n2 (default=none)  
 es=m1,m2 to exclude scans m1 and m2 (default=none)  
 lowtth=xx.xx to set min two theta to use (default=-30.0)  
 hightth=xx.xx to set max two theta to use (default=160.0)  
 step=x.xxxx to force a stepsize (if very small)  
 alp=xx for esd=sqrt(cts+alp), default is 0.5  
 mm=xx for to set minimum monitor counts threshold (default=5)  
 scalinp=xxx to rescale .inp files (default=10000)  
 gsas to output a .gsa files  
 spf to output a .spf files  
 pds to output a .pds files

Offsets are done by id31offsets. Get a good quality data file and run id31offsets in the same way as id31sum. It will make a temp.new file which hopefully has better offsets than before. Do this a few times, until it appears to converge.

Read the file binit.pdf to find out exactly what is happening, or send email an email to [wright@esrf.fr](mailto:wright@esrf.fr)

## 17 Pseudo-Random number generator

Get a re-producible stream of random numbers for doing randomised off-setting of scan start positions.

$\langle \text{random 136a} \rangle \equiv$

```

      real(kind=4) function rlcg( )
! This returns a float in the range 0 to 1 which is a
! deterministic sequence. Will be the same over many runs.
      integer(kind=4) :: a = 69069, c = 1327217885
      integer(kind=4), save :: x = 1
      real(kind=4), save :: m = 2147483647
      x = x*a + c
      rlcg = (real(x)/m+1.0)/2.0
!      write(*,'(I12,I1X,F12.10)')x, rlcg
      end function rlcg

```

◇

Fragment referenced in 28.

## 18 Sorting algorithm from netlib for medians

$\langle \text{dsort 136b} \rangle \equiv$

```

! *DECK DSORT
      SUBROUTINE DISORT (DX, DY, N, KFLAG)
!C***BEGIN PROLOGUE  DSORT
!C***PURPOSE  Sort an array and optionally make the same interchanges in
!C             an auxiliary array.  The array may be sorted in increasing
!C             or decreasing order.  A slightly modified QUICKSORT
!C             algorithm is used.
!C***LIBRARY  SLATEC
!C***CATEGORY  N6A2B
!C***TYPE  DOUBLE PRECISION (SSORT-S, DSORT-D, ISORT-I)
!C***KEYWORDS  SINGLETON QUICKSORT, SORT, SORTING
!C***AUTHOR  Jones, R. E., (SNLA)
!C            Wisniewski, J. A., (SNLA)
!C***DESCRIPTION
!C
!C  DSORT sorts array DX and optionally makes the same interchanges in
!C  array DY.  The array DX may be sorted in increasing order or
!C  decreasing order.  A slightly modified quicksort algorithm is used.
!C
!C  Description of Parameters
!C    DX - array of values to be sorted  (usually abscissas)
!C    DY - array to be (optionally) carried along
!C    N  - number of values in array DX to be sorted
!C    KFLAG - control parameter
!C           = 2  means sort DX in increasing order and carry DY along.
!C           = 1  means sort DX in increasing order (ignoring DY)
!C           = -1 means sort DX in decreasing order (ignoring DY)
!C           = -2 means sort DX in decreasing order and carry DY along.
!C
!C***REFERENCES  R. C. Singleton, Algorithm 347, An efficient algorithm
!C               for sorting with minimal storage, Communications of
!C               the ACM, 12, 3 (1969), pp. 185-187.
!C***ROUTINES CALLED  XERMSG
!C***REVISION HISTORY  (YYMMDD)
!C  761101  DATE WRITTEN

```

```

!C 761118 Modified to use the Singleton quicksort algorithm. (JAW)
!C 890531 Changed all specific intrinsics to generic. (WRB)
!C 890831 Modified array declarations. (WRB)
!C 891009 Removed unreferenced statement labels. (WRB)
!C 891024 Changed category. (WRB)
!C 891024 REVISION DATE from Version 3.2
!C 891214 Prologue converted to Version 4.0 format. (BAB)
!C 900315 CALLs to XERROR changed to CALLs to XERMSG. (THJ)
!C 901012 Declared all variables; changed X,Y to DX,DY; changed
!C         code to parallel SSORT. (M. McClain)
!C 920501 Reformatted the REFERENCES section. (DWL, WRB)
!C 920519 Clarified error messages. (DWL)
!C 920801 Declarations section rebuilt and code restructured to use
!C         IF-THEN-ELSE-ENDIF. (RWC, WRB)
!C***END PROLOGUE  DSORT
!C  .. Scalar Arguments ..
      INTEGER(kind=4) KFLAG
      INTEGER(kind=4) N
!C  .. Array Arguments ..
      DOUBLE PRECISION DX(*)
      INTEGER(kind=4) DY(*)
!C  .. Local Scalars ..
      DOUBLE PRECISION R, T, TT, TTY, TY
      INTEGER I, IJ, J, K, KK, L, M, NN
!C  .. Local Arrays ..
      INTEGER IL(21), IU(21)
!C  .. External Subroutines ..
!C      EXTERNAL XERMSG
!C  .. Intrinsic Functions ..
      INTRINSIC ABS, INT
!C***FIRST EXECUTABLE STATEMENT  DSORT
      NN = N
      IF (NN .LT. 1) THEN
!          CALL XERMSG ('SLATEC', 'DSORT',
!          +          'The number of values to be sorted is not positive.', 1, 1)
          RETURN
      ENDIF
!C
      KK = ABS(KFLAG)
      IF (KK.NE.1 .AND. KK.NE.2) THEN
!          CALL XERMSG ('SLATEC', 'DSORT',
!          +          'The sort control parameter, K, is not 2, 1, -1, or -2.', 2,
!          +          1)
          RETURN
      ENDIF
!C
!C      Alter array DX to get decreasing order if needed
!C
      IF (KFLAG .LE. -1) THEN
          DO 10 I=1,NN
              DX(I) = -DX(I)
10      CONTINUE
      ENDIF
!C
      IF (KK .EQ. 2) GO TO 100
!C
!C      Sort DX only
!C
      M = 1
      I = 1

```

```

      J = NN
      R = 0.375D0
!C
20 IF (I .EQ. J) GO TO 60
   IF (R .LE. 0.5898437D0) THEN
      R = R+3.90625D-2
   ELSE
      R = R-0.21875D0
   ENDIF
!C
30 K = I
!C
!C   Select a central element of the array and save it in location T
!C
      IJ = I + INT((J-I)*R)
      T = DX(IJ)
!C
!C   If first element of array is greater than T, interchange with T
!C
      IF (DX(I) .GT. T) THEN
         DX(IJ) = DX(I)
         DX(I) = T
         T = DX(IJ)
      ENDIF
      L = J
!C
!C   If last element of array is less than than T, interchange with T
!C
      IF (DX(J) .LT. T) THEN
         DX(IJ) = DX(J)
         DX(J) = T
         T = DX(IJ)
!C
!C   If first element of array is greater than T, interchange with T
!C
      IF (DX(I) .GT. T) THEN
         DX(IJ) = DX(I)
         DX(I) = T
         T = DX(IJ)
      ENDIF
!C
!C   Find an element in the second half of the array which is smaller
!C   than T
!C
40 L = L-1
   IF (DX(L) .GT. T) GO TO 40
!C
!C   Find an element in the first half of the array which is greater
!C   than T
!C
50 K = K+1
   IF (DX(K) .LT. T) GO TO 50
!C
!C   Interchange these elements
!C
      IF (K .LE. L) THEN
         TT = DX(L)
         DX(L) = DX(K)
         DX(K) = TT

```

```

        GO TO 40
    ENDIF
!C
!C    Save upper and lower subscripts of the array yet to be sorted
!C
    IF (L-I .GT. J-K) THEN
        IL(M) = I
        IU(M) = L
        I = K
        M = M+1
    ELSE
        IL(M) = K
        IU(M) = J
        J = L
        M = M+1
    ENDIF
    GO TO 70
!C
!C    Begin again on another portion of the unsorted array
!C
    60 M = M-1
        IF (M .EQ. 0) GO TO 190
        I = IL(M)
        J = IU(M)
!C
    70 IF (J-I .GE. 1) GO TO 30
        IF (I .EQ. 1) GO TO 20
        I = I-1
!C
    80 I = I+1
        IF (I .EQ. J) GO TO 60
        T = DX(I+1)
        IF (DX(I) .LE. T) GO TO 80
        K = I
!C
    90 DX(K+1) = DX(K)
        K = K-1
        IF (T .LT. DX(K)) GO TO 90
        DX(K+1) = T
        GO TO 80
!C
!C    Sort DX and carry DY along
!C
    100 M = 1
        I = 1
        J = NN
        R = 0.375D0
!C
    110 IF (I .EQ. J) GO TO 150
        IF (R .LE. 0.5898437D0) THEN
            R = R+3.90625D-2
        ELSE
            R = R-0.21875D0
        ENDIF
!C
    120 K = I
!C
!C    Select a central element of the array and save it in location T
!C
    IJ = I + INT((J-I)*R)

```

```

      T = DX(IJ)
      TY = DY(IJ)
!C
!C      If first element of array is greater than T, interchange with T
!C
      IF (DX(I) .GT. T) THEN
        DX(IJ) = DX(I)
        DX(I) = T
        T = DX(IJ)
        DY(IJ) = DY(I)
        DY(I) = TY
        TY = DY(IJ)
      ENDIF
      L = J
!C
!C      If last element of array is less than T, interchange with T
!C
      IF (DX(J) .LT. T) THEN
        DX(IJ) = DX(J)
        DX(J) = T
        T = DX(IJ)
        DY(IJ) = DY(J)
        DY(J) = TY
        TY = DY(IJ)
!C
!C      If first element of array is greater than T, interchange with T
!C
      IF (DX(I) .GT. T) THEN
        DX(IJ) = DX(I)
        DX(I) = T
        T = DX(IJ)
        DY(IJ) = DY(I)
        DY(I) = TY
        TY = DY(IJ)
      ENDIF
    ENDIF
!C
!C      Find an element in the second half of the array which is smaller
!C      than T
!C
    130 L = L-1
      IF (DX(L) .GT. T) GO TO 130
!C
!C      Find an element in the first half of the array which is greater
!C      than T
!C
    140 K = K+1
      IF (DX(K) .LT. T) GO TO 140
!C
!C      Interchange these elements
!C
      IF (K .LE. L) THEN
        TT = DX(L)
        DX(L) = DX(K)
        DX(K) = TT
        TTY = DY(L)
        DY(L) = DY(K)
        DY(K) = TTY
        GO TO 130
      ENDIF

```

```

!C
!C      Save upper and lower subscripts of the array yet to be sorted
!C
      IF (L-I .GT. J-K) THEN
        IL(M) = I
        IU(M) = L
        I = K
        M = M+1
      ELSE
        IL(M) = K
        IU(M) = J
        J = L
        M = M+1
      ENDIF
      GO TO 160
!C
!C      Begin again on another portion of the unsorted array
!C
      150 M = M-1
        IF (M .EQ. 0) GO TO 190
        I = IL(M)
        J = IU(M)
!C
      160 IF (J-I .GE. 1) GO TO 120
        IF (I .EQ. 1) GO TO 110
        I = I-1
!C
      170 I = I+1
        IF (I .EQ. J) GO TO 150
        T = DX(I+1)
        TY = DY(I+1)
        IF (DX(I) .LE. T) GO TO 170
        K = I
!C
      180 DX(K+1) = DX(K)
        DY(K+1) = DY(K)
        K = K-1
        IF (T .LT. DX(K)) GO TO 180
        DX(K+1) = T
        DY(K+1) = TY
        GO TO 170
!C
!C      Clean up
!C
      190 IF (KFLAG .LE. -1) THEN
        DO 200 I=1,NN
          DX(I) = -DX(I)
      200  CONTINUE
      ENDIF
      RETURN
      END          SUBROUTINE DISORT

```

◇

Fragment referenced in 52.